



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® TDX® SI Heavy Duty Power Wheelchair (Rehab Version) Group 2 Heavy Duty, 450 lb. Weight Capacity

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

POWER CHAIR TYPE

- TDXSI-HD TDX-SI HD with Rehab Seating ... \$5,970.00 K0824
 TDXSISEAT-HD TDX-SI HD Seat Only \$755.00

FRAME COLOR

- 61P Electric Red No Charge
 62P Electric Blue No Charge

USER WEIGHT LIMIT

- U450 User weight <450 lbs. Standard

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00

* Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

SEAT SELECTION

- ADJASBA Adj. Rehab Seat with Adj. Seat Tilt and Back Angle Standard
 SEMIRB Semi-Recline Van Seat with Solid Seat Pan \$70.00
 Seat-to-floor height is 17.5".

ADJUSTABLE ASBA REHAB SEAT

SEAT-TO-FLOOR HEIGHT

Due to variation in options the listed seat-to-floor heights have a tolerance of +/- 0.5"

- LOW 16 1/2" at 0° Standard
 MED 18 1/2" at 0° No Charge
 TALL 20 1/2" at 0° No Charge

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range Standard
 2024W 20" to 24" Seat Width Range No Charge

ASBA REHAB SEAT WIDTH

- W18 18" Seat Width No Charge
 W19 19" Seat Width No Charge
 W20 20" Seat Width Standard
 W21 21" Seat Width \$329.00
 W22 22" Seat Width \$329.00
 W23 23" Seat Width \$329.00
 W24 24" Seat Width \$469.00

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range Standard
 1922D 19" to 22" Seat Depth Range No Charge

ASBA REHAB SEAT DEPTH

- D16 16" Seat Depth No Charge
 D17 17" Seat Depth No Charge
 D18 18" Seat Depth No Charge
 D19 19" Seat Depth No Charge
 D20 20" Seat Depth Standard
 D21 21" Seat Depth \$529.00
 D22 22" Seat Depth \$529.00

ARM TYPE (PAIR)

- Clear
 28F Desk Length Fixed Ht Conv No Charge
 29F Desk Length Adj. Ht Conv \$230.00 E0973
 18F Full Length Fixed Ht Conv Standard
 19F Full Length Adj. Ht Conv \$230.00 E0973

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

**Invacare® TDX® SI Heavy Duty Power Wheelchair (Rehab Version)
Group 2 Heavy Duty, 450 lb. Weight Capacity**

ADJUSTABLE ASBA REHAB SEAT continued

MIXED ARMS (EACH)

- Right Arm Clear
- 19FR Right Arm, Full Length, Adj. Height \$115.00 E0973
- 29FR Right Arm, Desk Length, Adj. Height \$115.00 E0973
- Left Arm Clear
- 19FL Left Arm, Full Length, Adj. Height \$115.00 E0973
- 29FL Left Arm, Desk Length, Adj. Height \$115.00 E0973

SEAT BACK TYPE

- SECT Sectional Rear Uprights Standard

SEAT BACK ANGLE DEGREE

- 90BA 90° Back Angle Standard
- 95BA 95° Back Angle No Charge

ASBA BACK HEIGHT

- BH16 16" Back Height No Charge
- BH17 17" Back Height No Charge
- BH18 18" Back Height No Charge
- BH19 19" Back Height No Charge
- BH20 20" Back Height No Charge
- BH21 21" Back Height No Charge
- BH22 22" Back Height No Charge
- BH23 23" Back Height No Charge
- BH24 24" Back Height No Charge

SEAT POSITIONING STRAPS

- Clear
- 1515 Push Button Style Seat Pos. Strap No Charge
- 1311BK Airline Buckle Seat Pos. Strap \$50.00

CHEST POSITIONING STRAPS

- Clear
- 7311BK Airline Buckle Chest Pos. Strap \$100.00 E0960
- 9515 Push Button Style Chest Pos. Strap \$100.00 E0960
- 7321BK Hook & Loop Chest Pos. Strap \$100.00 E0960

VAN SEAT WITH SOLID SEAT PAN

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery Standard

VAN SEAT SIZE

- VW18 18" Wide Van Seat No Charge
- VD16 16" Seat Depth No Charge
- VD17 17" Seat Depth No Charge
- VD18 18" Seat Depth No Charge
- VW20 20" Wide Van Seat Standard
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth Standard
- VW22 22" Wide Van Seat No Charge
- VD18 18" Seat Depth No Charge
- VD19 19" Seat Depth No Charge
- VD20 20" Seat Depth No Charge

VAN ARM TYPE (PAIR)

- Clear
- Arm selections are flip-up and adjustable width.**
- VAH Full Length, Adj. Height \$230.00 E0973
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adj. Height \$230.00 E0973
- VFHD Desk Length, Fixed Height No Charge

VAN MIXED ARMS

- Arms are flip-up and adjustable width.**
- Right Arm Clear
- RVAHF Right Arm, Full Length, Adj. Height \$115.00 E0973
- RVFHF Right Arm, Full Length, Fixed Height No Charge
- RVAHD Right Arm, Desk Length, Adj. Height \$115.00 E0973
- RVFHD Right Arm, Desk Length, Fixed Height No Charge
- Left Arm Clear
- LVAHF Left Arm, Full Length, Adj. Height \$115.00 E0973
- LVFHF Left Arm, Full Length, Fixed Height No Charge
- LVAHD Left Arm, Desk Length, Adj. Height \$115.00 E0973
- LVFHD Left Arm, Desk Length, Fixed Height No Charge

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

**Invacare® TDX® SI Heavy Duty Power Wheelchair (Rehab Version)
Group 2 Heavy Duty, 450 lb. Weight Capacity**

FRONT RIGGING OPTIONS

CENTER PIVOT STYLE RIGGINGS Clear

- 70NHD 70° Swingaway Footrest 13 1/2" – 17 1/2" Standard
Standard with 1651 footplates.
- 70STAPER 70° Tapered Footrest 13" – 17" No Charge
Standard with Impact Guards and Flip-Up footplates (ATT5543)
- AT5544 Manual Elevating Legrest \$255.00 E0990

PIN STYLE RIGGINGS Clear

- PHW93 Swingaway Footrest 13 1/2" – 19" No Charge K0053
Available with 1651 footplates.
- PHAL4A Smartleg Articulating ELR \$225.00

FOOTBOARD

- 90FB Foot Platform \$150.00

FOOTPLATE OPTIONS Clear

- 1651 Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates \$115.00
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00
- CTMTFB Footboard for Center Mount \$150.00
One piece footboard for Center Mount Riggings (CTMT-F)

FOOTREST/LEGREST ACCESSORIES Clear

- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00
- A1337 Calf Strap \$27.00 K0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops \$45.00 E0951

CENTER MOUNT RIGGINGS

- CTMT-F Center Mount with Adjustable Knee Angle \$300.00
Allows 60 to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2 - 15".
- CTMT-FJ Junior Center Mount with Adj. Knee Angle \$300.00
Allows 60 to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" - 12".
- SM Small Footplates (5 3/4" x 7 1/2") No Charge
- MED Medium Footplates (5 1/2" x 9 1/2") No Charge
- LRG Large Footplates (5.90" x 11 1/2") No Charge
- CALF Calf Pads for Center Mount \$200.00

MISCELLANEOUS ACCESSORIES

INSTALLATION OF BATTERIES

- 22NFBATTERY 22NF Battery(On-Chair)(1 pair) \$300.00 E2361

BATTERY CHARGER OPTIONS

- 110CHARGERCTE 8-Amp Battery Charger Standard
- LESTR Lester 8-Amp Battery Charger \$320.0
- L8900 Omit Charger Less \$80.00

DRIVE WHEEL TYPE

- 1431-3 14" x 3" with Foam Filled Inserts Standard

ELECTRONICS

- MK690 MK6i 90-Amp Controller Standard
- SPJPLUS SPJ+ Single Drive Joystick Standard

JOYSTICK MOUNTING HARDWARE

- QLA Quad Link Retractable Joystick Mount SPJ+ \$279.00 E1028
- LEFT Left Hand Mounted Joystick No Charge

ACCESSORIES

- SEATO2 Oxygen Holder – ASBA Rehab Seat \$150.00 E2208
- VANO2 Oxygen Holder – Van Style Seat \$150.00 E2208
- HFO2-ML6 HomeFill O2 Holder ML6 Only \$150.00
Available for Van Seat only.
- 1813M6 MK6i Hand-Held Programmer \$475.00
- 1303 Wheel Locks \$45.00

Comments

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

©2013 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.