

INVACARE®  
**PERFECTO<sub>2</sub>™ V**  
OXYGEN CONCENTRATOR

---



---

The **Invacare® Perfecto<sub>2</sub>™ V Oxygen Concentrator** is an economical addition to the premier Invacare stationary concentrator line. This 5-liter concentrator is small, quiet and lightweight and comes standard with the **Invacare® SensO<sub>2</sub>® Oxygen Sensor**. The Perfecto<sub>2</sub> V unit incorporates many of the best features of the market-leading Invacare® Perfecto<sub>2</sub>™ Concentrator, including its contemporary, patient preferred external design and color.

---

**Features**

Top handle for ease of mobility

Pressure-based system is superior to timing-based systems

Self-diagnostic electronics minimize troubleshooting time

Alarms and safety systems: Battery-free power-loss alarm; Invacare® Sieve-GARD™ Monitor protects sieve from degradation; Invacare® SensO<sub>2</sub>® Oxygen Monitor; Compressor 35 PSI pressure-relief valve, Low flow alarm

Servicing parts compatibility: Compatible with many Invacare® Platinum™ 10 Oxygen Concentrator parts

Simplified preventative maintenance schedule: Replace preventative maintenance components as needed

**HCPCS Code: E1390**



Invacare® HomeFill® Oxygen System connection

HomeFill compressor and Ready-Rack have been updated to coordinate with the Perfecto<sub>2</sub> Concentrator. Use with Invacare cylinders

Easy to maintain filter



**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
Form No. 09-064  
Rev. 01/12  
111129

Model #	Dimensions	Weight	Shipping Weight	Quantity	Oxygen Flow	Flow Range	Maximum Outlet Pressure	Pressure Relief Mechanism
IRC5PO2V	23" H 13" W 11.5" D	39 lb.	45 lb.	1 Per Case 27 Per Pallet	5 LPM	0.5-5 LPM	5+/- 0.5 PSI	35 PSI +/- 3.5 PSI
	Oxygen Concentration 95.6%-87% At All Flow Rates	Filter Cabinet Compressor Inlet Outlet HEPA	Electrical Requirements 120 VAC +/- 10%, 60 Hz	Rated Current Input 3.0 A	Average Power Consumption 325 W @ 3 LPM	Sound Level 43 db.	Altitude 8,000 ft. 6,000 ft. w/ HomeFill System	Limited Warranty 3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

