

INVACARE®

Insignia™

Lightweight Wheelchair



Yes, you can.®

INVACARE®

Insignia™

Lightweight Wheelchair

The new Invacare® Insignia™ Lightweight Wheelchair will be the lightweight chair of choice for providers and patients alike. Fully adjustable arms, from desk to full length, allow providers to stock one chair instead of two, and standard features include a robust wheel lock, simple adjustable head tube design, sleek caster fork and adjustable angle back. Patients will benefit from having a stylish, lightweight chair that is easy to use and adjusts to their specific needs. The Insignia Wheelchair has it all in one package.



1



Convertible, height adjustable arms

In a single design, the arms offer 4" of height adjustment for maximum flexibility and convert from desk to full length for minimum inventory.

2



Angle adjustable caster journal

The simplicity of this design lends to no loose parts with a 1-2-3 quick and easy adjustment.

- 1) simply loosen the two inner fasteners,
- 2) adjust the two bottom fasteners using the rear edge of the caster fork as your guide,
- 3) re-tighten the inner fasteners.

3



Robust wheel lock

User, therapist and provider preferred, the robust wheel lock is easy to locate and engage with a "click" that confirms it is locked. It is also easy to maintain.

4



Large push button quick release axle

Comfort and convenience.



Standard Features

- Lightweight, aluminum frame
- Quick release axles
- Swingaway footrests
- Push-button, height adjustable arms that offer desk-to-full length changeability
- Height and angle adjustable back includes spring-loaded, flip-up anti-tippers
- Adjustable angle head tubes
- Stock models for quick ship convenience

5



Adjustable angle back

Industry proven bolt-through bracket adjusts from -5° to 20° in 5° increments.

6



Flip-back arm features a large paddle release

Ergonomically correct, the large paddle release is easy to operate and requires minimum dexterity and is removable.

7



Spring-loaded, flip-up anti-tippers

Allow the caregiver to pull and flip without removing for easy navigation over curbs and thresholds.

8



Height adjustable back

The back features a one piece, tool-free height adjustment knob that offers adjustment from 15" – 19" in 1" increments.



Insignia™ Wheelchair Specifications

| Overall | |
|-----------------------------|-----------------------|
| Width Open | 7.50" plus seat width |
| Width Closed | 12" |
| Depth (w/o Riggings) | 31" |
| Seat Widths | 16", 18", 20" |
| Seat Depths | 16", 18" |

Knee-to-Heel* (range adjusts in 1" increments)

| | Min. | Max. |
|------------------|-------------|-------------|
| 60° | 16" | 20" |
| Elevating | 17" | 20.5" |

* Measurement taken from top of seat rail to back of footplate

Arm Type**

| Space-Saver Dual-Point Desk-to-Full Length, Adjustable Height | Min. | Max. |
|--|-------------|-------------|
| | 9" | 13" |

**Measurement taken from top of seat rail to top of arm pad

| | |
|------------------------|--|
| Back Heights | 15" – 19" in 1" increments |
| Weight Capacity | 250 lb. |
| Seat-to-Floor | 19.5", 18.5", 17.5" or 16.5" with 22"/6" wheel kit |

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010

Invacare Corporation.

All rights reserved.

Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Form No. 09-101

Rev. 4/10

For more information
Visit www.invacare.com/insignia

