



**Yes, you can.®**

## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).


**Yes, you can.®**
**TRANSACTION TYPE:**
☐ Order ☐ Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

**Special Note**

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

**INVACARE LEO SCOOTERS**

The below models include base, seat, batteries, and lap belt.

<input type="checkbox"/>	LEO-4S Silver 4-Wheel Mid-size .....	\$2,495.00	K0807
<input type="checkbox"/>	LEO-4B Blue 4-Wheel Mid-size .....	\$2,495.00	K0807
<input type="checkbox"/>	LEO-3S Silver 3-Wheel Mid-size .....	\$2,395.00	
<input type="checkbox"/>	LEO-3B Blue 3-Wheel Mid-size .....	\$2,395.00	

**LEO SCOOTER OPTIONS AND ACCESSORIES**

Options and Accessories are for Leo Scooters only. Off-chair discounts apply to the selections below.

<input type="checkbox"/>	ACC140 Safety Flag .....	\$25.00	
	Only one of these accessories can be used at a time.		
<input type="checkbox"/>	ACC200 Crutch and Cane Holder .....	\$80.00	E2207
	Only one of these accessories can be used at a time.		
<input type="checkbox"/>	ACC210 Walker Holder .....	\$80.00	
	Only one of these accessories can be used at a time.		
<input type="checkbox"/>	ACC220 Oxygen Holder .....	\$80.00	E2208
	Only one of these accessories can be used at a time.		
<input type="checkbox"/>	ACC230 Rear Basket .....	\$80.00	
	Rear Basket dimensions are 18" L x 13" W x 8" H.		
	Only one of these accessories can be used at a time.		
<input type="checkbox"/>	ACC250 Rear View Mirror .....	\$25.00	
<input type="checkbox"/>	ACC260 Lockable Front Box, Silver .....	\$125.00	
	Front Box dimensions are 17" L x 9" W x 11.5" H.		
<input type="checkbox"/>	ACC270 Lockable Front Box, Blue .....	\$125.00	
	Front Box dimensions are 17" L x 9" W x 11.5" H.		
<input type="checkbox"/>	36AH-2 36 Amp Hr Batteries (pair) .....	\$232.00	

**Leo Specifications - Specifications are approximate**

	Leo 3-Wheel	Leo 4-Wheel
<b>Base Length</b>	47"	48"
<b>Base Width</b>	23"	23.5"
<b>Ground Clearance</b>	2.4"	2.5"
<b>Seat Width</b>	19"	19"
<b>Seat Depth</b>	16"	16"
<b>Back Height</b>	18"	18"
<b>Seat to Floorboard Range</b>	16.5" - 19.5"	16.5" - 19.5"
<b>Arm Width Range</b>	18" - 22"	18" - 22"
<b>Weight Limit</b>	350 lbs.	350 lbs.
<b>Speed Range</b>	0 - 5 mph	0 - 5 mph
<b>Incline Capability</b>	6° grade	8° grade
<b>Travel Range</b>	Est. 19 miles	Est. 15 miles
<b>Tires</b>	10" x 2"	10" x 2"
<b>Turning Radius</b>	47"	52"
<b>Weight</b>		
Front Base Section	N/A	60 lbs.
Rear Base Section	N/A	42 lbs.
Batteries	25 lbs. each	25 lbs. each
Seat	32 lbs.	32 lbs.
Total Weight	174 lbs.	184 lbs.
<b>Warranty</b>		
Frame, Electronics, Transaxle	1 year limited	1 year limited
Batteries	6 months	6 months

**The Invacare Corp. Leo Scooter is not intended for highway or street use.**

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).


**Invacare Corporation**

www.invacare.com

**USA**

One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

©2013 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.