



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® TDX® SI-2 Power Wheelchair Group 2, 300 lb. Weight Capacity



Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

STANDARD CONFIGURATIONS

For the best value, order the TDXSI-2 in Standard Configurations.
Each Standard Configuration is pre-configured with 22NF Gel Batteries, Adjustable Height Arms, Flat Free Tires, and a Footboard.
 Just choose the seat size and the frame color.

GROUP 2 - K0823 VAN SEAT

- TDXSIV-2-S \$6,500.00 ^{K0823/}_{K0011}
 Standard Features and Suggested List Prices:
 • Semi-Recline Van Seat, Gray Vinyl \$5970.00
 • 22NF Gel Batteries Installed on Chair \$300.00 E2361
 • Full Length, Adjustable Height Arms \$230.00 E0973
 • 14" Flat Free Tires
 • Foot Platform
 • Push-button Style Seat Positioning Strap
 • 2-Pole Motors
 • MK5 75 Amp Controller
 • 8 Amp Off-Board Charger

FRAME COLOR

- 61P Electric Red No Charge
 62P Electric Blue No Charge

VAN SEAT SIZE

- VW18 18" Wide Van Seat No Charge
 VD16 16" Seat Depth No Charge
 VD17 17" Seat Depth No Charge
 VD18 18" Seat Depth No Charge
 VW20 20" Wide Van Seat No Charge
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth No Charge

FOOTREST/FRONT RIGGINGS

- Clear
 90FB Foot Platform Standard
 70NHD Swingaway Footrests \$165.00
 AT5544 Manual Elevating Legrest \$255.00 E0990

GROUP 2 - K0822 VAN SEAT with SOLID PAN

- TDXSI-2-S \$6,500.00 ^{K0822/}_{K0011}
 Standard Features and Suggested List Prices:
 • Semi-Recline Van Seat with Solid Seat, Gray Vinyl \$5970.00
 • 22NF Gel Batteries Installed on Chair \$300.00 E2361
 • Full Length, Adjustable Height Arms \$230.00 E0973
 • 14" Flat Free Tires
 • Foot Platform
 • Push-button Style Seat Positioning Strap
 • 2-Pole Motors
 • MK5 75 Amp Controller
 • 8 Amp Off-Board Charger

FRAME COLOR

- 61P Electric Red No Charge
 62P Electric Blue No Charge

VAN SEAT SIZE

- VW18 18" Wide Van Seat No Charge
 VD16 16" Seat Depth No Charge
 VD17 17" Seat Depth No Charge
 VD18 18" Seat Depth No Charge
 VW20 20" Wide Van Seat No Charge
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth No Charge

FOOTREST/FRONT RIGGINGS

- Clear
 90FB Foot Platform Standard
 70NHD Swingaway Footrests \$165.00
 AT5544 Manual Elevating Legrest \$255.00 E0990

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare TDX SI-2 Power Wheelchair Group 2, 300 lbs. Weight Capacity

CUSTOM CONFIGURATIONS

For clients with more advanced seating and positioning needs, custom configurations allow you to customize the TDXSI-2 to meet the specific need of each client.

CUSTOM CONFIG - FRAME SELECTION

POWER CHAIR TYPE

- TDXSIV-2 TDX-SI with Van Seat \$5,970.00^{K0823/}_{K0011}
 TDXSISEAT TDX-SI Seat Only \$755.00

FRAME COLOR

- 61P Electric Red No Charge
 62P Electric Blue No Charge

USER WEIGHT LIMIT

- U300 User weight <300 lbs. Standard

* Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

VAN SEAT SELECTIONS

SEAT SELECTION

- SEMIR Semi-Recline Van Seat Standard

VAN UPHOLSTERY COLORS

- GV Gray Vinyl Upholstery Standard

VAN SEAT SIZE

- VW16 16" Wide Van Seat No Charge
 VD16 16" Seat Depth No Charge
 VD17 17" Seat Depth No Charge
 VD18 18" Seat Depth No Charge
 VW18 18" Wide Van Seat No Charge
 VD16 16" Seat Depth No Charge
 VD17 17" Seat Depth No Charge
 VD18 18" Seat Depth No Charge
 VW20 20" Wide Van Seat Standard
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth Standard
 VW22 22" Wide Van Seat No Charge
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth No Charge

VAN ARM TYPE (PAIR)

Clear

Arms are flip-up and adjustable width.

- VAH Full Length, Adj. Height \$230.00 E0973
 VFH Full Length, Fixed Height Standard
 VAHD Desk Length, Adj. Height \$230.00 E0973
 VFHD Desk Length, Fixed Height No Charge

VAN MIXED ARMS

Arms are flip-up and adjustable width.

- Right Arm Clear
 RVAHF Right Arm, Full Length, Adj. Height \$115.00 E0973
 RVFHF Right Arm, Full Length, Fixed Height No Charge
 RVAHD Right Arm, Desk Length, Adj. Height \$115.00 E0973
 RVFHD Right Arm, Desk Length, Fixed Height No Charge
 Left Arm Clear
 LVAHF Left Arm, Full Length, Adj. Height \$115.00 E0973
 LVFHF Left Arm, Full Length, Fixed Height No Charge
 LVAHD Left Arm, Desk Length, Adj. Height \$115.00 E0973
 LVFHD Left Arm, Desk Length, Fixed Height No Charge

SEAT POSITIONING STRAPS

Clear

- 1515 Push Button Style Seat Pos. Strap Standard
 1311BK Airline Buckle Seat Pos. Strap \$50.00

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare TDX SI-2 Power Wheelchair

Group 2, 300 lbs. Weight Capacity

MISCELLANEOUS ACCESSORIES

INSTALLATION OF BATTERIES

22NFBATTERY 22NF Battery(On-Chair)(1 pair) \$300.00 E2361

BATTERY CHARGER OPTIONS

- 110CHARGER CTE 8-Amp Battery Charger Standard
 LESTR Lester 8-Amp Battery Charger \$320.00
 L8900 Omit Charger Less \$80.00

DRIVE WHEEL TYPE

1431-3 14" x 3" with Foam Filled Inserts Standard

ELECTRONICS

- MK5 NX 75-Amp Controller Standard
 SPJPLUS SPJ+ Single Drive Joystick Standard

JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retractable Joystick Mount SPJ+ \$230.00 E1028
 LEFT Left Hand Mounted Joystick No Charge

ACCESSORIES

- VANO2 Oxygen Holder \$150.00 E2208
 HFO2-ML6 HomeFill O2 Holder ML6 Only \$150.00
 1813M6 MK6i Hand-Held Programmer \$475.00
 1303 Wheel Locks \$45.00

FRONT RIGGING OPTIONS

FOOTBOARD

90FB Foot Platform Standard

CENTER PIVOT STYLE RIGGINGS

- 70NHD 70° Swingaway Footrest 13 1/2" – 17 1/2" \$165.00
Standard with 1651 footplates.
 70STAPER 70° Tapered Footrest 13" – 17" \$165.00
Standard with impact guards (IMPACTG) and flip-up footplates (AT5543).
 AT5544 Manual Elevating Legrest \$255.00 E0990

PIN STYLE RIGGINGS

- PHW93 Swingaway Footrest 13 1/2" – 19" \$165.00
Standard with 1651 footplates.
 PHAL4A Smartleg Articulating ELR \$225.00 K0053

FOOTPLATE OPTIONS

- 1651 Composite Footplates No Charge
 1350 Extra Large Aluminum Footplates \$115.00
 AT5543 Adjustable Angle Flip-Up Footplates \$165.00
 CTMTFB Footboard for Center Mount \$150.00
One piece footboard for Center Mount Riggings (CTMT-F)

FOOTREST/LEGREST ACCESSORIES

- Clear
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
 ALPT Longer Pivot And Slide Tube \$84.00
 ALPT4 Longer Pivot And Slide Tube 4" \$84.00
 A1337 Calf Strap \$27.00 K0038
 IMPACTG Impact Guards \$51.00
 HLOOP Heel Loops (pair) \$45.00 E0951

CENTER MOUNT RIGGINGS

- CTMT-F Center Mount with Adjustable Knee Angle \$300.00
Allows 60 to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2" - 15".
 CTMT-FJ Junior Center Mount with Adj. Knee Angle \$300.00
Allows 60 to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" - 12".
 SM Small Footplates (5 3/4" x 7 1/2") No Charge
 MED Medium Footplates (5 1/2" x 9 1/2") No Charge
 LRG Large Footplates (5.90" x 11 1/2") No Charge
 CALF Calf Pads for Center Mount \$200.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare TDX SI-2 Power Wheelchair Group 2, 300 lbs. Weight Capacity

Specification	Van Seat	Van Seat w/ Solid Seat Pan
Weight capacity	300 lbs.	300 lbs.
Maximum speed	5 MPH	5 MPH
Incline capability	6°	6°
Per-charge range	Up to 13 mi.	Up to 13 mi.
Turning radius	23"	23"
Ground clearance	3"	3"
Motor Package	2-Pole	2-Pole
Standard Electronics	MK5 NX/SPJ+	MK5 NX/SPJ+
Overall Dimensions		
Base length w/o front riggings	35"	35"
Base length w/ front riggings	38 1/2"	38 1/2"
Height	38"	38"
with headrest	45" - 47"	45" - 47"
Seat Dimensions		
Width	16" - 22"	16" - 20"
Depth	16" - 20"	16" - 20"
Back height	18"	18"
with headrest	25" - 27"	25" - 27"
Back angle	85° - 115°	85° - 115°

Specification	Van Seat	Van Seat w/ Solid Seat Pan
Seat-to-floor height		
Low	21"	17.5"
Weight		
With batteries	267 lbs.	267 lbs.
Without batteries	191 lbs.	191 lbs.
With transport brackets	Add 10 lbs.	Add 10 lbs.
Battery options	22NF	22NF
Wheel and caster size		
Drive wheels	14" x 3" foam filled	14" x 3" foam filled
Front/rear casters	6" x 2"	6" x 2"
Limited warranties		
Base frame	60 months	60 months
Seat frame	36 months	36 months
Electronics	12 months	12 months
Motor/gearbox	12 months	12 months

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

©2013 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.