INVACARE® TOILET SAFETY





Invacare Toilet Safety Products, including raised toilet seats with and without arms and toilet safety frames, can help make the bathroom a safer environment.

Invacare® Raised Toilet Seat Features

Model # 1300RTS

Recommended for use with Toilet Safety Frame

Attaches without tools and fits on most round bowl configurations*

Rear seat wings for more secure attachment

Contoured seat for internal rotation of hips

Ample depth of front recess makes hygiene performance more convenient

HCPCS Code: E0243





- 1. Tool-less assembly
- 2. Retail packaging



Invacare® Clamp-On Raised Toilet Seat Features

Model # 1301RTS without arms

Clamp-On mechanism tightens easily with one hand

Comfortable, contoured seat is made of easy-to-clean, durable molded composite

Recommended for use with toilet safety frame

Fits on most round bowl configurations*

HCPCS Code: E0243





- 1. Tool-less assembly
- 2. Retail packaging



^{*} ALWAYS make sure that the Raised Toilet Seat is securely in place and stable before each use. Although the Raised Toilet Seat is designed for use on most popular toilets, it may not work on certain unusual, elongated or artistic bowls.

Invacare® Clamp-On Raised Toilet Seat Features

Model # 1302RTS with arms

Comfortable, contoured seat is made of easy-to-clean, durable molded composite

Recommended for use with toilet safety frame

Fits on most round bowl configurations*

New adjustable locking arms

HCPCS Code: E0243





- 1. Tool-less assembly
- 2. Retail packaging



Invacare® Toilet Safety Frame Features

Model #1392KD

Arms fold back for cleaning and are easily removed for transfers

Armrests provide secure hand support

Frame is made of anodized, non-corrosive 1" aluminum tubing

Aluminum bracket is easy to attach

HCPCS Code: E0243





- 1. Tool-less assembly
- 2. Retail packaging



^{*} ALWAYS make sure that the Raised Toilet Seat is securely in place and stable before each use. Although the Raised Toilet Seat is designed for use on most popular toilets, it may not work on certain unusual, elongated or artistic bowls.















Model #	Dimensions	Dimensions	s Seat	Seat-to Arm	Arm-to-Floor	Weight	Weight
	With Arms	Between Arr	ms Dimensions	Height	Height		Capacity
1300RTS	N/A	N/A	5"x14.75"x15.5"D	N/A	N/A	3 lb. 6 oz.	300 lb.
1301RTS	N/A	N/A	5"x16.1"x16.5"D	N/A	N/A	4.5 lb.	300 lb.
1302RTS	21.5"	19"	5"x16.1"x16.5"D	12"	26.5"	5.5 lb.	300 lb.
1392KD	N/A	20"-22"	N/A	N/A	25.5"-29.5"	2 lb. 5 oz.	250 lb.





Model #	Units	Limited	
	Per Carton	Warranty	
1300RTS	4	3 Years	
1301RTS	2	3 Years	
1302RTS	2	3 Years	
1392KD	2	3 Years	

Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2017

Invacare Corporation.

All rights reserved.

Trademarks are identified by the symbols ™ and ®.

All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form No. 10-152

Rev. 04/17

170543

