

INVACARE® BATH BOARD



The **Invacare® Bath Board** provides a foundation for all of your daily hygiene needs. The lightweight bath board features an easy to clean surface and stable handle for secure transfer, all while looking great in your bathroom.

Features

Model # H112-5

Attractively designed bathtub board

Non-slip seat surface

Hand grip provides additional safety

Practical shower head holder on both sides

Individually adjustable using adjustment elements



1. Non-slip seat
2. Handle
3. Shower head holder
4. Tool-less assembly

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 10-189
Rev. 03/12
111272



| Model # | Overall Dimensions | Width Handgrip | Adjustment Range, Inner | Product Weight | Weight Capacity | Limited Warranty |
|---------|---------------------------|----------------|-------------------------|----------------|-----------------|------------------|
| H112-5 | 27" D 13" W 6.25" H | 7" | 15"– 24.75" | 3.5 lb. | 330 lb. | 2 Years |