



INVACARE®
SEATING
FOR MANUAL WHEELCHAIRS



As the premier brand in the home health care industry, Invacare continues to be committed to meeting both the needs of providers and patients worldwide with Invacare's standard wheelchair products. Invacare general use seat cushions and backs are available on standard wheelchairs for all users, even without a qualifying diagnosis. Just Ask For It!



1. The **Invacare® Matrix® Genera™ Back (GB)** provides support and comfort at a price so affordable that no wheelchair user needs to make do with sling upholstery.
2. For individuals with mild to moderate seating needs, the **Invacare® Absolute™ Cushion (EC)** offers simple, but highly functional solutions. This cushion is designed to provide basic comfort using highly resilient foam and is intended for individuals who are able to do their own weight shift for pressure relief.
3. For individuals with mild to moderate seating needs, the **Invacare® Single Density Cushion (SD)** offers simple, but highly functional solutions. This cushion is designed to provide basic comfort using highly resilient foam for pediatrics to geriatrics who are able to do their own weight shift for pressure relief.
4. For individuals with mild to moderate seating needs, the **Invacare® Comfort-Mate Extra Cushion (CMEX)** goes one step further with a comfortable, molded design and a durable moisture-resistant coating. Ideal for individuals who are able to do their own weight shifts for pressure relief, who need comfort, stability, gentle positioning and some protection.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-128
Rev. 03/12
111249



Backs



Model #	Height	Width	Depth Adjustment	Angle Adjustment	Weight Capacity	Limited Warranty
GB	16", 20"	16", 18", 20"	N/A	N/A	300 lb.	1 Year

Cushions



Model #	Height	Width	Depth	Weight	Limited Warranty
EC	3.25" Adductor 2.25" Abductor 2" Rear	16", 18", 20"	16", 18"	2 lb.	1 Year
SD	2" Rear 1.5" Rear Center 2.5" Front	10"-24"	7"-20"	2 lb.	1 Year
CMEX	2.625" Adductor 3" Abductor 2.25 Leg Trough	14"-22"	14"-20"	2.3 lb.	1 Year