




INVACARE® BED RAILS



All **Invacare® Bed Rails** feature durable chrome-plated, welded steel construction and are designed to withstand heavy-duty use. These rails can be installed permanently or semi-permanently on most spring-fabric medical-style beds. When properly used, bed rails are some of the most important safety and accident prevention accessories used in the patient room today.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.




1. **Invacare® Reduced Gap Full-length Bed Rail**

HCPSC Code: E0310

Model # 6629

Durable 7/8" chrome-plated, welded steel construction

Self-adjusting spring-loaded crossbraces snap easily into place without tools for quick installation



2. **ProBasics® Reduced Gap Full-Length Rails**

HCPSC Code: E0310


Model # PB6034

Welded 7/8" steel construction with brown vein finish

Self adjusting spring loaded crossbraces

Four crossbars help reduce entrapment hazards

Compatible with VC5310, 5307IVC, 5310IVC and 5410IVC



3. **ProBasics® Clamp-On Half Rail**

HCPSC Code: E0305


Model # PB7035

Clamp on rail, no tools required

Welded steel construction with brown vein finish and chrome on the second and fifth rail for scratch-resistant movement of the adjustment bracket

Compatible with VC5210*, VC5310, 5307IVC, 5310IVC and 5410IVC

*Only compatible when used with casters



4. **Invacare® Assist Bed Rail**

HCPSC Code: E0305


Model # 6632

Heavy-duty 1" powder coated, welded steel construction

3 positions: Guard, Assist, Transfer

Moves easily out of the way with a quick release button

Compatible with 5307IVC, 5310IVC, 5410IVC and 5410LOW



5. **Invacare® Reduced Gap Heavy-Duty Half-Length Bed Rail**

HCPSC Code: E0305

Model # BAR6640IVC

Heavy-duty 1" chrome-plated, welded steel construction

Bolt-on hardware for added strength

For optimal safety and convenience, four half-length rails are recommended (one set for head area and one set for foot area)

Compatible with BAR600IVC bed only

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2015
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ®.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
ProBasics is a
registered trademark
of Garden City Medical
Inc. Specifications
are subject to change
without notification.
Form No. 11-174
Rev. 08/15
1506109

Model #	Rail Length*	Rail Height Above Deck Raised	Rail Height Above Deck Lowered	Width Between Rails	Shipping Weight	Limited Warranty
6629	55"	14.5"	4.25"	36"	27.5 lb.	1 Year
PB6034	55"	14"	4.5"	36"	27 lb.	2 Year
PB7035	37"	16"	5"	36"	23 lb.	2 Year
6632	27"	14.5"	N/A	36"	21 lb.	1 Year
BAR6640IVC	28.75"	15"	3.5"	42"	25 lb.	1 Year

*Measurement on flat bed with crossbraces installed as per instruction sheet.

Invacare and ProBasics beds, mattresses and rails, excluding powered support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

NOTE: Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

