

INVACARE®
HOMEFILL®
OXYGEN SYSTEM



The **Invacare® HomeFill® Oxygen System** has revolutionized ambulatory oxygen by allowing patients to fill their own cylinders from their Invacare 5-liter and 10-liter concentrators. With the Invacare HomeFill Oxygen System, providers reduce costs by virtually eliminating cylinder deliveries, while patients enjoy their freedom and independence with an unlimited supply of ambulatory oxygen. The HomeFill Oxygen System allows the patient to breathe at home from a continuous flow concentrator and then use a convenient ambulatory pack when they are “out and about.”

Features

Provides an unlimited refillable ambulatory oxygen supply

Patient breathes from a continuous flow oxygen concentrator at night

Patient protection ensures oxygen purity in cylinder is always greater than 90%

Easy cylinder connections with the Invacare® OxyLock™ Coupler

Provider Benefits

Significantly reduces operating costs by virtually eliminating cylinder deliveries, without sacrificing patient care

No extra equipment needed because patients can breathe off concentrator while filling cylinders

Encourages compliance with oxygen therapy as patients have no worries about running out of oxygen

System is compatible with 5 liter or 10 liter concentrators and a wide variety of cylinder options


Patient Benefits

Eliminating oxygen deliveries increases independence, convenience and comfort

Small lightweight cylinders allow greater mobility than larger, traditional oxygen cylinders

Using small, refillable cylinders indoors can help avoid trips or falls associated with using a concentrator and long oxygen tubing

HCPCS Code: K0738



1. The M2 cylinder is comparable, in size, to a water bottle
2. Ready Rack available
3. Compatible with Invacare® concentrators
4. Invacare® OxyLock™ Coupler

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900



Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2013

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-189
Rev 04/13
130341



Model #	Dimensions	Weight	Inlet Pressure Rating	Oxygen Output	Oxygen Input	Electrical Requirements	Rated Current Input
IOH200	15" H 20.25" W 16" D	33 lb.	14-21 PSI	0-3 LPM (5 liters) 0-5 LPM (10 liter)	0-2 LPM	120 VAC ± 10% @ 60 Hz	2 A
							
	Power Consumption	Limited Warranty					
	140 W	3 Years/ 5,000 hrs.					

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.