# INVACARE® PLATINUM™ 10 OXYGEN CONCENTRATOR





The Invacare® Platinum™ 10 Oxygen Concentrator provides an attractive alternative to liquid systems and is known for reliability, performance and value. The 10 liter unit features quiet operation, a top handle and bottom recesses for easy handling and a side filter access door. In addition, the concentrator is compatible with the Invacare® HomeFill® Oxygen System and is available with an Invacare® SensO2® Oxygen Purity Sensor.

# **Features**

Pressure-based system is superior to timing-based systems

Self diagnostic electronics minimize troubleshooting time

Top handle for ease of mobility

Battery-free power loss alarm, Invacare® Sieve-GARD™ Monitor protects sieve from degradation, compressor 35 psi pressure-relief valve and low flow alarm

# **Provider Benefits**

Provides continuous flow up to 10 liters per minute as an alternative to liquid systems

All Platinum 10 stationary concentrators include a HomeFill system port for inventory flexibility

Patient-accessible reset button helps minimize emergency service calls

# **Patient Benefits**

User friendly design with easy to see control panel, humidifier bottle holder and filter access

HomeFill system compatibility allows patients with flow rates up to 5LPM to enjoy life's experiences with small, refillable cylinders

Reliable concentrator provides unlimited, worry-free oxygen supply

# **HCPCS Code: E1390**







- 1. Invacare® HomeFill® Oxygen System connection
- 2. Easy to maintain filter
- 3. Also available with Invacare® SensO2® Oxygen Purity Sensor



# **Invacare Corporation**

www.invacare.com

# **USA**

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

# Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

©2013 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-190 Rev. 04/13 130376



Dimensions



Weight

54 lb.



Shipping

Weight

60 lb.



Quantity

1 Per Case

18 Per Pallet



Oxygen

Flow

10 LPM



Flow

Range

2-10 LPM





Maximum Pressure Relief Outlet Pressure Mechanism 9+/- 0.5 PSI 35 PSI +/- 3.5 PSI

IRC10LX 26 <sup>3</sup>/<sub>8</sub>" H IRC10LXO2 18 3/8 W 14 <sup>3</sup>/<sub>8</sub>" D

Model#













Sound Level	Limited Warranty
50 db.	3 Years
Average	



Oxygen

Purity

94%-87%

Rates



Filter Electrical Requirements Cabinet 120 VAC At All Flow Compressor Inlet, +/- 10%, 60 Hz

Outlet HEPA (132VAC/108VAC)

Rated Average Power Current Input Consumption 4.0 A 585 W



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.