

INVACARE® TRAPEZE EQUIPMENT



The **Invacare® Trapeze** unit is an important patient room accessory designed to help the users change positions while in bed and to aid in transferring from bed to chair with easy attendant assistance. Trapeze units have a wide range of height adjustments and handbar positions to maximize patient accessibility.

Invacare® Trapeze Floor Stand

Model # 7714P

Provides free-standing support for model # 7740P fixed trapeze unit
Folds flat for convenient storage and handling
Safe working load of 650 lb.

Invacare® Fixed Offset Trapeze Bar with Trapeze

Model # 7740P

Heavy-gauge, octagon steel construction; no tools required for easy assembly
Wall bumpers and padded mounting brackets protect bed ends and wall surfaces
For use on model # 7714P floor stand or on any Invacare Homecare bed head section
Safe working load of 650 lb.

Invacare® Bariatric Floor Stand with Trapeze

Model # BARTRAP

Three locking overhead boom angle positions for comfortable, secure use
Tool-free design simplifies assembly
Four piece sub-assembly significantly reduces setup time
Lightweight design for easy delivery
Free-standing structure permits use with Bariatric, Home Care or non-medical bed
Safe working load of 650 lb.

HCPCS Code: E0940/E0910



7714P: shown with 7740P
7740P: Includes bar, brackets and grab bar

HCPCS Code: E0912



BARTRAP

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2020
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-217
Rev. 8/20
200646



Model #	Height	Base	Reach	Diameter	Product Weight	Shipping Weight	Limited Warranty
7714P	38"	40" L 32" W	N/A	1.25"	24 lb.	N/A	3 Years
7740P	61"	N/A	34"	N/A	17 lb.	N/A	3 Years
BARTRAP	73"	50.5" L 40" W	35.5"	2.25"	72 lb.	80 lb.	3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

