

INVACARE®
PREVENTION OVERLAYS
THERAPEUTIC SUPPORT SURFACES



Invacare offers two static overlays to aid in the prevention of pressure ulcers for those patients who are at risk of developing a pressure ulcer. The **Invacare® Gel Overlay** and **Invacare® CareGuard™ Alternating Pressure Pad System** are designed to be placed overtop an existing mattress. The overlays are both easy to use and maintain.

Invacare® Gel Foam Mattress Overlay

24 horizontal bladders minimize the migration of use
 One piece construction requires no assembly
 Strong elastic corner straps secure overlay to mattress
 Meets California Technical Bulletin #117 for fire retardancy
 Non-toxic gel used

Invacare® CareGuard™ APP Alternating Pressure Pad System

Pump alternately inflates and deflates air cells
 Variable pressure setting for added patient comfort
 Pump can be mounted on the bed with built-in bracket
 Pad includes flaps that wrap around mattress for better stability
 Does not contain latex

HCPCS Code: E0185



HCPCS Code: E0181



Invacare Corporation
 www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2012
 Invacare Corporation.
 All rights reserved.
 Trademarks are identified
 by the symbols ™ and
 ®. All trademarks are
 owned by or licensed
 to Invacare Corporation
 unless otherwise noted.
 Specifications are subject
 to change without
 notification.
 Form No. 11-220
 111226



Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
IVCGFM02	79" L 34" W 3.5" H	275 lb.	16 CFR 1632 16 CFR 1633	6 Months (Mattress/ Cover)
CG9701	118" L 34" W 2.5" H	250 lb.	N/A	30 Days (Mattress) 1 Year (Power Unit)

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

