


INVACARE® CANES



Invacare® Canes are stylish and unique, and with a wide assortment of handgrip styles, you can fit virtually any of your clients with ease. Invacare® Quad Canes come in medium and large base varieties, have a 250 lb. weight capacity and allow for proper fit with right- or left-hand use. Invacare® Bariatric Canes are built with durable steel construction that provides maximum stability for individuals weighing up to 700 lb. All Invacare Canes are height adjustable for proper fit for the consumer and are designed with a twist locking ring which provides additional security during height adjustment.



1. The Invacare® Grip

The familiar triangular-shaped grip. The flat top helps reduce the stress a user's hand may experience through everyday cane use



2. Cane Tip and Collar

Locking adjustment collar secures cane at desired length. Tip features reinforced metal insert for longer wear

3. Cane Tips

Model # 40850-4, 40851-4, 40852-4, 40853-4
Each pair of cane tips is retail packaged

HCPCS Code: E0100 Standard Canes E0105 Quad Canes



1. Invacare® Standard Crook Canes

Model # 8916 – Black
Invacare Grip

2. Invacare® Offset Canes

Model # 40918-6 – Black
Model # 8917 – Silver
Invacare Grip with Strap

3. Medium Base Quad Cane

Model # 3920-2 – Medium Base Invacare Grip
Model # 7827 – Large Base Invacare Grip

Invacare® Orthopedic Grip Cane

4. Model # 8929 – Black (Right Hand)

This unique orthopedic grip is ergonomically shaped to spread pressure evenly across the palm

5. Bariatric Quad Canes

Model # 7831-2 – Small Base, 7835-2 – Large Base
700 lb. weight capacity
Right- or left-hand use

6. Bariatric Offset Cane

Model # 8940-6
700 lb. weight capacity
Height adjustable

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-235
111263



Model #	Adjustment Range	Fits Cane Diameter	Product Weight	Weight Capacity	Patient Height Range	Quantity Per Carton	Limited Warranty
8916	30"-39"	N/A	15 oz.	250 lb.	4'6"-6'6"	6	3 Years
40918-6	29"-38"	N/A	15 oz.	250 lb.	4'4"-6'4"	6	3 Years
8917	29"-38"	N/A	15 oz.	250 lb.	4'4"-6'4"	6	3 Years
3920-2	29.25"-38.25"	N/A	1 lb. 13 oz.	250 lb.	N/A	2	3 Years
7827	29.25"-38.5"	N/A	2 lb. 10 oz.	250 lb.	N/A	2	3 Years
8929	31"-41"	N/A	15 oz.	250 lb.	4'6"-6'6"	6	3 Years
7831-2	29.5"-38.5"	N/A	3 lb.	700 lb.	4'11"-6'5"	2	3 Years
7835-2	29.5"-38.5"	N/A	3 lb.	700 lb.	4'11"-6'5"	2	3 Years
8940-6	29"-38"	N/A	3 lb.	700 lb.	4'10"-6'4"	6	3 Years
40850-4	N/A	1"	N/A	N/A	N/A	8	N/A
40851-4	N/A	.875"	N/A	N/A	N/A	8	N/A
40852-4	N/A	.75"	N/A	N/A	N/A	8	N/A
40853-4	N/A	.625"	N/A	N/A	N/A	8	N/A

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

