

# INVACARE® CRUTCHES



**Invacare® Quick-Change Crutches** offer a solid one-piece crutch frame in lightweight anodized aluminum with push-button height adjustment and contoured handgrips. **Invacare® Quick-Adjust Crutches'** one-piece handgrip assembly eliminates the need for wing nuts and other assembly pieces. Plus, the push-button height adjustment and patient height guide allow the crutch to be fitted to the proper user height. **Invacare® Forearm Crutches** feature 1" solid composite cuffs, cushioned vinyl handgrips and extra large crutch tips. Strong and comfortable, **Invacare® Bariatric Crutches** range in size from adult to tall adult. Large, double push buttons allow for height adjustments in 1" increments.

---

**Invacare® Quick-Change Crutch Features**

Model # 8115-T: Tall Adult, 8115-A: Adult, 8115-J: Junior

Invacare quick-change crutches offer a solid one-piece crutch frame in lightweight anodized aluminum with push-button height adjustment and contoured handgrips

---

**HCPCS Code: E0114**



---

**Invacare® Quick-Adjust Crutches Features**

Model # 8120-T: Tall Adult, 8120-A: Adult, 8120-J: Junior

The one-piece handgrip assembly eliminates the need for wing nuts and other assembly pieces. Plus, the push-button height adjustment and patient height guide allow the crutch to be fitted to the proper user height. This patented model comes completely assembled for the consumer's convenience

---

**HCPCS Code: E0114**



---

## Invacare® Forearm Crutch Features

Model # 8153-T: Tall, 8153-A: Adult, 8153 -J: Junior

Forearm crutches feature 1" solid composite cuffs, cushioned vinyl handgrips and extra large crutch tips

---

**HCPCS Code: E0110**




### 1. Rubber Tips

Model # 6134 Rubber Tips – Large 1.75" O.D. Base  
Model # 6135 Rubber Tips – Extra Large 2" O.D. Base  
Help grip the surface and absorb shock  
Reinforced with metal inserts for longer-lasting wear



### 2. Arm Pads

Model # 6131 Arm Pads  
Add comfort to armrests



### 3. Handgrips

Model # 6125, 6126 Solid Handgrips  
Padded to keep hands and wrists comfortable  
Contoured grips that limit rotation

## Invacare® Bariatric Crutch Features

Model # 8130-A: Adult Bariatric, 8130-T: Tall Bariatric

Strong and comfortable, Invacare's bariatric crutches range in size from adult to tall adult. Large, double push buttons allow for height adjustments in 1" increments. Crutch frame offers durable steel construction

1000 lb. weight capacity

Push-button height adjustment

Patient height guide on crutch for easy adjustment

Contoured underarm design

Arm pads (7131), handgrips (7126) and 2" crutch tips (6135) also available as accessories


### HCPCS Code: E0114



#### 1. Rubber Tips

Model # 6135

Extra large, 2" outside diameter (O.D.) base



#### 2. Arm Pads

Model # 7131

Contoured underarm design



#### 3. Handgrips

Model # 7126

Solid handgrips

### Invacare Corporation

www.invacare.com

#### USA

One Invacare Way

Elyria, Ohio

44035-4190

(800) 333-6900

#### Canada

570 Matheson Blvd. E.,

Unit 8

Mississauga, Ontario

L4Z 4G4 Canada

(800) 668-5324

©2012

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols ™ and ®.

All trademarks are

owned by or licensed

to Invacare Corporation

unless otherwise noted.

Specifications are subject

to change without

notification.

Form No. 11-236

111264



Model #	Adjustment Range	Product Weight	Weight Capacity	Patient Height Range	Quantity Per Carton	Limited Warranty
8115-J	36"-44"	3 lb. 4 oz.	350 lb.	4'6"-5'2"	8 Pairs	3 Years
8115-A	44"-52"	3 lb. 8 oz.	350 lb.	5'5"-5'10"	8 Pairs	3 Years
8115-T	52"-60"	3 lb. 10 oz.	350 lb.	5'10"-6'6"	8 Pairs	3 Years
8120-J	36"-44"	3 lb. 7 oz.	250 lb.	4'10"-6'4"	8 Pairs	3 Years
8120-A	44"-52"	3 lb. 11 oz.	250 lb.	5'2"-5'10"	8 Pairs	3 Years
8120-T	52"-60"	4 lb.	250 lb.	5'10"-6'6"	8 Pairs	3 Years
8153-J	23"-32"	4 lb. 8 oz.	300 lb.	4'5"-5'2"	6 Pairs	3 Years
8153-A	29"-38"	4 lb. 8 oz.	300 lb.	5'2"-5'11"	6 Pairs	3 Years
8153-T	32"-41"	4 lb. 8 oz.	300 lb.	5'11"-6'8"	6 Pairs	3 Years
8130-A	44"-52"	4 lb. 6 oz.	1,000 lb.	5'2"-5'10"	1 Pair	3 Years
8130-T	52"-60"	5 lb.	1,000 lb.	5'10"-6'6"	1 Pair	3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

