

INVACARE® ASPIRATOR



The **Invacare® Aspirator** is a durable suction unit designed for most home care suction requirements.

Features

- High-performance vacuum specifications
- Adjustable vacuums from 0" to 23" Hg
- Durable piston vacuum pump
- Easy-to-read gauge
- Non-marring suction cup feet
- Single-use collection jars and connecting tubing for quick and easy replacement (1 set included)
- UL and CSA listed

HCPCS Code: E0600



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

- Suction cup feet
- Disposable collection jar, Model # IRC1140 (HCPCS Code: A7002)
- Tube and filter kit (not shown) Model # 1148804

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-252
111290



Model #	Dimensions	Weight	Oxygen flow	Vacuum range	Electrical Requirements	Limited Warranty
IRC1135	14.8" L 9.8"H 6.8" W	12.1 lb.	0-40 LPM	0-520mm Hg 0-23" Hg	115 V, 60Hz 1/10 hp Pole Motor	1 Year

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

