

INVACARE®  
**TDX® SC**  
POWER WHEELCHAIR



The **Invacare® TDX® SC Power Wheelchair** is an extension of the TDX SP family by offering the distinguishing features of Enhanced Invacare® SureStep® Suspension with Stability Lock, contemporary aesthetics and wider rear caster spacing for easy interfacing with aftermarket docking devices.


## Features

Invacare® Formula™ CG Powered Tilt and Powered Elevating Seat available  
 Adjustable rehab (ASBA) seat  
 Captain's seat available in 16", 18", 20" and 22" widths with a 21" seat-to-floor height  
 True center wheel drive for intuitive driving  
 Transport Ready Option (TRRO) available\*

## Invacare Core Technologies

Invacare® SureStep® Suspension with Stability Lock provides a smooth and stable transition  
 Invacare® G-Trac™ Technology provides precise tracking for imprecise environments  
 Invacare® MK6i™ Electronics are the benchmark in innovative technology  
 Motion Concepts power positioning systems can be interfaced with the Invacare TDX SC power base. Contact Motion Concepts for more information at 888-433-6818

## HCPCS Code: K0823, K0835, K0848



1. Enhanced Invacare® SureStep® Suspension with Stability Lock
2. Invacare® G-Trac™ Technology
3. Invacare® MK6i™ Electronics
4. Motion Concepts Powered Positioning System Compatibility

\*TRRO includes four factory-installed brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

**Invacare Corporation**  
 www.invacare.com

### USA

One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 333-6900

### Canada

570 Matheson Blvd. E.,  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

©2012

Invacare Corporation.  
 All rights reserved.  
 Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
 Form No. 12-238  
 120338

Model #	Overall Dimensions**	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
TDXSC, TDXSC2-CG, TDXSCV	39" L 25.8" W 36.5" H	16"-22"x 16"-20"	16.5"-21"	10.5"	3"	20.1"	9°
	Total Weight	Weight Capacity	Motor	Max Speed	Battery	Battery Range	Limited Warranty
	Up to 286 lb.	250 lb.	2 Pole	4.8 mph	U1	Up To 12.67 Miles	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/Motor/Gearbox

\*\*Rehab seat and 90° manual centermount  
 Do not operate on roads, streets or highways

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

