

INVACARE®



PRONTO®  
AIR


PERSONAL  
TRANSPORTER



INVACARE INVACARE

## Features | Highlights

The Invacare® Pronto® Air Personal Transporter redefines pure power mobility for the consumer. It is a culmination of years of refining Invacare's approach to mobility, with the penultimate delivery of innovation and user experience.



- (1) Lush, contoured Invacare® MyBody™ Seating
- (2) Ergonomic, modular armrests for natural rest and activity
- (3) Sleek 12" Drive Wheels with 4.5 mph top speed
- (4) Intuitive driver controls
- (5) Travel-ready design



## Innovative | Design

Inspiration is everywhere in the innovative design, but it is especially evident in **Invacare® MyBody® Seating**. Invacare designed MyBody Seating to promote better posture through contoured, pressure relieving foam and natural points of rest which may result in less fatigue at the end of a fulfilling, accomplished day.



## Technical Specifications

### Typical Product Parameters

All dimensions are  $\pm .50"$  unless otherwise indicated

#### Overall Dimensions

Base Width (without joystick)	23.5"
Base Length	37.8" (w/ footplate)
Overall Height	Base w/ seat: 42"-46" Base only: 15"-19"

#### Seat

Seat-to-Floor Height	21"-25" in 1" increments
Seat to Footplate	18"-22" (w/ footplate 3" from ground) 17"-21" (w/ footplate 4" from ground)
Seat Depth	18"-20"
Back Height	21" (without headrest)
Headrest	7" Height adjustment
Back Angle Adjustment	24° of Recline (infinite adjustment angle)
Armrest Adjustment	Width 3" (each arm) Height 3" (each arm) Angle 38° Depth 2.5" (each arm)

#### Wheels

Caster	6"x2" Urethane with Precision Sealed Bearings
Forks	Two sided
Drive Wheel/Tires	12"x2.5" Non-marking Urethane

#### Battery

Battery Requirements	17 AH x2
----------------------	----------

#### Footrests

Footboard	Height 3", 4" from floor Depth .5", 1", 1.5" from chassis
-----------	--

#### Weight

Weight (base without batteries)	75 lb.
Seat Weight (with Arms)	55 lb.
Batteries	30 lb. total (15 lb. each)
*Maximum Weight Limitations	Up to 250 lb.

\*Weight limitation is total weight (user weight plus any additional items that the user may require [backpack, etc]). Example: If weight limitation of the personal transporter is 250 lb. and additional items equals 25 lb., subtract 25 lb. from 250 lb., this means the maximum weight limitation of the user is 225lb.



**Invacare Corporation**  
www.invacare.com

#### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

#### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols <sup>TM</sup> and <sup>®</sup>. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
Form No. 12-489

