






INVACARE® MATRIX®

ELITE

(12"-20" HIGH MODELS) BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS								OVERALL DIMENSIONS (CONTACT AREA)									
	A WIDTH		B HEIGHT		C LATERAL CONTOUR HEIGHT		D HEIGHT		E INNER MAXIMUM WIDTH		F OUTER MAXIMUM WIDTH		G HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
PBE1212	9¾	25	12	31	5¾	14.5	2	5	8½	21.5	11	28	14	36	2.0	0.9	300	136
PBE1214	9¾	25	13½	34.5	8	20	1¼	3	8½	21.5	11	28	15	38	2.2	1.0	300	136
PBE1216	9¾	25	16	41	6½	16.5	5¾	14.5	8½	21.5	11	28	17	43	2.3	1.0	300	136
PBE1218	9¾	25	17	43	8¼	21	1½	4	8½	21.5	11	28	18	46	2.5	1.1	300	136
PBE1220	9¾	25	19	48	8¼	21	1½	4	8½	21.5	11	28	20	51	2.7	1.2	300	136
PBE1312	10¼	26.5	11¾	30	7	18	3	7.5	9½	24	12	31	14	36	2.2	1.0	300	136
PBE1314	10¼	26.5	13½	34.5	7	18	3	7.5	9½	24	12	31	15	38	2.3	1.0	300	136
PBE1316	10¼	26.5	16	41	8	20	3¼	8.5	9½	24	12	31	17	43	2.5	1.1	300	136
PBE1318	10¼	26.5	17	43	8	20	1½	4	9½	24	12	31	18	46	2.7	1.2	300	136
PBE1320	10¼	26.5	19	48	8	20	1½	4	9½	24	12	31	20	51	2.9	1.3	300	136
PBE1412	12½	32	12	31	5¼	13.5	4½	11.5	11½	29	14	36	13	33	2.4	1.1	300	136
PBE1414	12½	32	13½	34.5	7½	19	2¼	6	11½	29	14	36	15	38	2.6	1.2	300	136
PBE1416	12½	32	16	41	9¼	23.5	2¾	7	11½	29	14	36	17	43	2.8	1.3	300	136
PBE1418	12½	32	18	46	6	15.5	5¼	13.5	11½	29	14	36	19	48	3.0	1.4	300	136
PBE1420	12½	32	20	51	7¼	18.5	4¾	12	11½	29	14	36	21	53	3.2	1.5	300	136



INVACARE® MATRX®

ELITE

(12"-20" HIGH MODELS) BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS								OVERALL DIMENSIONS (CONTACT AREA)											
	A		B		C		D		E		F		G		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY			
	WIDTH	HEIGHT	WIDTH	HEIGHT	LATERAL CONTOUR HEIGHT	LATERAL CONTOUR HEIGHT	HEIGHT	HEIGHT	INNER MAXIMUM WIDTH	INNER MAXIMUM WIDTH	OUTER MAXIMUM WIDTH	OUTER MAXIMUM WIDTH	HEIGHT	HEIGHT	lb	kg	lb	kg		
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
PBE1512	13½	34.5	12	31	6	15.5	3½	9	12½	32	15	38	13	33	2.6	1.2	300	136		
PBE1516	13½	34.5	16	41	6	15.5	5¾	14.5	12½	32	15	38	17	43	3.0	1.4	300	136		
PBE1518	13½	34.5	18	46	6	15.5	5¼	13.5	12½	32	15	38	19	48	3.2	1.5	300	136		
PBE1520	13½	34.5	20	51	7¼	18.5	4¾	12	12½	32	15	38	21	53	3.4	1.5	300	136		
PBE1612	14½	37	12	31	6	15.5	3½	9	13½	34.5	16	41	13	33	2.8	1.3	300	136		
PBE1616	14½	37	16	41	6	15.5	5¾	14.5	13½	34.5	16	41	17	43	3.2	1.5	300	136		
PBE1618	14½	37	18	46	6	15.5	5¼	13.5	13½	34.5	16	41	19	48	3.4	1.5	300	136		
PBE1620	14½	37	20	51	7¼	18.5	4¾	12	13½	34.5	16	41	21	53	3.6	1.6	300	136		
PBE1712	15½	39.5	12	31	6	15.5	3½	9	14½	37	17	43	13	33	3.0	1.4	300	136		
PBE1716	15½	39.5	16	41	6	15.5	5¾	14.5	14½	37	17	43	17	43	3.3	1.5	300	136		
PBE1718	15½	39.5	18	46	6	15.5	5¼	13.5	14½	37	17	43	19	48	3.6	1.6	300	136		
PBE1720	15½	39.5	20	51	7¼	18.5	4¾	12	14½	37	17	43	21	53	3.8	1.7	300	136		
PBE1812	16½	42	12	31	6	15.5	3½	9	15½	39.5	18	46	13	33	3.2	1.5	300	136		
PBE1816	16½	42	16	41	6	15.5	5¾	14.5	15½	39.5	18	46	17	43	3.4	1.5	300	136		
PBE1818	16½	42	18	46	6¼	16	4¾	12	15½	39.5	18	46	19	48	3.7	1.7	300	136		
PBE1820	16½	42	20	51	7¼	18.5	4¾	12	15½	39.5	18	46	21	53	4.0	1.8	300	136		
PBE1912	17½	44.5	12	31	6	15.5	3½	9	16½	42	19	48	13	33	3.3	1.5	300	136		
PBE1916	17½	44.5	16	41	6	15.5	5½	13.5	16½	42	19	48	17	43	3.6	1.6	300	136		
PBE1918	17½	44.5	18	46	7¼	18.5	4¾	12	16½	42	19	48	19	48	3.9	1.8	300	136		
PBE1920	17½	44.5	20	51	7¾	18.5	4¾	12	16½	42	19	48	21	53	4.2	1.9	300	136		
PBE2012	18½	47	12	31	6	15.5	3½	9	17½	44.5	20	51	13	33	3.5	1.6	300	136		
PBE2016	18½	47	16	41	6	15.5	5¾	14.5	17½	44.5	20	51	17	43	3.8	1.7	300	136		
PBE2018	18½	47	18	46	7¼	18.5	4¾	12	17½	44.5	20	51	19	48	4.1	1.9	300	136		
PBE2020	18½	47	20	51	7¼	18.5	4¾	12	17½	44.5	20	51	21	53	4.4	2.0	300	136		

Elite Back Specifications Page 2 of 2

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form 12-571 Rev. 06/14 140442