

INVACARE®
CAREGUARD™
ALTERNATING PRESSURE PAD SYSTEM



Invacare® CareGuard™ Alternating Pressure Pad (APP) System is designed to help prevent pressure injuries for those patients who are at risk of developing a pressure injury. The APP System is designed to be placed overtop an existing mattress and is both easy to use and maintain.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Invacare® CareGuard™ Alternating Pressure Pad (APP) System

- Pump alternately inflates and deflates air cells
- Variable pressure setting for added patient comfort
- Pump can be mounted on the bed with built-in bracket
- Pad includes flaps that wrap around mattress for better stability
- Does not contain latex

HCPCS Code: E0181*



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2019
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ©.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 19-001
190101



Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
CG9701	80" L 34" W 2.5" H	250 lb.	N/A	30 Days (Mattress) 1 Year (Power Unit)

*All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

