

Unbeatable



Yes, you can.



Design. Technology. Performance.

AVIVA™ FX POWER WHEELCHAIR

Enjoy the ride

Invacare® AVIVA™ FX Power Wheelchairs are designed to help you do what you want with confidence, in style and in comfort. Whether navigating tight corners at home, out and about in unpredictable conditions or at work, you will be in control in the driver's seat...and enjoying our most comfortable front-wheel drive ride yet.



AVIVA™ FX chassis

Good to go



The **AVIVA FX** Power Wheelchair puts the power in power wheelchair. It is designed to handle slopes, curbs and uneven surfaces while you go about your day. It looks as good as it performs, too. You can personalize your **AVIVA FX** Power Wheelchair to match your style.

> Patent-pending Invacare® 4Sure™ Suspension

When it comes to traction, there's no substitute for all four wheels in contact with the ground at all times. That's why we designed the advanced 4Sure Suspension System with articulating front anti-tippers, shock absorbers and springs all working together. But all you need to know is that you remain in control, and in comfort, when you need it most.

> Invacare® G-TRAC® Technology

Patented digital gyroscopic technology comes standard on the **AVIVA FX** Power Wheelchair. It detects and corrects even the smallest deviation from your intended path, so you'll go exactly where you want to go without having to tweak the control, even at higher speeds. It's especially handy if you use any specialty controls with your wheelchair.

Navigating everyday life

Stability even on uneven surfaces

Because all four wheels are designed to always stay in contact with the ground, you'll enjoy a smooth, comfortable ride. And thanks to an advanced suspension system, obstacles aren't an obstacle for the **AVIVA FX** Power Wheelchair.

> Curbs happen

The 4Sure Suspension System is engineered to make short work of them so you can get to where you want to be. And because the best obstacle is neither felt nor heard, elastomer dampeners minimize annoying squeaking of springs in the 4Sure Suspension System.

> Tackle tight turns

The small turning radius and your position at the front of the chair helps you turn the corner in cramped spaces, ideal for getting around at home.

> Speed it up

If you feel the need for speed, the **AVIVA FX** Power Wheelchair is available with an optional high speed 7.5 mph motor.





Get closer to what matters

The **AVIVA FX** Power Wheelchair boasts a seat-to-floor height as low as only 16.75", even with a full-sized 14" drive wheel and powerful Group 24 batteries.

This means that not only do you enjoy improved traction and stability, but you can also get closer to tables...and what really matters.



Express yourself

It's yours. Make it that way

The sleek and modern design of the **AVIVA FX** Power Wheelchair is your canvas. Make it yours by personalizing the shrouds, fenders, casters and wheel hubs with vibrant, contemporary colors.

Mix shroud/fender and wheel colors for a unique combination that is distinctly you. And be prepared for the compliments you'll be receiving while you're out and about in your new **AVIVA FX** Power Wheelchair.

Available colors

	Obsidian Black		Ruby Red
	Denim Navy		Peacock Blue
	Purple Passion		Lime Twist
	Morning Mist		Orange Marmalade
	Pink Flamingo		Banana Split
	Arctic Blue		

Express yourself

Light up the night

See what you need to when you're out at night thanks to powerful optional LED lights, signals and tail lights. Because your day doesn't end when the sun goes down.





Elevate your expectations

Elevate up to 12”

for natural eye-level interaction or reaching for items on shelves.

> Discreet and dynamic

Thanks to a unique elevation mechanism, you can adjust the height of your seating discreetly, whenever you want.

> Moves with you (and your companion)

There's no need to stop before adjusting your height, the simply smart **AVIVA FX** Power Wheelchair allows you to raise or lower yourself while driving. Plus, the stability that the 4Sure Suspension System offers means you can drive at up to 3.2 mph while elevated so you can interact with others at a natural level while on the move.

The seating system

Ultra Low Maxx

We've paired the **AVIVA FX** Power Wheelchair's base with Ultra Low Maxx Power Positioning- a high end seating system that helps you feel comfortable in your chair.



170°

recline

50°

tilt

> Natural pivot points

Natural pivot points is the industry speak for our Ultra Low Maxx Power Positioning System's design. It means your head, shoulders, hips and feet stay in contact with the chair, mimicking your body's natural pivot points. It's a snug fit that will leave you feeling comfortable in your wheelchair.

> Supported at all times

No matter what dynamic positioning movements the seating system makes, you will remain comfortably supported at all times as the Ultra Low Maxx moves precisely with you.

> Built-in rail

This allows for accessories to be added to your power wheelchair depending on your needs and wants.

> A custom fit

Width: 16"-22"

Depth: 16"-23"

Backrest height: 14"-20"

Armrest height: 9"-18.75"



REM400 Remote

3.5"

color touchscreen

Built-in Bluetooth®

Connect to your smart phone, laptop and smart speaker

Direct Access

Jump straight into the menus - you have direct access to the profiles and functions

Touchscreen technology

Swipe or tap our color touchscreen and be in control of where you go with our REM400 remote

Ambidextrous

The controls are adaptable for both left and right-handed use. If a standard joystick doesn't work for you, you can customize your wheelchair with a specialty control

Clear speed gauge

On the top left of the control screen is a speedometer, indicating the speed of your chair as you are driving it. An odometer is also present to record the distance you have traveled





See for yourself

Driving is believing

Not all power wheelchairs drive alike. Take the **AVIVA FX** Power Wheelchair for a test drive and experience the performance, technology and design that we can only describe in one word. Unbeatable.

Ask your provider for more details.



Connected = Confident

> MyLiNX™ App

The MyLiNX App communicates with your **AVIVA FX** Power Wheelchair so you can see useful information on your smart phone about your power wheelchair, including battery charge remaining so you can go about your day in confidence.

For even more peace of mind and convenience, diagnostic data about the **AVIVA FX** Power Wheelchair is uploaded to your participating provider via the cloud for remote assistance in troubleshooting issues.

Unbeatable



Design. Technology. Performance.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: (905) 695-0138

INVACARE.COM/UNBEATABLE

©2021
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject to
change without notification.
Bluetooth is a registered
trademark of Bluetooth
SIG Inc.
Form No. 20-010 rev6/1
200103

