



**AQUATEC® OCEAN ERGO
ACCESSORIES AND ADJUSTMENTS**



Yes, you can.®

AQUATEC®

Introduction

This useful guide explains how the Ocean Ergo Family can be customized to suit a wide range of physical conditions by making simple adjustments or by adding some easy to fit, off-the-shelf accessories.

Meet the Ocean Ergo Family...



Ocean 24" Ergo XL ▲

An entry level, plus-size model with a maximum user weight of 400 lbs.

Ocean Ergo XL ▲

An entry level plus-size attendant propelled model, with a maximum user weight of 400 lbs.

Ocean VIP Ergo ▲

Designed for individuals who require slightly more assistance with postural control when seated

Ocean Ergo ▲

An entry level attendant propelled model

Ocean Dual VIP Ergo ▲

Designed for individuals with complex seating and positioning needs

Ocean 24" Ergo ▲

An entry level self-propelling model

Accessories

Seating and backrests ▼



Standard soft seat

Model # AP1603142

- ▶ With key shaped cut out.



Soft seat insert

Model # AP1603143

- ▶ Can be used with or without the soft seat when commode function is not required.



Lateral supports

Model # 1535077

- ▶ For safe positioning. Width and height adjustable – also for small users and young adults.



Small soft seat

Model # AP1603145


- ▶ With small, oval opening for smaller and young adults.



Small soft seat insert

Model # AP1603152

- ▶ Can be used with the small soft seat when the commode function is not required.



Soft backrest cushion

Model # 1558769

- ▶ For additional stability and comfort. Will also benefit smaller users and young adults.



Comfort soft backrest

Model # AP1603606 VIP Ergo

Model # AP1604136 Ocean Ergo

Model # AP1604016 Ergo XL

- ▶ Easy to clean Aquatec Ocean Comfort Soft Backrest for Ocean Ergo VIP Shower and Commode Chair model provides users with additional comfort.



Positioning ▼



Calf and amputee support

Model # AP1643573, Right

Model # AP1643574, Left

- ▶ Padded, both height and angle adjustable. For use with the standard plastic seat, standard soft seat and small soft seat.



Calf support with foot support

Model # AP1643543, Right

Model # AP1643572, Left

- ▶ Padded, both height and angle adjustable with additional foot support.



Self propelling wheel kit

Model # 16379

- ▶ Self propel wheel with quick release axle.

Safety ▼



Ocean Safety belts

Chest Model # 1470081

Pelvic Model # 1470082

- ▶ For safe positioning. Individually adjustable, high load capacity and easy to assemble. Machine washable.



Armrest lock

Model # 1470077

- ▶ Locks armrest into the lower or higher position.



Front and rear anti-tippers

Ergo VIP Model # AP1603992

Ergo Others Model # AP1601986

- ▶ For additional safety.

Footrests

- ▶ Height adjustable, swing away and detachable footrests.
- ▶ No tools required.

Seat Height

- ▶ Seat to ground height can be adjusted between 18.675" - 23.6"
- ▶ No tools required.

Toileting ▼



Splash guard

Model # AP1603155

- ▶ Clips on easily for use with the standard plastic seat, standard soft seat and small soft seat.



Toilet pan

Model # 10230

- ▶ Ergonomic handle allows caregiver to remove pan from holder easily.
- ▶ Scale inside the pan.
- ▶ Lid holder when lid is not required.



Pail Support Guide Rail

Model # AP1602078

- ▶ Ergonomic handle allows carer to remove pan from holder easily.
- ▶ Scale inside the pan.
- ▶ Lid holder when lid is not required.

Footplates

- ▶ Angle adjustable, flip-up and textured footplates with heel loops.
- ▶ No tools required.




Adjustments

All Ocean Ergo models can be adjusted to suit a range of user sizes.

Armrests

- ▶ Ergonomic, flip-up armrests can be adjusted to different heights.



Compatibility

- ✓ - Compatible
- ✗ - Not compatible
- - Included

Most accessories are cross compatible with all Ocean Ergo models.

	Product code	Ocean Ergo	Ocean 24" Ergo	Ocean Ergo XL	Ocean 24" Ergo XL	Ocean VIP Ergo	Ocean Dual VIP Ergo
Standard soft seat	AP1603142	✓	✓	✓	✓	✓	✓
Soft seat insert	AP1603143	✓	✓	✓	✓	✓	✓
Small soft seat	AP1603145	✓	✓	✓	✓	✓	✓
Small soft seat insert	AP1603152	✓	✓	✓	✓	✓	✓
Splash guard	AP1603155	✓	✓	✓	✓	✓	✓
Chest belt	1470081	✓	✓	✓	✓	✓	✓
Pelvic belt	1470082	✓	✓	✓	✓	✓	✓
Armrest locking device	1470077	✓	✓	✓	✓	✓	✓
Soft backrest cushion	1558769	✓	✓	✓	✓	✓	✓
Thigh support, Right	AP1643573	✓	✓	✓	✓	✓	✓
Thigh support, Left	AP1643574	✓	✓	✓	✓	✓	✓
Calf support with foot support, Right	AP1643543	✓	✓	✓	✓	✓	✓
Calf support with foot support, Left	AP1643572	✓	✓	✓	✓	✓	✓
24" self propel wheels (set)	16379	✓	-	✗	-	✗	✗
Lateral	1535077	✓	✓	✓	✓	✓	✓
Toilet pan	10230	-	-	-	-	-	-
Toilet pan holder	AP1602078	-	-	-	-	-	-
Front and back anti-tipper Ergo VIP models	AP1603992	✗	✗	✗	✗	✓	✓
Front anti-tipper Ocean Ergo models	AP1601986	✓	✓	✓	✓	✗	✗
Comfort soft backrest VIP Ergo models	AP1603606	✗	✗	✗	✗	✓	✓
Comfort soft backrest Ocean Ergo models	AP1604136	✓	✓	✗	✗	✗	✗
Comfort soft backrest Ocean Ergo XL models	AP1604016	✗	✗	✓	✓	✗	✗

*All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E., Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2020 Invacare Corporation.
All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted or its subsidiaries. Aquatec is owned by Invacare Int'l GmbH. Specifications are subject to change without notification.
Form No. 19-142 rev 0320 200318

