

The Invacare® AX3™ folding wheelchair raises the bar for what a lightweight wheelchair should offer. With standard options in abundance, the economical AX3 still maintains lightweight, performance driven quality while offering stability and independence for any consumer. The versatility is apparent in the numerous options available, but it is the functionality that makes it a winner. The AX3: easy to transport, easy to adjust, easy to function.



Standard models include:

- Square rear frame
- One-piece side frames
- Swing-in/swing-out front riggings
- Flip-down hideaway wheel locks
- 20 different rear wheel settings, multiple center-of-gravity (COG) positions

Other options include:

- Adjustable tension back/seat upholstery
- Push-button height adjustable stroller handles
- Elevating leg rests



Wheel-locks:
Versatile wheel-locks conveniently hideaway for easy safe transfers.



Axle plate:
Dependable axle plate offers no-nonsense rear wheel adjustability.



Optional adjustable stroller handles:
Easily adapt to any caregiver's ideal height preference.

Specifications

	Seat width (inches)	16", 18", 20"		Seat depth	16", 18"
	Overall width (inches)	Seat width plus approx. 8.5"		Backrest angle degree	0 - 30°
	Width of folded wheelchair (inches)	Approx. 13"		Knee to heel height (inches)	13.75" - 19"
	Total height (inches)	32" - 42.5"		Adjustable armrest height (inches)	6.5" - 10.5"
	Rear seat-to-floor (inches)	14" - 19"		Rear wheel diameter (inches)	22" or 24"
	Front seat-to-floor (inches)	16" - 20"		User's maximum weight (lb.)	275



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2007 Invacare Corporation.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Form No. 06-119



REHAB
AX3™

