



For 3D Image Systems: E-mail Order Form and Digitized File to:

ElyriaCustomSeating@invacare.com

Company Name:		Date:	
Account #:		Contact Name:	
Ship To Address:		E-mail Address:	
City/State/Zip:		Phone Number:	
Client Reference:		Purchase Order #:	
Custom Molded CLIP file name:	*Note: CLIP file name= 1st Initial, 1st 3 of last name & Account #		

Mobility Base Information			
Brand:		Model:	
Frame Width:		Frame Depth:	
Cushion Information			
Additional Back Height in Inches:		Seat to Back Angle:	
Extend Existing Shape of Lateral Depth in Inches:	Left	Right	

What we need from you

- Custom Molded Shape
- Cushion Width
- Purchase Order #

On Chair Information
Is This Order Going "ON CHAIR" (NEW Invacare Chair)

When ordering cushions "On Chair" email this order form and Invacare Chair order form to: ElyriaCustomSeating@invacare.com or provide quote # below:

Invacare quote # : _____

Is This Order Going "ON CHAIR" (NEW Freedom Designs Chair)

If a Freedom Designs chair quote is to be converted please contact at 1-800-331-8551. Cushions will be drop shipped to Freedom Designs.

Freedom Designs quote / PO # : _____

Is This Order Going "ON CHAIR" (NEW Motion Concepts Chair)

If a Motion Concepts chair quote is to be converted please contact at 1-888-433-6818. Cushions will be drop shipped to Motion Concepts.

Motion Concepts quote / PO # : _____

Included Back Cushion Options
HCPCS CODE E2617
Model # CM02B 4U1 Retail \$2499
<ul style="list-style-type: none"> • Wood Mount w/ t-nuts (4 Row Pattern) • J & L Mounting Brackets (Includes 1" and 7/8" brackets) • Drawstring Cover (Reverse Dartex) • Soft Foam • Rail Cut (1.5 x 2) • Aluminum Reinforcements • Loose Fit (1/4" overall) • 1/2" Visco Foam Overlay • Blue Urethane Moisture Resistant Coating

INSTANT QUOTE
Online at www.invacare.com or call 1-800-451-3553

	Mounted 2" between the canes or flush mount
--	--

	3/8" Baltic Birch wood base 1/4 - 20 threaded t-nuts 1" spacing / 4 Rows
--	--

	1/8" Thick Aluminum supporting all contact surfaces Fabricated to match the cushion contours
--	---