



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ Tracer® SX5 Recliner Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

PRE-CONFIGURED STOCK WHEELCHAIRS

STOCK CHAIRS

Quantity	Model ID	Description	Price
<input type="checkbox"/>	TRSXSRC6P	TRSXSRC/WD66/28/R/U67 (9153637779).....	\$1,246.09
<input type="checkbox"/>	TRSXSRC8P	TRSXSRC/WD86/28/R/U67 (9153637778)	\$1,246.09

STOCK FRONT RIGGINGS

Quantity	Model ID	Description	Price
<input type="checkbox"/>	T93HAP	Hemi Sw. Footrests Aluminum Footplates (9153637801)	\$99.13
<input type="checkbox"/>	T93HCP	Hemi Sw. Footrests Composite Footplates (9153637786)	\$99.13
<input type="checkbox"/>	T94HAP	Hemi Sw. Elev Legrests Aluminum Footplates (9153637804)	\$191.86
<input type="checkbox"/>	T94HCP	Hemi Sw. Elev Legrests Composite Footplates (9153637788)	\$191.86
<input type="checkbox"/>	AHL4A	Hemi Smart Leg Art. Legrest Aluminum (9153619911)	\$341.09
<input type="checkbox"/>	T94HEP	Hemi Legrest Composite w/Black Non Rem Calfpad (9153637919)	\$117.25
<input type="checkbox"/>	T93HEP	Hemi Footrest without Loop Composite Footplate (9153637918)	\$69.28

FRAMES

CHAIR TYPE

TRSX5RC Tracer SX5 Recliner Builder Wheelchair \$1,097.61

FRAME WIDTH AND DEPTH

- WD46 14" Wide x 16" Deep \$96.39 A9900
- WD48 14" Wide x 18" Deep \$126.38 A9900
- WD66 16" Wide x 16" Deep No Charge
- WD68 16" Wide x 18" Deep \$126.38 A9900
- WD86 18" Wide x 16" Deep No Charge
- WD88 18" Wide x 18" Deep \$126.38 A9900
- WD06 20" Wide x 16" Deep No Charge E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).
- WD08 20" Wide x 18" Deep \$272.03
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).
- WD26 22" Wide x 16" Deep \$152.08
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).
- WD28 22" Wide x 18" Deep \$272.03 E2201
20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).

ARMS

ARM TYPE

All arm types below are removable, conventional and nonflip-back.

- 18 Full Length Fixed Height Conv. 10" \$134.95 A9900
- 19 Full Length Adj. Height Conv. 10" - 14" \$139.23 E0973
- 28 Desk Length Fixed Height Conv. 10" No Charge A9900
- 29 Desk Length Adj. Height Conv. 10" - 14" \$139.23 E0973

FRAME TYPE

- R Recliner Frame \$439.11 E1226
90 degree to 180 degree of recline. 24" back height with 10" headrest extension.

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ Tracer® SX5 Recliner Wheelchair

COLOR OPTIONS

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery No Charge
- U550 Midnight Blue Vinyl Upholstery Standard
- U67 Black Vinyl No Charge
- U86 Ruby Red Vinyl No Charge

CASTERS

CASTER OPTIONS

- AT903 8" x 1.25" Composite Urethane Standard
- 1228 8" x 1.75" Semi Pneumatic Casters Lt Grey \$59.98 A9900
- 4252 8" x 1.25" Pneumatic Casters w/FF Inserts \$191.71

WHEELS

WHEEL OPTIONS

- U1763BK 20" Composite Urethane Tires \$121.02 A9900
- U1764BK 22" Composite Urethane Tires \$121.02 A9900
- 2244BK 24" Composite Pneumatic Tires \$149.94 A9900
- 4244BK 24" Composite Pneumatic with FF Inserts \$240.98 E2213
- U2222C 24" Composite Urethane Corded Tires Standard
- U2244BK 24" Composite Urethane Treaded Tires \$128.52 A9900

WHEEL LOCKS

- 9690B Pull To Lock Wheel Locks (Std on Hemi) No Charge
- 9691B Push To Lock Wheel Locks (Std on Adult) No Charge

HANDRIM OPTIONS

- COM Composite Handrims Standard
Available with 24" rear wheels only.
- AHDRM Aluminum Handrims No Charge
Not available with U2222C rear wheels.

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Positioning Strap \$53.35 A9900 E0978
- 1321BK Hook and Loop Seat Positioning Strap \$42.84 E0978

ACCESSORIES

BACK ACCESSORIES

- L8601 Omit Back Upholstery Less \$0.00
- STC Stencil Back Upholstery \$31.06 A9270
3/4 Block Letters.

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$0.00

MISCELLANEOUS ACCESSORIES

- 1322SXR IV Rod and Bracket, Telescopic, Recliner \$117.81
- MUPH Matching Upholstery No Charge
Applies to arm pads and calf pads.
- 1329BK Toggle Wheel Lock Extension (8") \$41.77

FRONT RIGGING OPTIONS

FRONT RIGGINGS

- T93HE Hemi Footrests without heel loop, Composite Ftpl .. \$69.28 A9900
- T93HA Hemi Footrests Aluminum Footplates \$99.13
Not available with 14" wide.
- T93HC Hemi Footrests \$99.13
Not available with 14" wide.
- T94HE Hemi Elevating Legrest Composite Ftpl Blk Plastic Calf Pad \$117.25 A9900
- T94HA Hemi Elev. Lgrt Aluminum Ftpl \$191.86 K0195 E0990
*Not available with 14" wide.
All calf pads are black.*
- T94HC Hemi Elevating Legrests \$191.86 K0195 E0990
*Not available with 14" wide.
All calf pads are black.*
- AHL4A Hemi Smartleg Articulating Lgrt Alum Ftpl \$341.09 K0053
*Not available with 14" wide.
All calf pads are black.*
- AHL4C Hemi Smartleg Art. Lgrt Comp Ftpl \$341.09 K0053
*Not available with 14" wide.
All calf pads are black.*
- AHS4A Hemi Smartleg Art. Lgrt-14" W Chairs \$341.09 K0053
All calf pads are black.
- T93HS Hemi Footrests Aluminum Footplates 14" W \$95.93 A9900
- T94HS Hemi Legrests Aluminum Footplates 14" W \$186.53 A9900
All calf pads are black.

FOOTPLATE TYPE

- 1350 Extra Large Footplates \$84.21 K0041
Not available on 14" width.
- AT5543 Adjustable Angle Flip-Up Footplates \$71.42

FOOTREST/LEGREST ACCESSORIES

- 1820 Non-Removable Device \$12.79 A9900
- 1346 Heel Loops \$41.57 K0951
Standard on swingaway footrests T93HA and T93HC.

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.


Invacare® IVC™ Tracer® SX5 Recliner Wheelchair

Seat-to-Floor 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	Recliner Axle Mounting Bracket
18 1/3"	Middle	22"	Recliner Axle Mounting Bracket
19 1/2"	Bottom	24"	Recliner Axle Mounting Bracket

Standard on SX5 Recliner

- Rear anti-tippers
- Patient only operated push-to-lock wheel lock

Weight Limitation 250 lbs., 300 lbs. 20" & 22" width
 Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.



Comments

All prices subject to change without notice.
NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation
 USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2020 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.