



Yes, you can.®

## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Transport Wheelchairs

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

Order  Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

## INVACARE IVC TRACER TRANSPORT WHEELCHAIR

### SILVER VEIN TRANSPORT CHAIRS, INCLUDES MATCHING FRONT RIGGINGS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	TRAN17FR	Painted Transport 17" Seat Width, Perm Arm, Fold-Down Back, Footrest with Heel Loop .....	\$395.00
<input type="checkbox"/>	TRAN19FR	Painted Transport 19" Seat Width, Perm Arm, Fold-Down Back, Footrest with Heel Loop .....	\$395.00

### REPLACEMENT FRONT RIGGING – TRAN17FR / TRAN19FR INCLUDE RIGGINGS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	T93HA_PTO_37801	Hemi Swingaway Footrest Aluminum Footplate .....	\$93.00
<input type="checkbox"/>	T93HC_PTO_37786	Hemi Swingaway Footrest Composite Footplate .....	\$93.00
<input type="checkbox"/>	T94HA_PTO_37804	Hemi Swingaway Legrest Aluminum Footplate .....	\$180.00
<input type="checkbox"/>	T94HC_PTO_37788	Hemi Swingaway Legrest Composite Footplate.....	\$180.00
<input type="checkbox"/>	T93HE_PTO_37918	Hemi Footrests without Heel Loop, Composite Footplate .....	\$65.00
<input type="checkbox"/>	T94HE_PTO_37919	Hemi Elev. Legrests. Composite Footplate Blk Plastic Calf Pad .....	\$110.00

### STOCK OPTIONS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	1514	Overhead Anti - Theft Device .....	\$46.00
<input type="checkbox"/>	A1576	Wheel Lock Extensions 6" .....	\$39.00

Chair Dimensions		
Model Number	TRAN17FR	TRAN19FR
Seat Width	17"	19"
Seat Depth	16"	16"
Seat - to - Floor Height	19.5"	19.5"
Back Height	19"	19"
Back Height Folded	10"	10"
Arm Height	10"	10"

Overall Dimensions		
Model Number	TRAN17FR	TRAN19FR
Overall Width (Open)	19.5"	23.5"
Overall Width (Folded)	7.25"	7.75"
Overall Height with Canes	39"	39"
Overall Height with Back Folded	29"	29"
Length without Riggings	24"	24"
Length with Footrests	35"	35"
Weight without Riggings	24.5 lbs.	25 lbs.
Weight with Riggings (Footrests)	28.5 lbs.	29 lbs.

**Weight Limitations 250 Pounds**  
 Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® Transport Wheelchairs

## INVACARE IVC ALUMINUM TRANSPORT WHEELCHAIRS

### ALUMINUM LIGHTWEIGHT TRANSPORT WHEELCHAIR - STOCK CHAIR CONFIGURATIONS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	ALR19HBFR	19" x 16", Permanent Arm, 8" Front Casters, 12" Rear Wheels, Footrest, Red .....	\$430.00
<input type="checkbox"/>	ALB19HBFR	19" x 16", Permanent Arm, 8" Front Casters, 12" Rear Wheels, Footrest, Blue .....	\$430.00
<input type="checkbox"/>	LTTR19FR	19" x 16", Permanent Arm, 8" Wheels, Footrest, Red .....	\$370.00
<input type="checkbox"/>	LTTB19FR	19" x 16", Permanent Arm, 8" Wheels, Footrest, Blue .....	\$370.00
<input type="checkbox"/>	LTTB17FR	17" x 16", Permanent Arm, 8" Wheels, Footrest, Blue .....	\$370.00

Invacare Aluminum Lightweight Transport Wheelchair Specifications*			
Model Numbers	ALR19HBFR ALB19HBFR	LTTR19FR LTTB19FR	LTTB17FR
<b>Chair Dimensions</b>			
Seat Width	19"	19"	17"
Seat Depth	16"	16"	16"
Back Height	17.5"	17.5"	17.5"
Seat-to-Floor Height	19.25"	19.25"	19.25"
<b>Overall Dimensions</b>			
Overall Width	23.5"	23.5"	23.5"
Overall Depth w/o Riggings	24.25"	24.25"	24.25"
Overall Height	37.50"	36.50"	36.50"
Weight w/o Riggings	19 lbs.	19 lbs.	19 lbs.
Weight Limitation	250 lbs.	250 lbs.	250 lbs.

\*Note: All specifications and dimensions are approximate.

Invacare Aluminum Lightweight Transport Wheelchair Standard Features
<ul style="list-style-type: none"> <li>• Lightweight Aluminum Frames in Dark Blue or Red.</li> <li>• Permanent Arms with Padded Armrests.</li> <li>• Black Nylon Upholstery is Flame Retardant.</li> <li>• Seat Positioning Strap Included.</li> <li>• Footrests with Composite Footplates.</li> <li>• Thirteen Month Warranty.</li> </ul>

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

# Invacare® Transport Wheelchairs

## INVACARE HEAVY DUTY TRANSPORT WHEELCHAIR

### HEAVY DUTY STOCK WHEELCHAIRS

Quantity	Model	Description	List Price
<input type="checkbox"/>	TRHD22FR	22" x 18", Removable Arms, Footrests, Black .....	\$455.00

Specifications	
Model Number	TRHD22FR
Seat Width	22"
Seat Depth	18"
Back Height	18"
Seat - to - Floor Height	19.5"
Overall Width	28"
Overall Depth without Riggings	26.5"
Overall Height	39.5"
Weight limitation 400 pounds.	

**Invacare Heavy - Duty Transport Wheelchair**  
**Standard Features**

- Heavy Duty Frames incorporate a dual crossbrace and gusseted side frames.
- 22" width.
- Removable desk arms.
- Padded armrests provides additional patient comfort.
- Footrests with Aluminum Footplates.
- Black nylon upholstery is flame retardent.
- Three Year Warranty

Comments

**NOTE: All specifications and dimensions are approximate.**  
*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*



**Invacare Corporation**  
**USA**  
 One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2012 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.