

Invacare® Rollite™ Rollator



Yes, you can.™



Yes, you can.™

The new Invacare® Rollite™ Rollator is unlike any other in the industry, addressing the needs of providers, clinicians and consumers at an affordable price. With its innovative lightweight design, the Rollite offers a flip-up seat with a built-in handle, flexible backrest, ergonomic dual paddle folding mechanism and handbrakes. The Rollite provides the comfort, durability and quality that you've come to expect from an Invacare product.



Model no. 65100. Color: Midnight Blue.

Folds into a compact unit for storage. Built-in seat handle allows for easy transport.

Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height
- Patent pending

Specifications

Depth opened at base	24"
Depth folded	7"
Width inside hand grips	18.25"
Width inside base of legs	22"
Width overall	27"
Height adjustment min.-max.	32" - 36"
Seat-to-floor height min.-max.	19" - 23"
Seat width & depth	17" x 9.5"
Wheels	5" solid
Product weight	13 lb.
Weight capacity	250 lb.
Quantity	1 per carton
Limited warranty	lifetime frame
brakes	6 years
other components	1 year



Pictured left: tray. Pictured right: basket.

Basket and Tray Accessories

Available accessories for Rollite Rollator*

- Tray (model no. 65408)
 - For food and beverage placement
 - Available for all Invacare 4-wheel Rollators
- Basket (model no. 65110)
 - 10 lb. weight capacity
 - Position basket underneath seat to secure valuables
 - Available only for Rollite Rollator

*model no. 65100



Flip-up seat allows user to step inside the frame, when rollator is in a stationary position, providing greater reach access.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare is a registered trademark and Rollite and Yes, you can. are trademarks of Invacare Corporation. Specifications are subject to change without notification.