



Yes, you can.™

The Invacare® Compass™ SPT™ Limited adult tilt-in-space wheelchair provides consumers with access to a durable, lightweight tilt-in-space option when funding options are limited. Based on the revolutionary Compass SPT –Self Propelling Tilt wheelchair, the SPT Limited provides a proven wheelchair base along with an appropriate number of the most important options required to meet the consumer’s needs.



Compass SPT Limited shown in midnight blue with 70° pin and saddle style swingaway front rigging, flip-up composite footplates, 6"x2" semi-pneumatic casters, 12" composite wheels with pneumatic tires, Personal Back™ 10, Infinity® seat cushion, flush mounted seat pan, pull-to-lock wheel locks and conventional fixed height full-length arms.

The Compass SPT Limited Base Chair

- Rigid, rear wheel tilt-in-space frame with 45° of posterior tilt*
- Seat widths (15"-19") and seat depths (15"-19")
- Conventional fixed height arms
- Fixed back height (20" or 24")
- Back upholstery
- Three frame finishes (Wet Black, Red or Midnight Blue)
- 12" composite wheels
- Front rigging – Pin and Saddle Style or Center Pivot Options
- Flip-up composite footplates
- Longer pivot and slide tubes
- Permanent axles
- Pull-to-lock wheel locks
- Anti-tippers
- Adjustable angle back
- 17" front seat-to-floor height

*SPT Limited is only available in rear wheel configuration.

Other Available Options

- Larger seat widths and depths (20", 21" and 22")
- Conventional and cantilever adjustable height arms
- Flush mounted seat pan
- 6" x 2" semi-pneumatic casters
- Pneumatic tires
- Flat-free inserts
- Elevating or articulating legrests
- Adjustable angle flip-up footplates
- Heel and toe loops
- Legrest accessories
- Invacare® seating and positioning products
- Positioning straps
- Headrests



Quick Set™ washer (left) snaps onto the frame – preventing loss of cove washers (right) or your temper.

Newly designed adjustable angle footplate provides a positive lock for dorsi or plantar flexion.

Specifications

Frame style	Tilt
Seat widths	15"–22"
Seat depths	15"–22"
Overall chair width	22.5"–30.5" (Seat width plus 8.5")
Overall chair length	33"–39" (without anti-tippers and front rigging)
Overall chair height range	35.5"–39.5"
Back angle adjustment	90° to 120° in 10° increments
Tubing diameter	1"
Tilt range	0° to 45°
Weight capacity	250 lb.
Weight of chair	35 lb.* / 61 lb.**
Limited warranty	Lifetime on frames and crossbraces

*18"x18", AT902, no rear wheels, no footrests, no seat pan, no back upholstery, spreader bar, no anti-tippers, no wheel locks, no arms

**18"x18", AT902, 22COM/Urethane, PHW93/AT5543, seat pan, back upholstery, spreader bar, anti-tippers, wheel locks, no arms

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare and Infinity are registered trademarks and Yes, you can., the blue medallion design, Compass, Personal Back, Quick Set and SPT are trademarks of Invacare Corporation. Specifications are subject to change without notification. Accessories and seating products shown in the photographs are not necessarily included in the base price of wheelchair. Refer to current price list located at www.invacare.com

© 2004
Invacare Corporation
Form No. 04-033



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."

Invacare®
Compass™
SPT™ Limited
Wheelchair



Invacare® **Compass™** **SPT™** Limited
Tilt-in-Space Wheelchair



Yes, you can.™