



**Yes, you can.®**

**Price List and Order Form**

Effective July 1, 2003 Revised April 9, 2008

For ease of ordering, contact Customer Service Toll Free at:  
 800-333-6900 Fax: 800-678-4682  
 To order Literature: 800-828-6282

www.invacare.com

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

**FRAME TYPE**

- COMET-S Standard Configuration (E1237) ..... 2,034.00 <sup>1</sup>
- COMET-R Reverse Configuration (E1237) ..... 2,034.00

**SEAT WIDTH**

- W10 10" Wide ..... No Charge
- W11 11" Wide ..... No Charge
- W12 12" Wide ..... No Charge
- W13 13" Wide ..... No Charge
- W14 14" Wide ..... No Charge <sup>2</sup>
- W15 15" Wide ..... No Charge <sup>2</sup>
- W16 16" Wide ..... No Charge <sup>2</sup>

**ADJUSTABLE SEAT DEPTH**

- D0713 Short Frame (7"-13") ..... No Charge
- D0916 Long Frame (9"-16") ..... No Charge

**SEAT DEPTH**

- D07 Seat Depth 7" ..... No Charge <sup>3</sup>
- D08 Seat Depth 8" ..... No Charge <sup>3</sup>
- D09 Seat Depth 9" ..... No Charge
- D10 Seat Depth 10" ..... No Charge
- D11 Seat Depth 11" ..... No Charge
- D12 Seat Depth 12" ..... No Charge
- D13 Seat Depth 13" ..... No Charge
- D14 Seat Depth 14" ..... No Charge <sup>2</sup>
- D15 Seat Depth 15" ..... No Charge <sup>2</sup>
- D16 Seat Depth 16" ..... No Charge <sup>2</sup>

**SEAT STYLE**

- RDGSEAT Rigid Seat Pan ..... **Standard** <sup>4</sup>
- L9004 Omit Rigid Seat Pan ..... Less 25.00

**ARM TYPE**

- AT55 Cantilever, Adjustable Length..... No Charge <sup>5</sup>
- TAF T-Arms, Full Adjustable Height (E0973) ..... 225.00 <sup>6</sup>
- TAD T-Arms, Desk Adjustable Height (E0973) ..... 225.00 <sup>6</sup>
- L8401 Omit Arms ..... Less 20.00

**FOOTNOTES**

- <sup>1</sup> Camber angle is not available with reverse configuration.
- <sup>2</sup> Only available with long frame.
- <sup>3</sup> Not available with long frame.
- <sup>4</sup> Omit rigid seat pan when ordering a seating system.
- <sup>5</sup> Comes standard with desk length arm pads. For full length arm pads refer to arm accessories category. Length is measured from the front of the back cane to the front of the arms in one inch increments.
- <sup>6</sup> Not available with reverse configuration frame. Not available with 10" or less seat depth.

Arm Heights	
<b>Cantilever</b>	<b>7" – 12"</b>
Add 1" with upholstered arm pads Range adjusts in 1" increments	
<b>T-Arm</b>	
Low Range adjusts in ½" increments. High Range adjusts in 1" increments. Measurement taken from top of seat rail to top of arm rest.	

**T-ARM HEIGHT**

- L6-10 Low Range 6"-10" ..... No Charge
- H8-14 High Range 8"-14" ..... No Charge

**ARM ATTACHMENT (CANTILEVER ONLY)**

- BOLT Bolt-On ..... No Charge
- CLAMP Clamp-On..... 40.00

**TILT OPTION**

- MTILT 20° Mini Tilt ..... 445.00 <sup>7</sup>

**BACK HEIGHT**

- 1014 Sportster Back w/o Push Handles(10"-14") ...No Charge
- S1016 20" Fixed Height Back(10"-16").....No Charge
- A1014 10° Bend Back w/Push Handles(10"-14").....No Charge

**WHEEL LOCKS <sup>8</sup>**

- A1598 Hill Holder ..... 77.00 <sup>9 10</sup>
- A1576 Wheel Lock Extensions (6") (E0961)..... 39.00
- A1572 Reverse Configuration Wheel Locks ..... No Charge <sup>11</sup>

**CASTER OPTIONS**

- AT900 3" Roller Blades ..... No Charge
- AT901 5" Composite Urethane Casters ..... No Charge
- 5ALUM 5" Aluminum Urethane Casters ..... 92.00
- AT902 6" Composite Urethane 3 Spoke..... No Charge <sup>12</sup>
- 6PNEU 6" X 1 1/4" Composite Pneumatic ..... 64.00 <sup>13</sup>
- L8840 Omit Casters ..... Less 25.00

**FOOTNOTES**

- <sup>7</sup> Seat-to-floor height will increase 1 ½" with the mini-tilt and the S1016 is the only back height that is available with mini-tilt option. Comes standard with adjustable angle footboard.
- <sup>8</sup> Pull-to-lock wheel locks are standard.
- <sup>9</sup> Not available with adjustable angle swing-under footboard.
- <sup>10</sup> Not available with urethane tires or 10" – 14" seat depths.
- <sup>11</sup> Option for reverse configuration.
- <sup>12</sup> Only available with 9" – 16" long frame.
- <sup>13</sup> 3" fork is standard.

# Invacare® Comet™

## Accessories and Options

### WHEEL OPTIONS

<input type="checkbox"/> 16COM	16" Composite .....	No Charge	14 15
<input type="checkbox"/> 18COM	18" Composite .....	No Charge	15
<input type="checkbox"/> 20COM	20" Composite .....	No Charge	16
<input type="checkbox"/> 22COM	22" Composite .....	No Charge	16
<input type="checkbox"/> 24COM7	24" Composite 7 Spoke .....	No Charge	17
<input type="checkbox"/> 20SP	20" Spoke.....	No Charge	15 18
<input type="checkbox"/> 22SP	22" Spoke.....	No Charge	15 18
<input type="checkbox"/> 24SP	24" Spoke.....	No Charge	18
<input type="checkbox"/> L8701	Omit Rear Wheels .....	Less 50.00	

Front and Rear Seat-to-Floor Heights		
Caster & Wheels	Forward Configuration	Reverse Configuration
3" & 16"	10 ¾", 11 ¾"*	10 ¾", 11 ¾"
3" & 18"	11 ¾", 13 ¾"*	11 ¾", 13 ¾"
3" & 20"	12 ¾", 14 ½"*	12 ¾", 13 ¾", 14 ½"
3" & 22"	13 ¾", 15 ½"*	14 ½", 15 ½"*
3" & 24"	14 ¾", 16 ½"*	16 ½", 16 ½"
5" & 16"	11 ¾"	
5" & 18"	11 ¾", 13 ¾"	11 ¾", 13 ¾"
5" & 20"	13 ¾", 14 ½", 15 ½"*	14 ½", 15 ½"*
5" & 22"	14 ½"	14 ½"
5" & 24"	15 ½", 16 ½"*	16 ½"
6" & 20"	16 ½"	N/A
6" & 24"	17 ½"*	16 ½", 17 ½"*

\* Denotes 1" spacer      \*\*Denotes 2"spacer  
 Seat-to-Floor heights of 10 ¾", 11 ¾", 12 ¾" apply to 7" – 13" adjustable seat depths only. Add 1 ½" to seat-to-floor heights with mini-tilt option.

### SEAT-TO-FLOOR HEIGHT

<input type="checkbox"/> STF1034	Seat-To-Floor 10 3/4".....	No Charge	19
<input type="checkbox"/> STF1134	Seat-To-Floor 11 3/4".....	No Charge	19
<input type="checkbox"/> STF1234	Seat-To-Floor 12 3/4".....	No Charge	19
<input type="checkbox"/> STF1334	Seat-To-Floor 13 3/4".....	No Charge	
<input type="checkbox"/> FSTF145	Seat-To-Floor 14 1/2".....	No Charge	
<input type="checkbox"/> FSTF155	Seat-To-Floor 15 1/2".....	No Charge	
<input type="checkbox"/> FSTF165	Seat-To-Floor 16 1/2".....	No Charge	
<input type="checkbox"/> FSTF175	Seat-To-Floor 17 1/2".....	No Charge	

### MOUNTING HANDRIMS (SPOKE WHEELS ONLY)

<input type="checkbox"/> IN	Mount Inside Handrim .....	No Charge	20
<input type="checkbox"/> OUT	Mount Outside Handrim.....	<b>Standard</b>	

### TIRE OPTIONS

<input type="checkbox"/> PNEU	Pneumatic .....	No Charge
<input type="checkbox"/> PNEUFF	Pneumatic-Flat Free (E2213) .....	125.00
<input type="checkbox"/> URETH	Urethane (Treaded).....	51.00
<input type="checkbox"/> PRIMO	Primo .....	59.00
<input type="checkbox"/> CONTL	Continental Tubular.....	166.00
<input type="checkbox"/> KIKB	KIK® Black.....	60.00

### FOOTNOTES

- <sup>14</sup> Not available with 9" - 16" long frame.  
<sup>15</sup> Only available with PNEU, PNEUFF tires.  
<sup>16</sup> Only available with PNEU, PNEUFF or URETH tires.  
<sup>17</sup> Only available with PNEU, PNEUFF, URETH and KIKB tires.  
<sup>18</sup> Spoke wheels are standard with dome, quick-release button axle.  
<sup>19</sup> Only available with short frame.  
<sup>20</sup> Decreases overall width by ¾".

### HANDRIM OPTIONS

<input type="checkbox"/> ALUMPC	Aluminum Anodized Handrims .....	No Charge	21
<input type="checkbox"/> A1512	Plastic Coated Handrims .....	72.00	21
<input type="checkbox"/> A1512K	Black Plastic Coated Handrims .....	72.00	23
<input type="checkbox"/> P1512	Black Wrinkle Coated .....	No Charge	24
<input type="checkbox"/> AT1224-8	Plastic Coated w/8 Oblique Projections .....	245.00	25
<input type="checkbox"/> AT1225-8	Plastic Coated w/8 Vertical Projections .....	245.00	
<input type="checkbox"/> L8702	Omit Handrims .....	Less 10.00	

Axle Position Matrix	
Recommended Axle Position	
<b>Short Frame 7" – 13"</b>	
Caster & Wheel	Axle Position
3" & 16"	All Positions
5" & 16"	All Positions
3" & 20"	5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
5" & 20"	5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
<b>Long Frame 9" – 16"</b>	
Caster & Wheel	Axle Position
3" & 20"	All Positions
5" & 20"	2 <sup>nd</sup> , 3 <sup>rd</sup> , 4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
6" & 20"	3 <sup>rd</sup> , 4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
3" & 22"	All Positions
5" & 22"	2 <sup>nd</sup> , 3 <sup>rd</sup> , 4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
6" & 22"	4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
3" & 24"	3 <sup>rd</sup> , 4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
5" & 24"	4 <sup>th</sup> , 5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole
6" & 24"	5 <sup>th</sup> , 6 <sup>th</sup> , Furthest Hole

**Reverse Configuration uses first two axle positions only.**

### LARGE WHEEL AXLE POSITION

<input type="checkbox"/> A01	Axle In Furthest Forward Position .....	No Charge
<input type="checkbox"/> A02	Axle In Second Hole.....	No Charge
<input type="checkbox"/> A03	Axle In Third Hole.....	No Charge
<input type="checkbox"/> A04	Axle In Fourth Hole.....	No Charge
<input type="checkbox"/> A05	Axle In Fifth Hole.....	No Charge
<input type="checkbox"/> A06	Axle In Sixth Hole .....	No Charge
<input type="checkbox"/> A0F	Axle In Furthest Rearward Position .....	No Charge

### CAMBER ANGLE <sup>26</sup>

<input type="checkbox"/> C00	0° Camber .....	<b>Standard</b>
<input type="checkbox"/> C03	3° Camber .....	No Charge
<input type="checkbox"/> C05	5° Camber .....	No Charge

### FOOTNOTES

- <sup>21</sup> Not available with 16COM<sup>1</sup> or 18COM composite wheels.  
<sup>22</sup> Only available with spoke wheels.  
<sup>23</sup> Only available with 16COM, 18COM or 20SP.  
<sup>24</sup> Only available with 24COM7 and 24SP.  
<sup>25</sup> Only available with 22COM, 24COM7 and 22SP.  
<sup>26</sup> Not available with reverse configuration.

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information go to [www.invacare.com](http://www.invacare.com).

# Invacare® Comet™

## Accessories and Options

### CONTEMPORARY COLORS

- 24P Wet Black ..... No Charge
- 71P Silver Vein ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 105P Midnight Blue..... No Charge

### TOO HOT COLORS

- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge

### WAY COOL COLORS

- 62P Electric Blue ..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue ..... No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green..... No Charge
- 104P Electric Teal..... No Charge
- Touch Up Paint..... 20.00 <sup>27</sup>

### MONOGRAMMING & EMBLEMS

- TEXT ..... 15.00 <sup>28</sup>
- EM00 Cometeer..... No Charge
- EM01 Butterfly Emblem ..... 26.00
- EM02 Deer Emblem ..... 26.00
- EM03 Rabbit Emblem..... 26.00
- EM06 Corvette Emblem..... 26.00
- EM07 F-14 Jet Fighter Emblem ..... 26.00
- EM08 Football Emblem ..... 26.00
- EM10 Basketball & Hoop Emblem ..... 26.00

### THREAD COLOR FOR EMBLEM

- TC02 Pink Thread ..... No Charge
- TC06 Teal Thread..... No Charge
- TC08 Purple Thread..... No Charge
- TC19 Green Thread ..... No Charge
- TC29 White Thread..... No Charge
- TC30 Yellow Thread..... No Charge
- TC33 Red Thread ..... No Charge

### SEAT POSITIONING STRAPS

- A1515C Auto Rainbow Seat Positioning Strap..... 50.00
- 1321BK Hook & Loop Seat Positioning Strap..... 40.00
- A1321 Hook & Loop Seat Positioning Strap..... 40.00

### CHEST POSITIONING STRAPS

- 7515C Auto Style Rainbow Chest Positioning Strap ..... 50.00
- 7321BK Hook & Loop Chest Positioning Strap ..... 40.00
- 7321C Hook & Loop Rainbow Chest Positioning Strap... 40.00

### BACK ACCESSORIES

- BP03 Comet Backpack ..... No Charge
- BP02 Back Pack ..... 50.00
- ADJTENS Adjustable Tension Back Upholstery..... 87.00
- 1483 Stroller Handles ..... 85.00
- 1485 Clamp-On Stroller Handles ..... 75.00
- L8601 Omit Back Upholstery ..... Less 10.00

### SEAT ACCESSORIES

- SP-01 Seat Pouch..... 35.00

### ARM ACCESSORIES

- A1571 Fabric Clothing Protectors(K0108)..... 46.00
- 1671-SM Rigid Side Guards, Small (K0108)..... 122.00 <sup>29</sup>
- 1671-SMH Quad Rigid Side Guards, Small (K0108)..... 122.00 <sup>29</sup>
- 1671-LG Rigid Side Guards, Large (K0108)..... 122.00 <sup>29</sup>
- 1671-LGH Quad Rigid Side Guards, Large (K0108)..... 122.00 <sup>29</sup>
- A1492 Cantilever Full Length Arm Pads ..... 35.00
- L8405 Omit Arm Pads ..... Less 10.00

### MISCELLANEOUS ACCESSORIES

- A1573 Quad Release Axles (K1080) ..... 89.00
- A1307 Spoke Protectors(For 20", 22", 24") (K0065)..... 90.00
- A1209 Tire Pump ..... 32.00
- 3421 3/4" To 1" Clamp Assembly ..... 21.00
- A1300 Utility Tray ..... 47.00
- PUSHHANDL Push Handles Bolt On..... 64.00
- A1569 Crutch & Cane Carrier..... 46.00

Knee-to-Heel Range							
90° SWUNDER		90° FTBD		90° INDFTPT		AT5046	
<b>Long Frame</b>							
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
4"	15 1/2"	0"	15 1/2"	0"	16"	12"	16"
<b>Short Frame</b>							
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
4"	15 1/2"	0"	15 1/2"	0"	16"	12"	16"
Range adjusts in 1" increments.							
Measurement taken from top of seat rail to back of footrest.							

### FOOTPLATE TYPE

- 90SWUNDER Adjustable Angle Swing Under Ft. Brd ..... No Charge
- 90FTBD Adjustable Angle Footboard..... No Charge
- 90INDFTPT Adjustable Angle Indiv. Footplates..... No Charge
- AT5046 Rigid Tubular Footrests ..... No Charge

### FOOTREST/LEGREST ACCESSORIES

- A1337 Calf Strap (K0038) ..... 34.00
- A1600BK Heel Loops With Ankle Straps (E0951)..... 51.00
- 1346 Heel Loops ..... No Charge

### FOOTNOTES

<sup>27</sup> Touch-up paint is not available in all colors.  
<sup>28</sup> No thread color choice available.

### FOOTNOTES

<sup>29</sup> Not available with T-Arms. Not available in reverse configuration. Only available with 11" or greater seat depth in standard configuration.

**Invacare® Comet™**  
Accessories and Options

**Specifications**

• Frame Style	Rigid
• Seat Widths	10" – 16"
• Adjustable Height Back	10" – 16" (See Back Height Category)
• Adjustable Angle Back	-4°, 0°, 4°, 8°, 12°, 16°
• Camber Angles	0°, 3°, 5°
• Tubing O.D.	1" Aluminum
• Weight Limit 7" – 13" Adj. Seat	100 lbs
• Weight Limit 9" – 16" Adj. Seat	150 lbs
• Weight of Chair	22.5 lbs
• Warranty	Lifetime on frame and crossbraces

**Minimum Overall Width**

Seat width plus 9"  
Measurement taken from outside of handrim to outside of handrim with 0° camber.

**Weight Limitation 100/150 Pounds**

Invacare recommends all wheel chairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**INVACARE REHAB SEATING AND POSITIONING PRODUCTS  
PRICE LIST FORM # 07-027**



**Invacare Corporation**

**USA**  
One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

[www.invacare.com](http://www.invacare.com)

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

KIK is a registered trademark of KIK Technology, Inc.  
Specifications and prices are subject to change without notification.  
©2007 Invacare Corporation

Form No. 95-80 Rev 0408