

Invacare® Patriot™ - Lightweight Folding Wheelchair



Price List and Order Form

Effective July 1, 2003 Revised April 2, 2007

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Yes, you can.®

Company Name _____ Account #. _____ Date _____

Phone # _____ P.O.# _____

Ship To Name & Address _____

FRAME SELECTION

1 FRAME TYPE

K0004 ■ PATRIOT Folding Frame w/Swingaway Front End..... 991.00

FRAME WIDTH

2 W14 14" Wide..... 100.00
 W16 16" Wide No Charge
 W17 17" Wide No Charge
 W18 18" Wide No Charge
 E2201 W20 20" Wide..... 300.00

SEAT DEPTH

D16 16" Seat Depth No Charge
 D18 18" Seat Depth No Charge

3 SEAT-TO-FLOOR HEIGHT

K0056 SHE 15 1/2" Front Seat-To-Floor Height..... 100.00
 HEMI 17 1/2" Front Seat-To-Floor Height..... No Charge
 ADULT 19 1/2" Front Seat-To-Floor Height..... No Charge

ARM & BACK SELECTION

ARM TYPE

L8401 Omit Arms.....Less 75.00
 4 78 SpaceSaver Fixed Ht Full Lgth Arms (pair).....No Charge
 4 88 SpaceSaver Fixed Ht Desk Lgth Arms (pair).....No Charge
 18 Conventional Fixed Ht Full Lgth Arms (pair).....No Charge
 28 Conventional Fixed Ht Desk Lgth Arms (pair).....No Charge
 E0973 4 79 SpaceSaver Adjust Ht Full Lgth Arms (pair).....195.00
 E0973 4 89 SpaceSaver Adjust Ht Desk Lgth Arms (pair).....195.00
 E0973 19 Conventional Adjust Ht Full Lgth Arms (pair).....195.00
 E0973 29 Conventional Adj Ht Desk Lgth Arms (pair).....195.00
 E0973 4 FULLADJLWSS SpaceSaver Low Adj Ht Full Lgth Arms (pair).....195.00
 E0973 4 DESKADJLWSS SpaceSaver Low Adj Desk, Lgth Arms (pair).....195.00

BACK TYPE

5 ADJ Adjustable Back w/10° Bend.....No Charge
 6 AT20 10° Bend w/Push Handles 20"..... 75.00
 5 S1519 Straight Back w/Push Handles (15"-19").....No Charge
 7 ANGADJ Angle Adjustable Back w/10° Bend.....165.00
 6 ANGAT20 Angle Adjustable 10° Bend w/Push Handles 20".....165.00
 7 ANGS1519 Angle Adj. Straight Back w/Push Handles (15"-19").....165.00

CASTER & WHEEL SELECTION

CASTER OPTIONS

AT902 6" Composite Urethane 3-Spoke Casters (pair) ..No Charge
 6PNEU 6" X 1 1/4" Composite Pneumatic Casters (pair).....
 AT903 8" X 1 1/4" Composite Urethane Casters (pair)**Standard**
 K0074 8 AT1292 8" X 1 1/4" Pneumatic Casters (pair).....60.00
 K0064 7 1492 Flat Free Inserts for 8" X 1 1/4" Pneumatic Casters (pair)64.00
 L8840 Omit CastersLess 25.00
 L8790 Omit Caster And Fork AssembliesLess 50.00

WHEEL OPTIONS

1763BK 20" Composite Pneumatic Tires (pair).....171.00
 4763BK 20" Composite Pneumatic Flat Free Inserts (pair).....232.00
 U1763BK 20" Composite Urethane Tires (pair)113.00
 1764BK 22" Composite Pneumatic Tires (pair).....117.00
 4764BK 22" Composite Pneumatic Flat Free Inserts (pair)232.00
 U1764BK 22" Composite Urethane Tires (pair)113.00
 U2222 24" Composite Urethane Tires Low Profile (pair)..... 50.00
 2244BK 24" Composite Pneumatic Tires (pair).....140.00
 4244BK 24" Composite Pneumatic w/FF Inserts (pair).....225.00
 9 U2244BK 24" Composite Urethane Tires (pair)**Standard**
 L8701 Omit Rear WheelsLess 50.00

HANDRIM OPTIONS

ALUMPC Aluminum Anodized Handrims (pair).....**Standard**
 A1512 Plastic Coated Handrims (pair)..... 72.00 K0059
 COM Composite Handrim (pair)No Charge
 11 P1512 Black Wrinkle Coated (pair)No Charge
 11 AT1224-8 Plastic Coated w/8 Oblique Projections (pair)240.00 E0967
 11 AT1225-8 Plastic Coated w/8 Vertical Projections (pair).....240.00 E0967
 11 AT1225-12 Plastic Coated w/12 Vertical Projections (pair).....240.00 E0967
 11 L8702 Omit Handrims.....Less 15.00

CLD ONE ARM DRIVE

LDR96RI CLD Lever Drive RH 6" Install.....797.00
 LDR96LI CLD Lever Drive LH 6" Install.....797.00
 LDR98RI CLD Lever Drive RH 8" Install.....797.00
 LDR98LI CLD Lever Drive LH 8" Install.....797.00

WHEEL LOCK

11,12 9695 Hill Holder Wheel Locks - Bolt Through (pair).....170.00 E0974
 12 9690 Pull-To-Lock Wheel Lock (Std. On Hemi) (pair) ..No Charge
 12 9691 Push-To-Lock Wheel Lock (Std. On Adult) (pair) ..No Charge
 A1567 Pull-To-Lock Clamp On Wheel Locks (pair)..... 76.00
 1721 Hi Push-To-Lock Clamp On Wheel Locks (pair).....71.00
 A1576 Wheel Lock Extensions (6") (pair) 65.00 E0961

FOOTNOTES

1. Seat angle is set at 3 degrees. See service manual for details.
2. Not available with 18" Seat Depth.
3. Refer to Seat-to-Floor Matrix on bottom of page 2 to determine proper caster and rear wheel size selections for desired seat-to-floor height.
4. Not available with adjustable angle back type.
5. Back height is set at 16".
6. AT20 back height is 20".
7. Rear anti-tippers come standard with this selection. Back angle adjusts -5 degree to 20 degree in 5 degree increments.
8. P70 or P70T front rigging are not available with 8" pneumatic or pneumatic flat-free casters if a 16" deep frame and a hemi seat-to-floor is chosen.
9. High profile treaded urethane tires.
10. Only available w/24" rear wheels.
11. Not available with urethane tires.
12. Bolt through wheel locks. Available in push or pull-to-lock.

Note: All specifications and dimensions are approximate.

Invacare® Patriot™ - Lightweight Folding Wheelchair

Accessories and Options

FRAME & UPHOLSTERY SELECTION

FRONT RIGGING SELECTION

FRAME FINISH

- 24P Wet BlackNo Charge
- 71P Silver VeinNo Charge
- 60P Silver Metallic.....No Charge
- 115P Black Prism.....No Charge
- 105P Midnight Blue.....No Charge
- 30P Sunny YellowNo Charge
- 61P Electric RedNo Charge
- 119P Bubblegum Pink.....No Charge
- 120P Grape MadnessNo Charge
- 62P Electric Blue.....No Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic BlueNo Charge
- 125P Emerald Green.....No Charge
- 127P Grasshopper Green.....No Charge
- 104P Electric Teal.....No Charge

FRONT RIGGING

- P70 70° Swingaway Footrests (pair).....**Standard**
- P60 60° Swingaway Footrests (pair).....No Charge
- P70T 70° Tapered Footrest (pair)No Charge
- P30E Elevating Legrests (pair)255.00 E0990
- ATL4A Center Pivot Articulating Legrests (pair)320.00 K0053
- FT70 70° Footrest - Composite Footplates w/ChairNo Charge
- AT5044 Elevating Legrests.....255.00

UPHOLSTERY COLOR

- U240 Black Nylon Upholstery.....No Charge
- U245 Black Sail Cloth Adjustable Tension Seat**Standard**

FOOTPLATE TYPE

- A1562 Aluminum Footplates (pair)No Charge
- AT5543 Adjustable Angle Flip-Up Footplates (pair)165.00 K0040
- 1651 Flip-Up Composite Footplates (pair)**Standard**
- A1564 Extra Large Aluminum Footplates (pair)59.00 K0041

TOUCH-UP COLOR

- Touch-Up Paint Colors20.00

FOOTREST/LEGREST ACCESSORIES

- A1337 Calf Strap (each).....24.00 K0038
- IMPACTG Impact Guards (pair).....51.00
- ALPT Longer Pivot And Slide Tube (pair)84.00
- A9342 "H" Strap (each).....51.00 K0039
- A1600BK Heel Loops With Ankle Straps (pair)51.00 E0951
- HLOOP Heel Loops Without Ankle Straps (pair)51.00 E0951
- L8515 Omit Calf Pads (pair)Less 10.00

POSITIONING STRAPS

SEAT POSITIONING

- E0978 A1515 Push Button Seat Positioning Strap (each).....50.00
- E0978 A1321 Hook & Loop Seat Positioning Strap (each)50.00

CHEST POSITIONING

- E0960 9515 Push Button Chest Positioning Strap (each)100.00
- E0960 7321BK Hook & Loop Chest Positioning Strap (each)100.00

ACCESSORIES

ARM ACCESSORIES

- 1474W Full/Desk Uph. Waterfall Arm Pad (pair)25.00
- NRAS Non-Removable Arm with Screws (pair)No Charge
- L8405 Omit Arm Pads OnlyLess 10.00
- A1571 Fabric Clothing Protectors (pair)46.00
- 1300 Utility Tray Full Length.....157.00

BACK ACCESSORIES

- L8601 Omit Back Upholstery.....Less 10.00
- ADJTENS Adjustable Tension Back Upholstery125.00
- L8651 Omit Seat UpholsteryLess 10.00
- BP02 Back Pack50.00

MISCELLANEOUS ACCESSORIES

- E0971 P1580 Anti-Tippers For 8" Casters (pair)150.00
- A1577 Quick Release Axles (pair).....60.00
- K0102 A1569 Crutch & Cane Carrier46.00
- E0959 9357L Amputee Attachments (pair)95.00
- K0105 2218 IV Rod, Telescopic148.00
- E0971 P1580A Anti-Tippers For 6" Casters - 17 1/2" STF (pair)150.00
- E0971 P1580SH Anti-Tippers For 6" Casters - 15 1/2" STF (pair)150.00

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

FOOTNOTES

1. Not available with adjustable angle back type.
2. For 0° and 3° seat angles, standard w/adjustable angle back type. See service manual for details.
3. P70 or P70T front riggings are not available with 8" pneumatic or pneumatic flat-free casters if a 16" deep frame and a hemi seat-to-floor is chosen.
4. Not available on 14" wide or w/extra large footplates.
5. Must use pull-to-lock wheel locks, due to interference.
6. Not available on 14" wide.
7. Only available with aluminum footplates.
8. Standard on P70T footrests.
9. Not available in all colors.

Note: All specifications and dimensions are approximate.

Invacare® Patriot™ - Lightweight Folding Wheelchair Accessories and Options

Basic Chair Includes:	
♦ Folding Frame w/Swingaway Front End	
♦ Seat Width	(W16 - W18)
♦ Seat Depth	(D16 - D18)
♦ Seat-To-Floor Height	(HEMI, ADULT)
♦ Arm Type	(78, 88, 18, 28 or NRAS)
♦ Back Type Selection	(ADJ or S1519)
♦ Frame Finish	
♦ Caster Options	(AT902 or AT903)
♦ Wheel Options	(U2244BK w/URETH)
♦ Handrim Options	(ALUMPC, COM or P1512)
♦ Seat Upholstery	
♦ Front Rigging	(P70, P60 or P70T)
♦ Footplate Type	(A1562 or 1651)
♦ Slip on Back Upholstery	
♦ Wheel Locks	(9690 or 9691)

Knee-to-Heel Range							
70° Tapered		70°		60°		Elevating	
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
14 3/4"	18 3/4"	12 1/2"	17 1/2"	13 3/4"	17 1/2"	14 1/4"	20 1/4"
Range adjusts in 1" increments. Measurement taken from top of seat rail to back of footrest.							

Seat-to-Floor Matrix 6" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Bottom	24"	Middle (upper)
17 1/2"	Bottom	22"	Middle (lower)
17 1/2"	Bottom	20"	Bottom
15 1/2"	Top	20"	Middle (upper)
For a complete list of Seat-to-Floor heights and angles that are configurable by the dealer, refer to service manual.			

Seat-to-Floor Matrix 8" Casters			
Seat-to-Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
19 1/2"	Bottom	24"	Bottom
17 1/2"	Top	24"	Middle (upper)
17 1/2"	Top	22"	Middle (lower)
17 1/2"	Top	20"	Bottom
For a complete list of Seat-to-Floor heights and angles that are configurable by the dealer, refer to service manual.			

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications and prices are subject to change without notification.
©2007 Invacare Corporation

Form No. 03-010 Rev 04/07
P2003_010F_Patriot