



**Yes, you can.®**

**Price List and Order Form**

Effective September 15, 2003 **Revised July 17, 2008**

For ease of ordering, contact Customer Service Toll Free at:  
 800-333-6900 Fax: 800-678-4682  
 To order Literature: 800-828-6282

[www.invacare.com](http://www.invacare.com)

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

**Special Note**

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

**FRAME SELECTION**

**FRAME TYPE**

- SPRXT Spree XT Folding Tilt (E1234) ..... 2,895.00 <sup>1</sup>

**FRAME WIDTH**

- W10 10" Wide (Growable to 12")..... No Charge
- W11 11" Wide (Growable to 12")..... No Charge
- W12 12" Wide (Growable to 14")..... No Charge
- W13 13" Wide (Growable to 14")..... No Charge
- W14 14" Wide (Growable to 16")..... No Charge
- W15 15" Wide (Growable to 16")..... No Charge
- W16 16" Wide (Growable to 18")..... No Charge
- W17 17" Wide (Growable to 18")..... No Charge
- W18 18" Wide ..... No Charge

**SEAT DEPTH RANGE**

- FT1214 Short Seat Rail (12" - 14" Deep) ..... No Charge
- FT1218 Long Seat Rail (12" - 18" Deep) ..... No Charge

**SEAT RAIL EXTENSION**

- SRE2 2" Seat Rail Extension Installed ..... 25.00

**SEAT DEPTH SETTING <sup>2</sup>**

- D12 Seat Depth 12"..... No Charge
- D13 Seat Depth 13" ..... No Charge
- D14 Seat Depth 14"..... No Charge
- D15 Seat Depth 15"..... No Charge <sup>3</sup>
- D16 Seat Depth 16"..... No Charge <sup>3</sup>
- D17 Seat Depth 17"..... No Charge
- D18 Seat Depth 18"..... No Charge
- D19 Seat Depth 19"..... No Charge <sup>4</sup>
- D20 Seat Depth 20" (E2201)..... 300.00 <sup>4</sup>

**SEAT-TO-FLOOR HEIGHT**

(MUST WRITE IN SEAT TO FLOOR TO COMPLETE ORDER)

- Seat-To-Floor ..... No Charge

Seat-To-Floor Heights							
*Measured From Floor To Top of Seat Rail							
Rear Wheel	3"	4"	5"	6"	6"	8"	Tilt Range
	AT900	4SPLT 4SOFT	AT901 5SOFT	AT902 6PNEU	1259	AT903	
12"	16.00 16.50 17.00	15.50 16.00 16.50 17.50	16.00 16.50 17.00 17.50	17.00 17.50	17.50	N/A	45°
20"	N/A	N/A	17.50	17.50 18.00 18.50	17.50 19.00	18.50	45°
22"	16.50 17.00	16.50 17.00	16.50 17.00 17.50	17.00 17.50 18.00	17.50	N/A	30°
24"	N/A	N/A	17.50	17.50 18.00 18.50	17.50 19.00	18.50	30°

All Dimensions Are Approximate

**FRAME FINISHES**

**CONTEMPORARY COLORS**

- 24P Wet Black ..... No Charge
- 71P Silver Vein ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 105P Midnight Blue ..... No Charge

**TOO HOT COLORS**

- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge

**WAY COOL COLORS**

- 62P Electric Blue ..... No Charge
- 121P Lolly Pop Blue ..... No Charge
- 122P Cosmic Blue ..... No Charge
- 125P Emerald Green ..... No Charge
- 127P Grasshopper Green ..... No Charge
- 104P Electric Teal ..... No Charge
- Touch Up Paint ..... 20.00

**FOOTNOTES**

- <sup>1</sup> A seating system must be purchased for use with this wheelchair.
- <sup>2</sup> No seat pan option is offered with this wheelchair.
- <sup>3</sup> Setting requires selection of FT1214 and SRE2.
- <sup>4</sup> Setting requires selection of FT1218 and SRE2.

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information go to [www.invacare.com](http://www.invacare.com).

# Invacare® Spree XT Accessories and Options

## TRANSPORT READY OPTION

### TRANSPORT TIE DOWN <sup>5</sup>

- TRRO Transport Ready Option ..... No Charge
- TRBKTS Transport Brackets & Pelvic Belt ..... No Charge
- TROMIT Less Transport Ready Option..... Less 150.00

#### \*Transport Tie Down

TRRO includes four factory-installed transport brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lb to 209 lb.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

**\* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.**

## ARM & BACK SELECTIONS

### ARM TYPE (PAIR)

- AT55D Cantilever Locking Desk  
Length. Bolt On Arms ..... 195.00
- AT55F Cantilever Locking Full  
Length Bolt On Arms ..... 195.00
- AT99N Cantilever Non-Locking Bolt On Arms ..... No Charge
- TAD T-Arms, Desk Adjustable Height (E0973) ..... 225.00
- TAF T-Arms, Full, Adjustable Height (E0973) ..... 225.00
- L8401 Omit Arms ..... Less 20.00
- L8405 Omit Arms Pads ..... Less 10.00

### T-ARM HEIGHT

- L8-12 Low Range..... No Charge
- H10-16 High Range..... No Charge

### ARM ATTACHMENT – CANTILEVER ARMS

- BOLT Bolt-On ..... No Charge
- CLAMP Clamp-On ..... 40.00 <sup>6</sup>

Arm Heights	
AT99N	6" – 9"
Range adjusts in 1" increments	
AT55D	9" – 12"
AT55F	9" – 12"
Range adjusts in 1" increments	
TAD	Low Range: 8" – 12"
TAF	High Range: 10" – 16"
All measurements from top of seat rail to top of armrest pad.	

### FIXED BACK HEIGHT <sup>7</sup>

- BH20 20" Height..... No Charge
- BH24 24" Height..... No Charge

### FOOTNOTES

<sup>5</sup> Transport Tie Down options not retrofittable to existing models, not field serviceable. **TRRO requires KSS seating system to be ordered.**

<sup>6</sup> Not available with AT99N.

<sup>7</sup> Back does not come with upholstery

## CASTERS & WHEEL LOCKS

### WHEEL LOCKS

- A1567 Pull To Lock Clamp On Wheel Locks (pair)..... No Charge
- A1598 Hill Holder Clamp On Wheel Locks (pair)..... 77.00 <sup>8</sup>
- 9691 Push to Lock Clamp On Wheel Locks (pair)..... No Charge <sup>8</sup>
- ATHUB Cane Mount Hublock (attendant) (E0961)..... 50.00 <sup>9</sup>
- USEHUB Seat Rail Mount Hublock ..... 50.00 <sup>9</sup>

### CASTER OPTIONS (PAIR)

- AT900 3" Roller Casters ..... 38.00
- 4SOFT 4" Semi Pneumatic Casters ..... 60.00
- 4SPLT 4" Comp. Ureth. Light Up Casters ..... No Charge
- 5SOFT 5" Semi Pneumatic Casters ..... 60.00
- AT901 5" Composite Urethane Casters ..... 64.00
- AT902 6" Composite Urethane 3 Spoke Casters..... No Charge
- 6PNEU 6" Composite Pneu. 6 Spoke Casters ..... 64.00
- 1259 6" X 2" Semi Pneumatic Casters ..... 60.00
- AT903 8" X 1.25" Composite Urethane Casters ..... 78.00
- AT1292FF Flat Free Inserts for 8" X 1.25" Composite  
Pneumatic Casters ..... 60.00
- 1292 8" X 1.25" Composite Pneu. Casters ..... 60.00

## REAR WHEELS

### WHEEL OPTIONS (PAIR)

- 12COM 12" Composite Wheels ..... No Charge
- 20COM 20" Composite Wheels ..... No Charge
- 22COM 22" Composite Wheels ..... No Charge
- 24COM7 24" Composite 7 Spoke Wheels ..... No Charge
- 20SP 20" Spoke Wheels..... No Charge
- 22SP 22" Spoke Wheels..... No Charge
- 24SP 24" Spoke Wheels..... No Charge

### TIRE OPTIONS (PAIR)

- PNEU Pneumatic Tires ..... 60.00
- PNEUFF Flat Free Inserts (E2213)..... 65.00 <sup>10</sup>
- URETH Urethane (Treaded) Tires..... No Charge <sup>11</sup>
- PRIMO Primo Tires ..... 60.00 <sup>11</sup>
- KIKB KIK® Black Tires ..... 60.00 <sup>12</sup>

### AXLE TYPE (PAIR)

- A1577 Quick Release Axles..... No Charge
- PAXLE Permanent Axles..... No Charge

### HANDRIM OPTIONS (PAIR)

- ALUMPC Aluminum Anodized Handrims ..... No Charge <sup>13</sup>
- A1512K Black Plastic Coated Handrims ..... 72.00 <sup>13</sup>
- AT1224-8 Plastic Coated w/8 Oblique Projection..... 245.00 <sup>12,13</sup>
- AT1225-8 Plastic Coated w/8 Vertical Projection..... 245.00 <sup>12,13</sup>
- AT1225-12 Plastic Coated w/12 Vertical Projection..... 245.00 <sup>13,14</sup>
- L8702 Omit Handrims ..... Less 15.00

### MOUNTING HANDRIMS (SPOKE ONLY)

- IN Mounting Inside Handrim ..... No Charge <sup>15</sup>
- OUT Mounting Outside Handrim..... No Charge

### FOOTNOTES

<sup>8</sup> Not available with 22" and 24".

<sup>9</sup> Only available with 12COM & 20COM rear wheels.

<sup>10</sup> Not available on 12COM, 20SP, 22SP and 24SP.

<sup>11</sup> Only available on 24SP.

<sup>12</sup> Only available with 24COM7 and 24SP.

<sup>13</sup> Not available with 12COM.

<sup>14</sup> Only available on 24COM7.

<sup>15</sup> Decreases overall width by .75".

# Invacare® Spree XT Accessories and Options

## ACCESSORIES

### BACK ACCESSORIES

- BP02 Back Pack ..... 50.00
- 1483 Non-Adjustable Rigidizing Handle Bar ..... **Standard**
- AASH Adj. Angle Stroller Handles (Outside Mount) ..... 235.00 <sup>16</sup>
- AASHIN Adj. Angle Stroller Handles (Inside Mount) ..... 235.00 <sup>16</sup>
  
- RIGB Flip Down Rigidizer Bar ..... 85.00
- SHOMIT Omit Non Adjust Rigidizer Handle Bar ..... Less 85.00

### SEAT ACCESSORIES

- SP-01 Seat Pouch ..... 35.00

### ARM ACCESSORIES (PAIR) <sup>17</sup>

- 1671-SM Rigid Side Guard, Small (K0108) ..... 122.00
- 1671-LG Rigid Side Guard, Large (K0108) ..... 122.00
- 1671-SMH Quad Rigid Side Guard Small (K0108) ..... 122.00
- 1671-LGH Quad Rigid Side Guard Large (K0108) ..... 122.00

### MISCELLANEOUS ACCESSORIES

- 3421 ¾" To 1" Clamp Assembly (each) ..... 21.00
- A1307 Spoke Protectors (For 20", 22", 24", 26") (K0065) ..... 90.00 <sup>18</sup>
- TLC Tilt Lock Out Collars (pair) ..... 35.00 <sup>19</sup>
- A1576 Wheel Lock Extensions (6") (pair) (E0961) ..... 39.00

#### POSITIONING STRAPS (EACH)

- 1515 Push Btn. Buckle Seat Pos. Strap (E0978) ..... 50.00
- 1321BK Hook & Loop Seat Pos. Strap (E0978) ..... 40.00
- 9515 Push Btn. Chest Pos. Strap (E0960) ..... 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960) ..... 100.00

## FRONT RIGGINGS

### FRONT RIGGING (PAIR)

- SWING1319 60° Swingaway Footrests ..... No Charge <sup>20</sup>
- SWING914 60° Swingaway Footrests MFX ..... No Charge <sup>20</sup>
- 70NTAPER 70° Swingaway Footrests ..... No Charge <sup>20</sup>
- 70MFX 70° Swingaway Footrests MFX ..... No Charge <sup>20 21</sup>
- 90SWING 90° Swingaway Footrests ..... No Charge <sup>22</sup>
- AT5044 Elevating Legrests (E0990) ..... 255.00 <sup>20 23</sup>
- 70STAPER 70° Tapered Footrests ..... No Charge
- 70LIFTNT 70° Lift off Non-Taper Footrests ..... 95.00
- 90LIFTMFX 90° Lift off MFX Footrests ..... 95.00

### FOOTNOTES

- <sup>16</sup> AASH release buttons are mounted on the outside and AASHIN release buttons are mounted on the inside.
- <sup>17</sup> Not available with TAD and TAF.
- <sup>18</sup> Not available with projection handrims.
- <sup>19</sup> Additional pair of tilt lock out collars, one pair standard with chair.
- <sup>20</sup> Not available with ALPT2or ALPT3.
- <sup>21</sup> 16" – 18" widths only. Standard with AT5543 Footplates.
- <sup>22</sup> Standard with AT5543 Footplates.
- <sup>23</sup> 14" – 18" widths only.

### FOOTPLATE TYPE (PAIR)

- AT5543 Adjustable Angle Flip-Up Footplates (K0040) ..... 165.00 <sup>24</sup>
- 1651 Flip-Up Composite Footplates ..... No Charge
- A1562 Aluminum Footplates ..... 59.00 <sup>25</sup>

### HEEL LOOP/TOE LOOP (PAIR)

- A1600BK Heel Loops w/ Ankle Straps (E0951) ..... 51.00
- HLOOP Heel Loops w/o Ankle Straps (E0951) ..... 51.00
- 1345 Toe Loops ..... 39.00

### FOOTREST/LEGREST ACCESSORIES (PAIR)

- 1337 Calf Strap (K0038) ..... 24.00
- IMPACTG Impact Guards ..... 51.00 <sup>26</sup>
- ALPT Longer Pivot And Slide Tubes ..... No Charge <sup>27</sup>
- ALPT2 Foot Support Extension 2" for 90° ..... 25.00
- ALPT3 Foot Support Extension 3" for 90° ..... 31.00
- A9342 "H" Strap (K0039) ..... 51.00

Knee-To-Heel Range		
Front Rigging Type	Min.	Max
60°	13.50	17.50
60° Max	9.00	12.50
70°	13.50	16.50
70° Tapered	12.75	16.25
70° MFX	5.75	11.75
90°	4.50	7.50
Elevating	14.00	20.00

Range adjust in 1: increments  
Measurement taken from top of seat rail to back of footrests.

## SEATING

### SEATING CUSHION SELECTION <sup>28</sup>

- EC Absolute Cushion (E2601) ..... 104.00 <sup>29</sup>
- ITSG InTouch™ Stabilite™ Gentle Contour Cushion (E2607) ..... 340.00
- ITSM InTouch Stabilite Max Contour Cushion (E2605) 340.00
- ITFG InTouch Flovair™ Gentle Contour Cushion (E2607) ..... 360.00
- ITFM InTouch Flovair Max Contour Cushion (E2607). 360.00

### SEATING BACK SELECTION <sup>28</sup>

- PBR Personal Back 10™ Reg. Height (E2615) ..... 565.00
- PBT Personal Back 10 Tall Height (E2615) ..... 565.00
- PPR Personal Back 10 Plus Reg. Height (E2615) ..... 645.00
- PPT Personal Back 10 Plus Tall Height (E2615) ..... 645.00
- ITPR InTouch Propel™ Back (E2611) ..... 438.00

**FOR ADDITIONAL SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**INVACARE REHAB SEATING AND POSITIONING PRODUCTS  
PRICE LIST FORM # 07-027.**

### FOOTNOTES

- <sup>24</sup> Dorsi/Plantar angle adjustable.
- <sup>25</sup> Not available on 90SWING, 70MFX and 90LIFTMFX.
- <sup>26</sup> Not available with 90SWING, 70MFX, SWING914 and AT5044.
- <sup>27</sup> 2" longer slide tube used on composite and adjustable angle footplates only. Only available with 60°, 70°, and 70° Tapered.
- <sup>28</sup> For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.
- <sup>29</sup> Available in even sizes only. See form 07-027 for additional information.

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).**

# Invacare® Spree XT

## Basic Chair Includes:

- Spree XT Folding Tilt
- TRRO or TRBKTS Transport Tie Down Option
- Seat Width Selection
- Seat Depth Setting (D12 - D19)
- Arm Type (AT99N)
- Bolt On Arm Attachment
- Fixed Back Height Selection
- Frame Finish
- Caster Options (2SPLT or AT902)
- Seat-to-Floor Height Selection
- Wheel Options
- URETH - Urethane (Treaded) Tires
- Axle Type (A1577 or PAXLE)
- ALUMPC - Aluminum Anodized Handrims
- Front Rigging Selection - SWING1319, SWING914, 70NTAPER, 70STAPER, 70MFX or 90SWING
- 1651- Flip-Up Composite Footplates (pair)
- ALPT - Longer Pivot and Side Tubes
- 1483 - Non-Adjustable Rigidizing Handle Bar
- Wheel Lock Selection (9691 or A1567)
- Anti-Tippers (pair)
- Tilt Trigger Lock-Out Levers
- Tilt Lock Out Collars (1 pair)

## Weight Limitation 150 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

## Specifications

- |   |  |
|---|--|
| • Frame Style                                 | Folding Tilt   |
| • Seat Widths                                 | 10" - 18"  |
| • Seat Depths                                 | 12" - 20" (with seat rail extension)                             |
| • Back Angle Adj.                             | -4° to 16°   |
| • Tubing Size                                 | 1"   |
| • Tilt Range                                  | 0° - 45°   |
| • Overall Chair Width                         | Seat width plus 7.50" (12, 20)<br>Seat width plus 8.75" (22, 24) |
| • Warranty Lifetime on frames and crossbraces |  |

## Specials Department

For your customization requests, please contact Customer Service for details and prices.

## Growth Kit Available



## Invacare Corporation

USA

One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

[www.invacare.com](http://www.invacare.com)

All rights reserved. Trademarks are identified by the symbols <sup>TM</sup> and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Otto Bock is a trademark of Otto Bock Orthopaedische Industries GMBH & Co.

Freedom Designs is a registered trademark of Freedom Designs, Inc.

Specifications and prices are subject to change without notification.

©2007 Invacare Corporation

Form No. 07-024 Rev 0708