

Invacare® Compass XE™ Adult Positioning Folding Chair



Yes, you can.®

Price List and Order Form

Effective January 15, 2008 Revised April 9, 2008

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

FRAME TYPE

CXE Folding, Depth Adj Extended Frame (K0005) 2,236.00

CXE CROSSBRACE TYPE

CONV Modular Crossbrace No Charge

LSR Less Seat Rails No Charge

TRANSPORT TIE DOWN ¹

TRRO Transport Ready Option No Charge

TRBKTS Transport Brackets No Charge

TROMIT Less Transport Ready Option Less 150.00

Transport Ready Option

*TRRO Must be ordered with conventional crossbraces, seat cushion and PERSONAL BACK

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

FRAME DEPTH

MEDM Medium Side Frames (12" - 16") No Charge

LONG Long Side Frames (14" - 18") No Charge

XLONG Extra Long Side Frames (16" - 20") No Charge

FOOTNOTES

¹ Transport Tie Down options not retrofittable to existing models, not field serviceable. TRRO and TRBKTS not available with LSR (Less Seat Rail).

SEAT SELECTIONS

SEAT WIDTH

W12 12" Wide No Charge

W13 13" Wide No Charge

W14 14" Wide No Charge

W15 15" Wide No Charge

W16 16" Wide No Charge

W17 17" Wide No Charge

W18 18" Wide No Charge

W19 19" Wide No Charge

W20 20" Wide (E2201) 300.00

Available Seat Widths/Depths for Compass XE

	D12	D13	D14	D15	D16	D17	D18	D19	D20
W12	X	X	X	X	S	S	S	S	S
W13	X	X	X	X	X	S	S	S	S
W14	X	X	X	X	X	X	X	S	S
W15	S	X	X	X	X	X	X	S	S
W16	S	X	X	X	X	X	X	S	S
W17	S	X	X	X	X	X	X	S	S
W18	S	S	X	X	X	X	X	X	X
W19	S	S	S	S	X	X	X	X	X
W20	S	S	S	S	X	X	X	X	X

S = Call Customer Service for availability of Special size.

SEAT DEPTH

D12 12" Seat Depth No Charge

D13 13" Seat Depth No Charge

D14 14" Seat Depth No Charge

D15 15" Seat Depth No Charge

D16 16" Seat Depth No Charge

D17 17" Seat Depth No Charge

D18 18" Seat Depth No Charge

D19 19" Seat Depth No Charge

D20 20" Seat Depth (E2203) 300.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare Compass XE Adult Positioning Folding Chair Accessories and Options

Arm Heights	
AT99N Range adjusts in 1" increments	5" - 10"
AT55D Range adjusts in 1" increments	5" - 13"
AT55F Range adjusts in 1" increments	5" - 13"
DESKADJCV, FULLADJCV Range adjusts in 1" increments	10" - 14"
TAD Low Range Range adjusts in 1/2" increments	6" - 10"
TAF High Range Range adjusts in 1" increments	8" - 14"
All measurements from top of seat rail to top of armrest. On AT55D and AT55F, add 1" for upholstered armrests. For less seat rail chairs, add 1-1/2" to above measurements (top of upper frame rail to top of armrest).	

ARM SELECTIONS

ARM TYPE (PAIR)

- AT99N Cantilever Non Locking Bolt On Arms No Charge ²
- FXFULL Conventional Fixed Ht.Full Lgth. Arms No Charge ³
- FXDESK Conventional Fixed Ht.Desk Lgth. Arms No Charge ³
- DESKADJCV Conventional Adjustable Ht Desk Lgth Arms(E0973).... 195.00 ³
- FULLADJCV Conventional Adj. Ht.Full Lgth Arms (E0973).. 195.00 ³
- AT55D Cantilever Locking Desk Length
Bolt On Arms (E0973) 195.00 ²
- AT55F Cantilever Locking Full Length
Bolt On Arms (E0973) 195.00 ²
- TAD T-Arm Adjustable Height, Desk Length (E0973)225.00 ⁴
- TAF T-Arm Adjustable Height, Full Length (E0973) 225.00 ⁴

ARM SELECTIONS (CONTINUED)

T-ARM HEIGHT

- L6-10 Low Range 6" - 10" (T-Arm Height) No Charge
- H8-14 High Range 8" - 14" (T-Arm Height) No Charge

BACK SELECTIONS

BACK TYPE

- A1418 10° Bend Back w/Push Handles (14"-18")..... No Charge
- 1014 Sportster Back w/o Push Handles (10"-14") ..No Charge ⁵
- S1216 Sportster Back w/o Push Handles (12"-16") ..No Charge ⁵
- A1216 Straight Back w/ Push Handles (12"-16") 30.00
- S1418 Straight Back w/ Push Handles (14"-18") 30.00
- FBH18 Fixed Back Height 18" No Charge
- FBH20 Fixed Back Height 20" No Charge
- RECLINER Recliner Frame No Charge

BACK UPHOLSTERY

- SLIPBACK Slip On Back Upholstery No Charge
- ADJTENS Adjustable Tension Back Upholstery (K0108)... 125.00
- L8601 Omit Back Upholstery Less 10.00

RECLINER BACK ⁶

- REC18 Reclining Back 18" High (E1226) 570.00
- REC20 Reclining Back 20" High (E1226) 570.00
- REC24 Reclining Back 24" High (E1226) 570.00
- RECEXT Reclining Headrest Pillow and Extension (E0966).....75.00 ⁷

FOOTNOTES

- ² Only available with FBH18 and FBH20 fixed back height.
- ³ Only available with CONV, Conventional Crossbraces and Seat Depths greater than 15".
- ⁴ Bolt-On Arm Attachment.
- ⁵ For push handles see Back Accessories section.
- ⁶ 1" diameter anodized straight canes with push handles comes standard with back upholstery. Must select when RECLINER (Recliner Frame) is selected.

CASTER & WHEEL SELECTIONS ⁸

Front and Rear Seat-To-Floor Heights *Measured From Floor To Top Of Seat Rail					
Fork	3" AT900	5" 5ALUM AT901	6" AT902 6SKYWAY	6" 1259	8" 1228 AT903
A1103 3" Fork	15.25" 15.75" 16.50" 17.00"	16.25" 16.75" 17.50" 18.00"	N/A	N/A	N/A
A1105 5" Fork	16.25" 16.75" 17.25" 17.50" 18.00" 18.50"	16.75" 17.25" 17.75" 18.00" 18.25" 18.50" 19.00" 19.50" 19.50"	17.25" 17.75" 18.25" 18.50" 18.75" 19.00" 19.50" 20.00"	N/A	N/A
A1108 8" Fork	N/A	N/A	N/A	18.00" 18.50" 19.00" 19.25" 19.50" 19.75" 20.25" 20.75"	19.00" 19.50" 20.00" 20.25" 20.50" 20.75" 21.25" 21.75"
All Dimensions Are Approximate					

FRONT SEAT-TO-FLOOR HEIGHT ⁹

- FSTF1525 15.25" Front Seat-To-Floor Height (K0056)..... 100.00
- FSTF1575 15.75" Front Seat-To-Floor Height (K0056)..... 100.00
- FSTF1625 16.25" Front Seat-To-Floor Height (K0056)..... 100.00
- FSTF165 16.5" Front Seat-To-Floor Height (K0056)..... 100.00
- FSTF1675 16.75" Front Seat-To-Floor Height (K0056)..... 100.00
- FSTF17 17" Front Seat-To-Floor Height No Charge
- FSTF1725 17.25" Front Seat-To-Floor Height No Charge
- FSTF175 17.5" Front Seat-To-Floor Height No Charge
- FSTF1775 17.75" Front Seat-To-Floor Height No Charge
- FSTF18 18" Front Seat-To-Floor Height No Charge
- FSTF1825 18.25" Front Seat-To-Floor Height No Charge
- FSTF185 18.5" Front Seat-To-Floor Height No Charge
- FSTF1875 18.75" Front Seat-To-Floor Height No Charge
- FSTF19 19" Front Seat-To-Floor Height No Charge
- FSTF1925 19.25" Front Seat-To-Floor Height No Charge
- FSTF1950 19.50" Front Seat-To-Floor Height No Charge
- FSTF1975 19.75" Front Seat-To-Floor Height No Charge
- FSTF20 20" Front Seat-To-Floor Height No Charge
- FSTF2025 20.25" Front Seat-To-Floor Height No Charge
- FSTF2050 20.50" Front Seat-To-Floor Height No Charge
- FSTF2075 20.75" Front Seat-To-Floor Height No Charge
- FSTF2125 21.25" Front Seat-To-Floor Height (K0056)..... 100.00

FOOTNOTES

- ⁷ RECEXT must be chosen with reclining back selection.
- ⁸ For folding frame without seat rails, subtract 1-1/2" from the seat-to-floor height to obtain distance from floor to top of side frame.
- ⁹ Rear seat-to-floor height is set approximately 1" lower than the front seat-to-floor height.

Invacare Compass XE Adult Positioning Folding Chair

Accessories and Options

CASTER OPTIONS (PAIR)

<input type="checkbox"/> AT900	3" Roller Blade Casters	38.00
<input type="checkbox"/> AT901	5" Composite Urethane Casters.....	64.00
<input type="checkbox"/> 5ALUM	5" Aluminum Urethane Casters.....	92.00
<input type="checkbox"/> AT902	6" Composite Urethane 3 Spoke Casters	No Charge
<input type="checkbox"/> 6SKYWAY	6" Composite Urethane 6 Spoke Casters	No Charge
<input type="checkbox"/> 1259	6" X 2" Semi Pneumatic Casters (K0075)	60.00
<input type="checkbox"/> 1228	8" X 1.75" Semi Pneumatic Casters Lt Grey	60.00
<input type="checkbox"/> AT903	8" X 1.25" Composite Urethane Casters	No Charge

FORK TYPE

<input type="checkbox"/> A1103	3" Fork (Single Hole Fork).....	No Charge
<input type="checkbox"/> A1105	5" Fork (3 Hole Fork).....	No Charge
<input type="checkbox"/> A1108	8" Fork (2 Hole Fork).....	No Charge

WHEEL OPTIONS (PAIR) ¹⁰

<input type="checkbox"/> 20COM	20" Composite Wheels	No Charge
<input type="checkbox"/> 22COM	22" Composite Wheels	No Charge
<input type="checkbox"/> 24COM7	24" Composite 7-Spoke Wheels.....	No Charge
<input type="checkbox"/> 20SP	20" Spoke Wheels.....	No Charge
<input type="checkbox"/> 22SP	22" Spoke Wheels.....	No Charge
<input type="checkbox"/> 24SP	24" Spoke Wheels.....	No Charge
<input type="checkbox"/> L8701	Omit Rear Wheels	Less 50.00

TIRE OPTIONS (PAIR)

<input type="checkbox"/> PNEU	Pneumatic Tires	60.00
<input type="checkbox"/> PNEUFF	Flat Free Inserts (E2213).....	65.00
<input type="checkbox"/> URETH	Urethane (Treaded) Tires	No Charge
<input type="checkbox"/> PRIMO	Primo Tires	60.00
<input type="checkbox"/> KIKB	KIK [®] Black Tires	60.00

NEW OPTION

<input type="checkbox"/> CRF123	Schwalbe 1 in Blk Flat Resistant Non Marking.....	15
	Tire.....	90.00

MOUNTING HANDRIMS (SPOKE WHEELS ONLY)

<input type="checkbox"/> IN	Mounting Inside Handrims (pair).....	No Charge
<input type="checkbox"/> OUT	Mounting Outside Handrims (pair).....	No Charge

HANDRIM OPTIONS (PAIR)

<input type="checkbox"/> ALUMPC	Aluminum Anodized Handrims	No Charge
<input type="checkbox"/> A1512	Black Plastic Coated Handrims.....	72.00
<input type="checkbox"/> AT1225-8	Plastic Coated w/8 Vertical Projections	240.00
<input type="checkbox"/> AT1224-8	Plastic Coated w/8 Oblique Projections	240.00
<input type="checkbox"/> AT1225-12	Plastic Coated w/12 Vertical Projections.....	240.00
<input type="checkbox"/> NFITSD	Natural Fit Handrims Standard Grip	295.00
<input type="checkbox"/> NFITSP	Natural Fit Handrims Super Grip.....	340.00
<input type="checkbox"/> L8702	Omit Handrims	Less 15.00

AXLE & CAMBER SELECTIONS

AXLE PLATE OPTIONS (PAIR)

<input type="checkbox"/> STDPT	Standard Axle Plates	No Charge
<input type="checkbox"/> MPOSAP	Multi-Position Offset Axle Plates	No Charge
<input type="checkbox"/> EXTPT	Extended/Amputee Axle Plates (E0959)	72.00

FOOTNOTES

- ¹⁰ Spoke wheels are standard with dome button axle.
- ¹¹ Not available on medium frame set at 14" deep.
- ¹² Not available with 20SP, 22SP, or 24SP.
- ¹³ Only available on 24SP.
- ¹⁴ Only available on 24COM7 and 24SP.
- ¹⁵ The Schwalbe tire is a black high-pressure pneumatic tire that is non-marking, and includes a Smart Guard[®] puncture resistant layer. Schwalbe tire is only available with: 24SP,26SP,24HPS,26HPS,SPIN24,SPIN25,SPIN26.
- ¹⁶ Decrease overall width by ¾".
- ¹⁷ Only available on 22COM, 24COM7 and 24SP.

AXLES (PAIR)

<input type="checkbox"/> A1577	Quick Release Axles.....	No Charge
<input type="checkbox"/> A1573	Quad Release Axles (K0108)	89.00
<input type="checkbox"/> PAXLE	Permanent Axles	Standard

CAMBER ANGLE

<input type="checkbox"/> C00	0° Camber.....	No Charge
<input type="checkbox"/> C03	3° Camber.....	No Charge
<input type="checkbox"/> C06	6° Camber.....	No Charge

FRAME FINISH SELECTIONS

CONTEMPORARY COLORS

<input type="checkbox"/> 24P	Wet Black	No Charge
<input type="checkbox"/> 71P	Silver Vein	No Charge
<input type="checkbox"/> 60P	Silver Metallic	No Charge
<input type="checkbox"/> 115P	Black Prism.....	No Charge
<input type="checkbox"/> 105P	Midnight Blue.....	No Charge

TOO HOT COLORS

<input type="checkbox"/> 30P	Sunny Yellow	No Charge
<input type="checkbox"/> 61P	Electric Red.....	No Charge
<input type="checkbox"/> 119P	Bubblegum Pink	No Charge
<input type="checkbox"/> 120P	Grape Madness	No Charge

WAY COOL COLORS

<input type="checkbox"/> 62P	Electric Blue	No Charge
<input type="checkbox"/> 121P	Lolly Pop Blue.....	No Charge
<input type="checkbox"/> 122P	Cosmic Blue.....	No Charge
<input type="checkbox"/> 125P	Emerald Green	No Charge
<input type="checkbox"/> 127P	Grasshopper Green	No Charge
<input type="checkbox"/> 104P	Electric Teal.....	No Charge
<input type="checkbox"/>	Touch Up Paint.....	20.00

POSITIONING STRAPS

SEAT POSITIONING STRAPS (EACH)

<input type="checkbox"/> A1515	Push Button Buckle Seat Pos. Strap (E0978).....	50.00
<input type="checkbox"/> A1321	Hook & Loop Seat Positioning Strap (E0978).....	50.00

CHEST POSITIONING STRAPS (EACH)

<input type="checkbox"/> 9515	Push Button Chest Positioning Strap (E0960)	100.00
<input type="checkbox"/> 7321	Hook & Loop Chest Positioning Strap (E0960).....	100.00

ACCESSORIES

BACK ACCESSORIES

<input type="checkbox"/> BP02	Back Pack	50.00
<input type="checkbox"/> 1483	Non-Adjustable Rigidizing Handle Bar	85.00
<input type="checkbox"/> AASH	Adjustable Angle Rigidizing Handle Bar.....	235.00
<input type="checkbox"/> RIGB	Flip Down Rigidizer Bar	85.00

SEAT ACCESSORIES

<input type="checkbox"/> SP-01	Seat Pouch.....	35.00
<input type="checkbox"/> L8651	Omit Seat Upholstery	Less 10.00

ARM ACCESSORIES (PAIR) ²⁰

<input type="checkbox"/> A1571	Fabric Clothing Protector (K0108).....	46.00
<input type="checkbox"/> 1671-SM	Rigid Side Guard, Small (K0108)	122.00
<input type="checkbox"/> 1671-LG	Rigid Side Guard, Large (K0108)	122.00
<input type="checkbox"/> 1671-SMH	Quad Rigid Side Guard, Small (K0108)	122.00
<input type="checkbox"/> 1671-LGH	Quad Rigid Side Guard, Large (K0108)	122.00

FOOTNOTES

- ¹⁸ Not available with Reclining Back Option.
- ¹⁹ Not available on 19" or 20" widths. Not available with Reclining Back Option.
- ²⁰ Only available with AT99N, AT55D, or AT55F.

Invacare Compass XE Adult Positioning Folding Chair

Accessories and Options

WHEEL LOCKS (PAIR)

- A1598 Hill Holder Clamp On Wheel Locks (E0974) 170.00 ²¹
- A1567 Pull to Lock Clamp On Wheel Locks No Charge
- HUMW Hideaway Undermount Wheel Locks 65.00 ²²
- 9691 Push to Lock Clamp On Wheel Locks No Charge

MISCELLANEOUS ACCESSORIES

- 3421 ¾ " To 1" Clamp Assembly 21.00
- A1541 Caster Swivel Locks (pair) (K0073) 80.00 ²³
- A1569 Crutch & Cane Holder (E2207) 46.00
- A1307 Spoke Protectors (for 20", 22", 24", 26") (K0065)..... 90.00
- A1209 Tire Pump 32.00
- A1580 Rear Anti-Tippers (pair) (E0971)..... 150.00
- A1576 Wheel Lock Extensions (6") (pair) (E0961) 65.00

Knee-to-Heel Range							
For Folding Frames with seat rails.							
70° Tapered		70°		70° MFX		Elevating	
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
16.25"	20.25"	14"	19"	8.25"	14"	15.75"	21.75"
90°		60°		60° MFX			
MIN	MAX	MIN	MAX	MIN	MAX		
5.75"	10.75"	15.25"	19"	11.5"	15.5"		

Range adjusts in 1" increments.

FRONT RIGGING SELECTIONS

FRONT RIGGING (PAIR) ²⁴

- SWING1319 60° Swingaway Footrests No Charge
- SWING914 60° Swingaway Footrests MFX No Charge ²⁵
- 70MFX 70° Swingaway Footrests MFX No Charge
- 70NTAPER 70° Swingaway Footrests No Charge
- 70STAPER 70° Tapered Footrests..... No Charge ²⁶
- 90SWING 90° Swingaway Footrests No Charge
- 70LIFTNT 70° Lift Off Non-Taper Footrests 95.00
- 70LIFTMFX 70° Lift Off MFX Footrests 95.00
- 90LIFTMFX 90° Lift Off MFX Footrests 95.00
- AT5044 Elevating Legrests (E0990)..... 255.00
- ATL4A Center Pivot Articulating Legrests (K0053)..... 320.00

FOOTPLATE OPTION (PAIR)

- 1651 Flip-Up Composite Footplates No Charge ²⁷
- AT5543 Adjustable Angle Flip-Up Footplates(K0040)..... 165.00 ²⁸
- A1564 Extra Large Aluminum Footplates59.00 ²⁹
- TUBULAR Individual Tubular Footplates.....59.00

FOOTNOTES

- ²¹ Not available with URETH or KIKB tires.
- ²² Not available with ATL4A Center Pivot Articulating Legrests.
- ²³ Only available on tubular footplates AT900, 5ALUM, AT901, AT902, 6SKYWAY.
- ²⁴ All front riggings are center pivot style. Footplates are attached directly to rigging tube when MFX style is selected.
- ²⁵ Standard with AT5543 adjustable angle footplate.
- ²⁶ 16" – 20" wide frames only. Standard with AT5543 adjustable angle footplate.
- ²⁷ Not available with 70MFX, 90SWING Swingaway footrests or 12" wide.
- ²⁸ AT5543 footplates are dorsi-plantar adjustable.
- ²⁹ Not available with 70STAPER or 12" or 15" wide.

MAX KNEE TO HEEL LENGTH

- KTH14 14" Knee to Heel.....No Charge
- KTH1516 15"-16" Knee to Heel.....No Charge
- KTH1618 16"-18" Knee to Heel.....No Charge

FOOTREST/LEGREST ACCESSORIES (PAIR)

- A1337 Calf Strap (K0038) 24.00 ³⁰
- IMPACTG Impact Guards 51.00 ³¹
- ALPT Longer Pivot And Slide Tubes 84.00
- A1600BK Heel Loops with Ankle Straps (E0951) 51.00
- HLOOP Heel Loops w/o Ankle Straps (E0951)..... 51.00
- 9342 "H" Strap (Each) (K0039)..... 51.00

SEATING CUSHION SELECTION ³²

- EC Absolute Cushion (E2601) 101.00 ³³
- ITSG InTouch™ Stabilite™ Gentle Contour Cushion (E2607) 330.00
- ITSM InTouch Stabilite Max Contour Cushion (E2605) 330.00
- ITFG InTouch Flovair™ Gentle Contour Cushion (E2607)..... 350.00
- ITFM InTouch Flovair Max Contour Cushion (E2607). 350.00

SEATING BACK SELECTION ³²

- PBR Personal Back 10™ Reg. Height (E2615)..... 550.00
- PBT Personal Back 10 Tall Height (E2615)..... 550.00
- PPR Personal Back 10 Plus Reg. Height (E2615)..... 625.00
- PPT Personal Back 10 Plus Tall Height (E2615)..... 625.00
- ITPR InTouch Propel™ Back (E2611)..... 425.00

**FOR ADDITIONAL SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027**

FOOTNOTES

- ³⁰ Not available with elevating legrests.
- ³¹ 2" longer slide tube used on composite, adjustable angle or tubular footplates only, not available with 70MFX, 90SWING, 90LIFTMFX, SWING914, and ATL4A.
- ³² **For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.**
- ³³ Available in even sizes only. See form 07-027 for additional information.

Invacare Compass XE Adult Positioning Folding Chair

Basic Chair Includes:

- Folding, Depth Adjustable Extended Frame
- Crossbrace Type Selection
- Transport Tie Down Options
- Seat Width (W12 - W20)
- Frame Depth Selection
- Seat Depth (D12 - D20)
- Arm Type (AT99N, FXDESK, FXFULL)
- Back Type Selection (A1418, 1014, S1216, FBH18, FBH20)
- Slip-On Back Upholstery
- Frame Finish
- Caster Options (AT902, AT903 or 6SKYWAY)
- Seat to Floor Height Selection (FSTF1725 – FSTF2075)
- Wheel Options
- Tire Options (URETH)
- Handrims Options (ALUMPC)
- Axle Type (A1577 or PAXLE)
- Axle Brackets (STDPT or MPOSAP)
- Camber Angle Selection
- Front Rigging (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER or 90SWING).
- Flip-Up Composite Footplates (1651)
- Seat Upholstery
- Wheel Locks (9691 or A1567)

Specifications

Frame Style	Folding
Seat Widths	12" – 20"
Seat Depths	12" – 20"
Overall Chair Width (closed)	Seat width plus 11.5"*
Overall Chair Width (open)	Seat width plus 8.5"*
Adjustable Height Back	14" – 18"
Back Angle Adjustment	80° to 110° in 5° increments
Camber Angles	0°, 3°, 6°
Tubing Outside Diameter	1"
Weight of Chair	19 lbs**/36lbs***
Warranty	Lifetime on frames and crossbraces
*0° Angles	
**18" x 16", AT902, no rear wheels, no footrests, no anti-tippers, no wheel locks, no arms.	
***18" x 16", AT902, 24SP, PNEU, 70NTAPER, seat upholstery, back upholstery, wheel locks, no arms.	

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.



Invacare Corporation

USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

All rights reserved. Trademarks are identified by the symbols TM and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

KIK is a registered trademark of KIK Technology, Inc.
 Specifications and prices are subject to change without notification.
 ©2008 Invacare Corporation