

Invacare® Compass XE™ Adult Positioning Folding Chair



Yes, you can.®

Price List and Order Form

Effective January 15, 2008 Revised June 10, 2009

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

FRAME TYPE

CXE Folding, Depth Adj Extended Frame (K0005)2,236.00

CXE CROSSBRACE TYPE

CONV Modular Crossbrace..... No Charge

LSR Less Seat Rails..... No Charge

TRANSPORT TIE DOWN ¹

TRRO Transport Ready Option..... No Charge

TRBKTS Transport Brackets No Charge

TROMIT Less Transport Ready OptionLess 150.00

Transport Ready Option

*TRRO Must be ordered with conventional crossbraces, seat cushion and PERSONAL BACK

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

SEAT SELECTIONS

SEAT WIDTH

W12 12" Wide No Charge

W13 13" Wide No Charge

W14 14" Wide No Charge

W15 15" Wide No Charge

W16 16" Wide No Charge

W17 17" Wide No Charge

W18 18" Wide No Charge

W19 19" Wide No Charge

W20 20" Wide (E2201) 300.00

Available Seat Widths/Depths for Compass XE

	D12	D13	D14	D15	D16	D17	D18	D19	D20
W12	X	X	X	X	S	S	S	S	S
W13	X	X	X	X	X	S	S	S	S
W14	X	X	X	X	X	X	X	S	S
W15	S	X	X	X	X	X	X	S	S
W16	S	X	X	X	X	X	X	S	S
W17	S	X	X	X	X	X	X	S	S
W18	S	S	X	X	X	X	X	X	X
W19	S	S	S	S	X	X	X	X	X
W20	S	S	S	S	X	X	X	X	X

S = Call Customer Service for availability of Special size.

SEAT DEPTH

D12 12" Seat Depth No Charge

D13 13" Seat Depth No Charge

D14 14" Seat Depth No Charge

D15 15" Seat Depth No Charge

D16 16" Seat Depth No Charge

D17 17" Seat Depth No Charge

D18 18" Seat Depth No Charge

D19 19" Seat Depth No Charge

D20 20" Seat Depth (E2203)..... 300.00

FRAME DEPTH

MEDM Medium Side Frames (12" - 16") No Charge

LONG Long Side Frames (14" - 18") No Charge

XLONG Extra Long Side Frames (16" - 20")..... No Charge

FOOTNOTES

¹ Transport Tie Down options not retrofittable to existing models, not field serviceable. TRRO and TRBKTS not available with LSR (Less Seat Rail).

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare Compass XE Adult Positioning Folding Chair

Accessories and Options

Arm Heights	
AT99N	5" - 10"
Range adjusts in 1" increments	
AT55D	5" - 13"
AT55F	5" - 13"
Range adjusts in 1" increments	
DESKADJCV, FULLADJCV	10" - 14"
Range adjusts in 1" increments	
TAD	Low Range
Range adjusts in 1/2" increments	
TAF	High Range
Range adjusts in 1" increments	
All measurements from top of seat rail to top of armrest. On AT55D and AT55F, add 1" for upholstered armrests. For less seat rail chairs, add 1-1/2" to above measurements (top of upper frame rail to top of armrest).	

ARM SELECTIONS

ARM TYPE (PAIR)

- AT99N Cantilever Non Locking Bolt On Arms..... No Charge ²
- FXFULL Conventional Fixed Ht.Full Lgth. Arms No Charge ³
- FXDESK Conventional Fixed Ht.Desk Lgth. Arms No Charge ³
- DESKADJCV Conventional Adjustable Ht Desk Lgth Arms(E0973)..... 195.00 ³
- FULLADJCV Conventional Adj. Ht.Full Lgth Arms (E0973)... 195.00 ³
- AT55D Cantilever Locking Desk Length
Bolt On Arms (E0973)..... 195.00 ²
- AT55F Cantilever Locking Full Length
Bolt On Arms (E0973)..... 195.00 ²
- TAD T-Arm Adjustable Height, Desk Length (E0973)225.00 ⁴
- TAF T-Arm Adjustable Height, Full Length (E0973).225.00 ⁴

T-ARM HEIGHT

- L6-10 Low Range 6" - 10" (T-Arm Height)..... No Charge
- H8-14 High Range 8" - 14" (T-Arm Height) No Charge

BACK SELECTIONS

BACK TYPE

- A1418 10° Bend Back w/Push Handles (14"-18").... No Charge
- 1014 Sportster Back w/o Push Handles (10"-14").. No Charge ⁵
- S1216 Sportster Back w/o Push Handles (12"-16").. No Charge ⁵
- A1216 Straight Back w/ Push Handles (12"-16")30.00
- S1418 Straight Back w/ Push Handles (14"-18")30.00
- FBH18 Fixed Back Height 18"..... No Charge
- FBH20 Fixed Back Height 20"..... No Charge
- RECLINER Recliner Frame No Charge ⁶

BACK UPHOLSTERY

- SLIPBACK Slip On Back Upholstery No Charge
- ADJTENS Adjustable Tension Back Upholstery (K0108) ...125.00
- L8601 Omit Back Upholstery Less 10.00

RECLINER BACK ^{7,8,9}

- REC18 Reclining Back 18" High (E1226)570.00
- REC20 Reclining Back 20" High (E1226)570.00
- REC24 Reclining Back 24" High (E1226)570.00
- RECEXT Reclining Headrest Pillow and Extension (E0966).....75.00 ⁹

FOOTNOTES

- ² Only available with FBH18 and FBH20 fixed back height.
- ³ Only available with CONV, Conventional Crossbraces and Seat Depths greater than 15".
- ⁴ Bolt-On Arm Attachment.
- ⁵ For push handles see Back Accessories section.
- ⁶ Must select Recliner Back options with this selection.
- ⁷ Must select when Recliner Frame (RECLINER) is selected.
- ⁸ 1" diameter anodized straight canes with push handles comes standard with back upholstery.

CASTER & WHEEL SELECTIONS ¹⁰

Front and Rear Seat-To-Floor Heights *Measured From Floor To Top Of Seat Rail					
	3"	5" 5ALUM AT901	6" AT902 6SKYWAY	6"	8" 1228 AT903
Fork	AT900			1259	
A1103 3" Fork	15.25" 15.75" 16.50" 17.00"	16.25" 16.75" 17.50" 18.00"	N/A	N/A	N/A
A1105 5" Fork	16.25" 16.75" 17.25" 17.50" 18.00" 18.50"	16.75" 17.25" 17.75" 18.00" 18.25" 18.50" 18.50" 19.00" 19.50" 19.50"	17.25" 17.75" 18.25" 18.50" 18.75" 19.00" 19.50" 20.00"	N/A	N/A
A1108 8" Fork	N/A	N/A	N/A	18.00" 18.50" 19.00" 19.25" 19.50" 19.75" 20.25" 20.75"	19.00" 19.50" 20.00" 20.25" 20.50" 20.75" 21.25" 21.75"
All Dimensions Are Approximate					

FRONT SEAT-TO-FLOOR HEIGHT ¹¹

- FSTF1525 15.25" Front Seat-To-Floor Height (K0056) 100.00
- FSTF1575 15.75" Front Seat-To-Floor Height (K0056) 100.00
- FSTF1625 16.25" Front Seat-To-Floor Height (K0056) 100.00
- FSTF165 16.5" Front Seat-To-Floor Height (K0056) 100.00
- FSTF1675 16.75" Front Seat-To-Floor Height (K0056) 100.00
- FSTF17 17" Front Seat-To-Floor Height..... No Charge
- FSTF1725 17.25" Front Seat-To-Floor Height..... No Charge
- FSTF175 17.5" Front Seat-To-Floor Height..... No Charge
- FSTF1775 17.75" Front Seat-To-Floor Height..... No Charge
- FSTF18 18" Front Seat-To-Floor Height..... No Charge
- FSTF1825 18.25" Front Seat-To-Floor Height..... No Charge
- FSTF185 18.5" Front Seat-To-Floor Height..... No Charge
- FSTF1875 18.75" Front Seat-To-Floor Height..... No Charge
- FSTF19 19" Front Seat-To-Floor Height..... No Charge
- FSTF1925 19.25" Front Seat-To-Floor Height..... No Charge
- FSTF1950 19.50" Front Seat-To-Floor Height..... No Charge
- FSTF1975 19.75" Front Seat-To-Floor Height..... No Charge
- FSTF20 20" Front Seat-To-Floor Height..... No Charge
- FSTF2025 20.25" Front Seat-To-Floor Height..... No Charge
- FSTF2050 20.50" Front Seat-To-Floor Height..... No Charge
- FSTF2075 20.75" Front Seat-To-Floor Height..... No Charge
- FSTF2125 21.25" Front Seat-To-Floor Height (K0056) 100.00

FOOTNOTES

- ⁹ RECEXT must be chosen with other Recliner Back option.
- ¹⁰ For folding frame without seat rails, subtract 1-1/2" from the seat-to-floor height to obtain distance from floor to top of side frame.
- ¹¹ Rear seat-to-floor height is set approximately 1" lower than the front seat-to-floor height.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare Compass XE Adult Positioning Folding Chair

Accessories and Options

PAXLE Permanent Axles Standard

CASTER OPTIONS (PAIR)

AT900 3" Roller Blade Casters38.00
 AT901 5" Composite Urethane Casters.....64.00
 5ALUM 5" Aluminum Urethane Casters.....92.00
 AT902 6" Composite Urethane 3 Spoke Casters..... No Charge
 6SKYWAY 6" Composite Urethane 6 Spoke Casters..... No Charge
 1259 6" X 2" Semi Pneumatic Casters (K0075).....60.00
 1228 8" X 1.75" Semi Pneumatic Casters Lt Grey60.00
 AT903 8" X 1.25" Composite Urethane Casters No Charge

FORK TYPE

A1103 3" Fork (Single Hole Fork) No Charge
 A1105 5" Fork (3 Hole Fork) No Charge
 A1108 8" Fork (2 Hole Fork) No Charge

WHEEL OPTIONS (PAIR) ¹²

20COM 20" Composite Wheels..... No Charge
 22COM 22" Composite Wheels..... No Charge
 24COM7 24" Composite 7-Spoke Wheels..... No Charge ¹³
 20SP 20" Spoke Wheels..... No Charge
 22SP 22" Spoke Wheels..... No Charge ¹³
 24SP 24" Spoke Wheels..... No Charge
 L8701 Omit Rear Wheels..... Less 50.00

TIRE OPTIONS (PAIR)

PNEU Pneumatic Tires60.00
 PNEUFF Flat Free Inserts (E2213)65.00 ¹⁴
 URETH Urethane (Treaded) Tires No Charge ¹⁵
 PRIMO Primo Tires60.00 ¹⁶
 KIKB KIK® Black Tires60.00 ¹⁶
 CRF123 Schwalbe 1 in Blk Flat Resistant Non Marking ¹⁷
Tire90.00

MOUNTING HANDRIMS (SPOKE WHEELS ONLY)

IN Mounting Inside Handrims (pair)..... No Charge ¹⁸
 OUT Mounting Outside Handrims (pair)..... No Charge

HANDRIM OPTIONS (PAIR)

ALUMPC Aluminum Anodized Handrims No Charge
 A1512 Black Plastic Coated Handrims.....72.00
 AT1225-8 Plastic Coated w/8 Vertical Projections240.00 ¹⁹
 AT1224-8 Plastic Coated w/8 Oblique Projections240.00
 AT1225-12 Plastic Coated w/12 Vertical Projections240.00 ¹⁶
 NFISTD Natural Fit Handrims Standard Grip295.00 ¹⁶
 NFITSP Natural Fit Handrims Super Grip.....340.00 ¹⁶
 L8702 Omit HandrimsLess 15.00 ¹⁶

AXLE & CAMBER SELECTIONS

AXLE PLATE OPTIONS (PAIR)

STDPT Standard Axle Plates No Charge
 MPOSAP Multi-Position Offset Axle Plates No Charge
 EXTPT Extended/Amputee Axle Plates (E0959).....72.00

AXLES (PAIR)

A1577 Quick Release Axles No Charge
 A1573 Quad Release Axles (K0108).....89.00

FOOTNOTES

- ¹² Spoke wheels are standard with dome button axle.
¹³ Not available on medium frame set at 14" deep.
¹⁴ Not available with 20SP, 22SP, or 24SP.
¹⁵ Only available on 24SP.
¹⁶ Only available on 24COM7 and 24SP.
¹⁷ The Schwalbe tire is a black high-pressure pneumatic tire that is non-marking, and includes a Smart Guard® puncture resistant layer. Schwalbe tire is only available with: 24SP.
¹⁸ Decrease overall width by ¾".
¹⁹ Only available on 22COM, 24COM7 and 24SP.

CAMBER ANGLE

C00 0° Camber.....No Charge
 C03 3° Camber.....No Charge
 C06 6° Camber.....No Charge

FRAME FINISH SELECTIONS

CONTEMPORARY COLORS

24P Wet Black..... No Charge
 71P Silver Vein No Charge
 60P Silver Metallic..... No Charge
 115P Black Prism..... No Charge
 105P Midnight Blue No Charge

TOO HOT COLORS

30P Sunny Yellow..... No Charge
 61P Electric Red..... No Charge
 119P Bubblegum Pink..... No Charge
 120P Grape Madness..... No Charge

WAY COOL COLORS

62P Electric Blue..... No Charge
 121P Lolly Pop Blue No Charge
 122P Cosmic Blue..... No Charge
 125P Emerald Green No Charge
 127P Grasshopper Green No Charge
 104P Electric Teal No Charge
 Touch Up Paint 20.00

POSITIONING STRAPS

SEAT POSITIONING STRAPS (EACH)

A1515 Push Button Buckle Seat Pos. Strap (E0978) 50.00
 A1321 Hook & Loop Seat Positioning Strap (E0978)..... 50.00

CHEST POSITIONING STRAPS (EACH)

9515 Push Button Chest Positioning Strap (E0960).... 100.00
 7321 Hook & Loop Chest Positioning Strap (E0960) . 100.00

ACCESSORIES

BACK ACCESSORIES

BP02 Back Pack..... 50.00 ²⁰
 1483 Non-Adjustable Rigidizing Handle Bar..... 85.00 ²¹
 RIGB Flip Down Rigidizer Bar 85.00 ²¹
 AASH Adj. Angle Stroller Handles (Outside Mount)..... 235.00 ²²
 AASHIN Adj. Angle Stroller Handles (Inside Mount)..... 235.00 ²²

SEAT ACCESSORIES

SP-01 Seat Pouch..... 35.00
 L8651 Omit Seat Upholstery..... Less 10.00

ARM ACCESSORIES (PAIR) ²³

A1571 Fabric Clothing Protector (K0108)..... 46.00
 1671-SM Rigid Side Guard, Small (K0108) 122.00
 1671-LG Rigid Side Guard, Large (K0108) 122.00
 1671-SMH Quad Rigid Side Guard, Small (K0108)..... 122.00
 1671-LGH Quad Rigid Side Guard, Large (K0108)..... 122.00

FOOTNOTES

- ²⁰ Not available with Reclining Back Option.
²¹ Not available on 19" or 20" widths. Not available with Reclining Back Option.
²² AASH release buttons are mounted on the outside and AASHIN release buttons are mounted on the inside.
²³ Only available with AT99N, AT55D, or AT55F.

Invacare Compass XE Adult Positioning Folding Chair

Accessories and Options

WHEEL LOCKS (PAIR)

- A1598 Hill Holder Clamp On Wheel Locks (E0974).....170.00 ²⁴
- A1567 Pull to Lock Clamp On Wheel Locks No Charge
- HUMW Hideaway Undermount Wheel Locks65.00 ²⁵
- 9691 Push to Lock Clamp On Wheel Locks No Charge

MISCELLANEOUS ACCESSORIES

- 3421 ¾ " To 1" Clamp Assembly.....21.00
- A1541 Caster Swivel Locks (pair) (K0073)80.00 ²⁶
- A1569 Crutch & Cane Holder (E2207)46.00
- A1307 Spoke Protectors (for 20", 22", 24", 26") (K0065).....90.00
- A1209 Tire Pump.....32.00
- A1580 Rear Anti-Tippers (pair) (E0971).....150.00
- A1576 Wheel Lock Extensions (6") (pair) (E0961).....65.00

MAX KNEE TO HEEL LENGTH

- KTH14 14" Knee to Heel No Charge
- KTH1516 15"-16" Knee to Heel No Charge
- KTH1618 16"-18" Knee to Heel No Charge

FOOTREST/LEGREST ACCESSORIES (PAIR)

- A1337 Calf Strap (K0038) 24.00
- FRMPRBK ChairArmor Black Neoprene Frame Protectors 65.00 ³³
- FRMPRB ChairArmor Blue Neoprene Frame Protectors..... 65.00 ³³
- FRMPRR ChairArmor Red Neoprene Frame Protectors..... 65.00 ³³
- FRMPRS ChairArmor Silver Neoprene Frame Protectors..... 65.00 ³³
- ALPT Longer Pivot And Slide Tubes 84.00 ³⁴
- A1600BK Heel Loops with Ankle Straps (E0951)..... 51.00
- HLOOP Heel Loops w/o Ankle Straps (E0951)..... 51.00
- 9342 "H" Strap (Each) (K0039)..... 51.00

Knee-to-Heel Range							
For Folding Frames with seat rails.							
70° Tapered		70°		70° MFX		Elevating	
MIN	MAX	MIN	MAX	MIN	MAX	MIN	MAX
16.25"	20.25"	14"	19"	8.25"	14"	15.75"	21.75"
90°		60°		60° MFX			
MIN	MAX	MIN	MAX	MIN	MAX		
5.75"	10.75"	15.25"	19"	11.5"	15.5"		

Range adjusts in 1" increments.

FRONT RIGGING SELECTIONS

FRONT RIGGING (PAIR) ²⁷

- SWING1319 60° Swingaway Footrests..... No Charge
- SWING914 60° Swingaway Footrests MFX No Charge
- 70MFX 70° Swingaway Footrests MFX No Charge ²⁸
- 70NTAPER 70° Swingaway Footrests..... No Charge
- 70STAPER 70° Tapered Footrests No Charge ²⁹
- 90SWING 90° Swingaway Footrests..... No Charge
- 70LIFTNT 70° Lift Off Non-Taper Footrests95.00
- 70LIFTMFX 70° Lift Off MFX Footrests95.00
- 90LIFTMFX 90° Lift Off MFX Footrests95.00
- AT5044 Elevating Legrests (E0990).....255.00
- ATL4A Center Pivot Articulating Legrests (K0053).....320.00

FOOTPLATE OPTION (PAIR)

- 1651 Flip-Up Composite Footplates No Charge ³⁰
- AT5543 Adjustable Angle Flip-Up Footplates(K0040)..... 165.00 ³¹
- A1564 Extra Large Aluminum Footplates59.00 ³²
- TUBULAR Individual Tubular Footplates.....59.00

FOOTNOTES

- ²⁴ Not available with URETH or KIKB tires.
- ²⁵ Not available with ATL4A Center Pivot Articulating Legrests.
- ²⁶ Only available on tubular footplates AT900, 5ALUM, AT901, AT902, 6SKYWAY.
- ²⁷ All front riggings are center pivot style. Footplates are attached directly to rigging tube when MFX style is selected.
- ²⁸ Standard with AT5543 adjustable angle footplate.
- ²⁹ 16" – 20" wide frames only. Standard with AT5543 adjustable angle footplate.
- ³⁰ Not available with 70MFX, 90SWING Swingaway footrests or 12" wide.
- ³¹ AT5543 footplates are dorsi-plantar adjustable.
- ³² Not available with 70STAPER or 12" or 15" wide.

SEATING CUSHION SELECTION ³⁵

- EC Absolute Cushion (E2601) 104.00 ³⁶
- ITSG InTouch™ Stabilite™ Gentle Contour Cushion (E2607)..... 340.00
- ITSM InTouch Stabilite Max Contour Cushion (E2605) 340.00
- ITFG InTouch Flovair™ Gentle Contour Cushion (E2607) 360.00
- ITFM InTouch Flovair Max Contour Cushion (E2607). 360.00

SEATING BACK SELECTION ³³

- PBR Personal Back 10™ Reg. Height (E2615)..... 565.00
- PBT Personal Back 10 Tall Height (E2615)..... 565.00
- PPR Personal Back 10 Plus Reg. Height (E2615) 645.00
- PPT Personal Back 10 Plus Tall Height (E2615) 645.00
- ITPR InTouch Propel™ Back (E2611) 438.00

**FOR ADDITIONAL SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027**

FOOTNOTES

- ³³ Not available with elevating legrests.
- ³⁴ 2" longer slide tube used on composite, adjustable angle or tubular footplates only, not available with 70MFX, 90SWING, 90LIFTMFX, SWING914, and ATL4A.
- ³⁵ For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.
- ³⁶ Available in even sizes only. See form 07-027 for additional information.

Invacare Compass XE Adult Positioning Folding Chair

Basic Chair Includes:

- Folding, Depth Adjustable Extended Frame
- Crossbrace Type Selection
- Transport Tie Down Options
- Seat Width (W12 - W20)
- Frame Depth Selection
- Seat Depth (D12 - D20)
- Arm Type (AT99N, FXDESK, FXFULL)
- Back Type Selection (A1418, 1014, S1216, FBH18, FBH20)
- Slip-On Back Upholstery
- Frame Finish
- Caster Options (AT902, AT903 or 6SKYWAY)
- Seat to Floor Height Selection (FSTF1725 – FSTF2075)
- Wheel Options
- Tire Options (URETH)
- Handrims Options (ALUMPC)
- Axle Type (A1577 or PAXLE)
- Axle Brackets (STDPT or MPOSAP)
- Camber Angle Selection
- Front Rigging (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER or 90SWING).
- Flip-Up Composite Footplates (1651)
- Seat Upholstery
- Wheel Locks (9691 or A1567)

Specifications

Frame Style	Folding
Seat Widths	12" – 20"
Seat Depths	12" – 20"
Overall Chair Width (closed)	Seat width plus 11.5"*
Overall Chair Width (open)	Seat width plus 8.5"*
Adjustable Height Back	14" – 18"
Back Angle Adjustment	80° to 110° in 5° increments
Camber Angles	0°, 3°, 6°
Tubing Outside Diameter	1"
Weight of Chair	19 lbs**/36lbs***
Warranty	Lifetime on frames and crossbraces
*0° Angles	
**18" x 16", AT902, no rear wheels, no footrests, no anti-tippers, no wheel locks, no arms.	
***18" x 16", AT902, 24SP, PNEU, 70NTAPER, seat upholstery, back upholstery, wheel locks, no arms.	

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.



Yes, you can.®

Invacare Corporation

USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

All rights reserved. Trademarks are identified by the symbols TM and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

KIK is a registered trademark of KIK Technology, Inc.

Specifications and prices are subject to change without notification.

©2009 Invacare Corporation

Form No. 04-094 Rev 0609