



Price List and Order Form

Effective January 15, 2008

For ease of ordering, contact customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Yes, you can.®

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

FRAME TYPE

PROX4S	X4 Folding Frame w/Sw Front End (K0005) ..	2,043.00	
PROX4F70	70 Degree Fixed Front Frame	2,314.00	¹
PROX4F80	80 Degree Fixed Front Frame	2,314.00	¹

FRAME WIDTH

W14	14" Wide	No Charge	²
W15	15" Wide	No Charge	²
W16	16" Wide	No Charge	²
W17	17" Wide	No Charge	²
W18	18" Wide	No Charge	
W19	19" Wide	No Charge	
W20	20" Wide (E2201).....	300.00	
W22	22" Wide (E2201).....	300.00	³

SEAT DEPTH

D14	14" Seat Depth	No Charge	²
D15	15" Seat Depth	No Charge	²
D16	16" Seat Depth	No Charge	²
D17	17" Seat Depth	No Charge	²
D18	18" Seat Depth	No Charge	
D19	19" Seat Depth	No Charge	
D20	20" Seat Depth (E2203).....	300.00	
D21	21" Seat Depth (E2203).....	300.00	³
D22	22" Seat Depth (E2204).....	300.00	³

HEAVY DUTY OPTION ⁴

HD Heavy Duty Option (K0108).....	400.00
-----------------------------------	--------

FRONT SEAT-TO-FLOOR

FSTF155	Front STF Height 15.5" (K0056)	95.00
FSTF16	Front STF Height 16" (K0056)	95.00
FSTF165	Front STF Height 16.5" (K0056)	95.00
FSTF17	Front STF Height 17"	No Charge
FSTF175	Front STF Height 17.5"	No Charge
FSTF18	Front STF Height 18"	No Charge
FSTF185	Front STF Height 18.5"	No Charge
FSTF19	Front STF Height 19"	No Charge
FSTF195	Front STF Height 19.5"	No Charge
FSTF20	Front STF Height 20"	No Charge
FSTF205	Front STF Height 20.5"	No Charge
FSTF21	Front STF Height 21"	No Charge
FSTF215	Front STF Height 21.5" (K0056)	95.00

Width & Depth Matrix									
Depth	14"	15"	16"	17"	18"	19"	20"	21"	22"
Width									
14"	X	X	X	X	S	S	S		
15"	X	X	X	X	S	S	S		
16"	X	X	X	X	S	S	S		
17"	X	X	X	X	S	S	S		
18"	X	X	X	X	X*	X*	X*	S	S
19"	S	S	X	X	X*	X*	X*	S	S
20"	S	S	X	X	X*	X*	X*	*	*
22"			S	S	S	S	S	S	*

X- Sizes available S- Sizes available through Specials
 *-Sizes available with HD

Front Seat-To-Floor Heights					
Fork	Spacer	CASTER			
		3" AT900	4" 4SPLT 4ALUM2	5" 5SOFT AT901 5ALUM2	6" AT902 6ALUM2
BRG1 (Short)	None	15.5 16.0	15.5 16.0 16.5	16.5 17.0	17.5
	1"	16.5 17.0	16.5 17.0 17.5	17.5 18.0	18.5
	2"	17.5 18.0	17.5 18.0 18.5	18.5 19.0	19.5
BRG2 (Long)	None	17.5 18.0	17.5 18.0 18.5	18.0 18.5 19.0	18.5 19.0 19.5
	1"	18.5 19.0	18.5 19.0 19.5	19.0 19.5 20.0	19.5 20.0
	2"	19.5 20.0	19.5 20.0 20.5	20.0 20.5 21.0	20.5 21.0 21.5
SOFR (Suspension)	None	16.5 17.0	16.5 17.0 17.5	17.5 18.0	18.5
	1"	17.5 18.0	17.5 18.0 18.5	18.5 19.0	19.5
	2"	NA	NA	NA	NA

FOOTNOTES

- ¹ Not available with HD - Heavy Duty, TRRO or TRBKTS.
- ² Not available with HD - Heavy Duty Option.
- ³ Must order HD - Heavy Duty Option.
- ⁴ Only available on the PROX4S, X4 Folding Frame w/Swingaway Front Frame.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

SUSPENSION SELECTIONS

TRANSPORT TIE DOWN

TRRO	Transport Ready Option.....	300.00
TRBKTS	Wheelchair Transport Brackets.....	250.00

FRONT SUSPENSION

SOFR	Frog Legs™ Front Suspension Forks (E1015).....	315.00
------	--	--------

FRONT SUSPENSION SETTING

SET160	Soft Setting (Max. User Weight 160 lbs.).....	No Charge
SET300	Firm Setting (Max. User Weight 250 lbs.).....	Standard

TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*** Transport Tie Down option is not retrofittable to existing models and is not field serviceable.**

ARM SELECTIONS

ARM TYPE (PAIR)

AT99	Swingaway Padded Arms	No Charge	⁵
CAD	Conventional Adj Hgt Desk Arms (E0973).....	195.00	
CAF	Conventional Adj Hgt Full Arms (E0973).....	195.00	
TAD	T-Arms, Adj. Height Desk Length (E0973)	225.00	⁶
TAF	T-Arms, Adj. Height Full Length (E0973)	225.00	⁶
TAF	T-Arms, Adj. Height Full Length	225.00	
LESSARMS	Less Arms	No Charge	

T-ARM HEIGHT

L6-10	Low Range 6" – 10" (T-Arm Height).....	No Charge
H8-14	High Range 8" – 14" (T-Arm Height).....	No Charge

Arm Heights

AT99		8"-10"
Range adjusts in 1" increments		
CAD		10"-14"
CAF		10"-14"
Range adjusts in 1" increments		
TAD	Low Range	6"-10"
Range adjusts in 1/2" increments		
TAF	High Range	8"-14"
Range adjusts in 1" increments		
All measurements from top of seat rail to top of armrest. On AT99 add 1" for upholstered arm pads. For less seat rail chairs, add 1/2" to above measurements (top of upper frame rail to top of armrest).		

FOOTNOTES

⁵ Not available with Heavy Duty Option or Angle Adjustable Backs (AAB1216 & AAB1720)

⁶ Rigid sideguards included with the T-Arms. Bolt-on attachment.

BACK SELECTIONS

BACK TYPE

BACK1014	Adj. Back w/o Push Handles (10" -14").....	No Charge
BACK1214	Adj. Back w/o Push Handles (12" -14").....	No Charge
BACK1418	Adj. Back w/o Push Handles (14" -18").....	No Charge
BACK1620	Adj. Back w/o Push Handles (16" -20").....	No Charge
PH1216	Adj. Back w/ Push Handles (12" -16").....	30.00
PH1620	Adj. Back w/ Push Handles (16" -20").....	30.00
A1418	Adj. Hgt 10° Bendback w/Push Hndls (14"-18")...	30.00
A1620	Adj. Hgt 10° Bendback w/Push Hndls (16"-20")...	30.00
AAB1216	Adj. Angle Back w/Push Handles (12" -16").....	165.00 ⁷
AAB1720	Adj. Angle Back w/Push Handles (17" -20").....	165.00 ⁷

BACK HEIGHT

ABH10	10" High.....	No Charge
ABH11	11" High.....	No Charge
ABH12	12" High.....	No Charge
ABH13	13" High.....	No Charge
ABH14	14" High.....	No Charge
ABH15	15" High.....	No Charge
ABH16	16" High.....	No Charge
ABH17	17" High.....	No Charge
ABH18	18" High.....	No Charge
ABH19	19" High.....	No Charge
ABH20	20" High.....	No Charge

BACK UPHOLSTERY TYPE

SLIPBACK	Slip on Back Upholstery.....	No Charge
ADJTENS	Adjustable Tension Back Upholstery.....	125.00
L8601	Omit Back Upholstery.....	Less 10.00

FRAME FINISH SELECTIONS

CONTEMPORARY COLORS

24P Wet Black	No Charge
71P Silver Vein	No Charge
60P Silver Metallic	No Charge
115P	Black Prism..... No Charge
105P	Midnight Blue..... No Charge

TOO HOT COLOR

30P Sunny Yellow	No Charge
61P Electric Red	No Charge
119P	Bubblegum Pink..... No Charge
120P	Grape Madness..... No Charge

WAY COOL COLOR

62P Electric Blue	No Charge
121P	Lolly Pop Blue..... No Charge
122P	Cosmic Blue..... No Charge
125P	Emerald Green..... No Charge
127P	Grasshopper Green..... No Charge
104P	Electric Teal..... No Charge

FOOTNOTES

⁷ Not available with HD-Heavy Duty Option, TRRO, TRBKTS or ADJTENS.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

CASTER SELECTIONS

CASTER OPTIONS (PAIR)

AT900	3" Roller Casters.....	38.00	⁸
4ALUM2	4" Frog Legs® Alum Soft Roll Casters.....	92.00	
4SPLT	4" Comp Urethane Light Up Casters.....	No Charge	⁸
AT901	5" Composite Urethane Casters.....	No Charge	⁸
5SOFT	5" Semi Pneumatic Casters.....	60.00	⁸
5ALUM2	5" Frog Legs® Aluminum Soft Roll Casters.....	92.00	
AT902	6" Comp Urethane 3 Spoke Casters.....	No Charge	⁸
6ALUM2	6" Frog Legs® Aluminum Soft Roll Casters.....	92.00	

FORK TYPE (NON-SUSPENSION) (PAIR)

BRG1	Fork (3-Hole Short Fork.....	No Charge
BRG2	Fork (3-Hole Long Fork).....	No Charge

FORK SPACER

SPCR0	No Fork Spacer.....	No Charge
SPCR1	1" Fork Spacer.....	No Charge
SPCR2	2" Fork Spacer.....	No Charge

REAR WHEEL SELECTIONS ¹⁰

Rear Seat-To-Floor Ranges	
SEAT-TO-FLOOR HEIGHT	REAR WHEEL
13.5"-18"	20COM, 20SP
14.5"-19"	22COM, 22SP
15.5"-20"	24COM7, 24SP, 24HSP, 24SPIN
16"-20.5"	25SPIN, 559TS, 26SP, 26HPS, 26SPIN

WHEEL OPTIONS (PAIR) ¹¹

20COM	20" Composite Wheels.....	No Charge	¹²
22COM	22" Composite Wheels.....	No Charge	¹²
24COM7	24" Composite-7 Spoke Wheels.....	No Charge	¹³
20SP	20" Spoke Wheels.....	No Charge	^{8,14}
22SP	22" Spoke Wheels.....	No Charge	^{8,14}
24SP	24" Spoke Wheels.....	No Charge	^{8,15}
26SP	24" Spoke Wheels.....	48.00	^{8,16}
559TS	559 Turbo Spoke.....	117.00	^{8,17}
24HPS	24" High Performance Spoke Wheels.....	117.00	^{8,15}
26HPS	26" High Performance Spoke Wheels.....	117.00	^{8,16}
SPIN24	24" Spinergy® Spoke Wheels.....	695.00	^{8,18}
SPIN25	25" Spinergy® Spoke Wheels.....	695.00	^{8,17}
SPIN26	26" Spinergy® Spoke Wheels.....	695.00	^{8,18}
L8701	Omit Rear Wheels.....	Less 50.00	^{8,19}

FOOTNOTES

- ⁸ Not available with HD-Heavy Duty Option.
⁹ Not available with SOFR-F.
¹⁰ Rear Seat-To-Floor is factory set at 1"-1.5" lower than Front Seat-To-Floor.
¹¹ Spoke wheels come standard with dome, quick-release button axle.
¹² Available with PNEU, PNEUFF or URETH tires.
¹³ Available with PNEU, PNEUFF, URETH or KIKB tires.
¹⁴ Available with PNEU or PNEUFF tires.
¹⁵ Available with PNEU, PNEUFF, PRIMO, or KIKB tires.
¹⁶ Available with PNEU or PRIMO tires.
¹⁷ Available with PNEU or PRIMOTRD tires.
¹⁸ Available with PNEU tires.

TIRE OPTIONS (PAIR)

PNEU	Pneumatic Tires.....	60.00
PNEUFF	flat free Inserts (K0064).....	65.00
URETH	Urethane (treaded) Tires.....	No Charge
PRIMO	Primo Tires.....	60.00
PRIMOTRD	Primo 559 Treaded Pneumatic.....	60.00
KIKB	KIK® Black Tires.....	60.00

HANDRIM OPTIONS (PAIR)

ALUMPC	Aluminum Anodized Handrim.....	No Charge
TIPC	Titanium Handrims.....	400.00
A1512	Black Plastic Coated Handrims (K0059).....	72.00
NFITSD	Natural Fit Handrims Standard Grip.....	295.00
NFITSP	Natural Fit Handrims Super Grip.....	340.00
L8702	Omit Handrims.....	Less 15.00

MOUNTING HANDRIMS (SPOKE WHEELS ONLY) ⁸

IN	Mount Inside Handrim.....	No Charge
OUT	Mount Outside Handrim.....	No Charge

AXLES (PAIR)

A1577	Quick Release Axles.....	No Charge
A1573	Quad Release Axles.....	89.00
PAXLE	Permanent Axles.....	No Charge

CAMBER ANGLE

C000°	Camber	No Charge
C033°	Camber	No Charge
C066°	Camber	No Charge

POSITIONING STRAPS/ACCESSORIES

POSITIONING STRAPS (EACH)

A1515	Push Button Buckle Seat Pos. Strap (E0978).....	50.00
A1321BK	Hook & Loop Seat Pos. Strap (E0978).....	40.00

BACK ACCESSORIES

BP02	Back Pack.....	50.00
------	----------------	-------

SEAT ACCESSORIES

SP01	Seat Pouch.....	35.00
CR	Cushion Rigidizer (for use w/sling seat) (E2618).....	37.00

ARM ACCESSORIES (PAIR) ²²

A1571	Fabric Clothing Protector.....	46.00
1671-SM	Rigid Side Guard, Small.....	122.00
1671-LG	Rigid Side Guard, Large.....	122.00
1671-SMH	Quad Rigid Side Guard, Small.....	122.00
1671-LGH	Quad Rigid Side Guard, Large.....	122.00

PERFORMANCE DETAIL PACKAGES

SPORT	ProSpin Sport Performance Detail Package.....	65.00
-------	---	-------

FOOTNOTES

- ¹⁹ When choosing omit rear wheels the frame will be constructed for 24" rear wheels.
²⁰ Not available with 20COM, 22COM, 24COM7, 20SP, 22SP, 559TS.
²¹ Not available with odd seat widths.
²² Available with AT99 arms only.
²³ Sport Package includes impact guards, rear frame protectors and clear coat back cane brackets.

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information go to www.invacare.com.**

MISCELLANEOUS ACCESSORIES

A1307	Spoke Protector (For 20", 22", 24", 26")(pr)	90.00
1209	Tire Pump	35.00
A1580	Rear Anti-Tippers (pair).....	60.00
A1541	Caster Swivel Locks	80.00
EXTPT	Extended Amputee Axle Plate	150.00 ⁸
6WLEXT	6" Wheel Lock Extensions	39.00

WHEEL LOCK SELECTIONS

WHEEL LOCKS

A1598	Hill Holder Clamp-On Wheel Locks	77.00 ²⁴
WLPULL	Pull To Lock Clamp-On Wheel Locks	No Charge
WLPUSH	Push To Lock Wheel Locks	No Charge
HUMV	Hiway Undermount Wheel Locks	65.00 ^{8,24,25}
L8704	Omit Wheel Locks.....	Less 10.00

FRONT RIGGING SELECTIONS

FRONT RIGGING (PAIR)

70SSNTAPER	70° Non-Tapered Swing In/Swing Out Ftrsts	No Charge
60SSNTAPER	60° Swing In/Swing Out Footrests.....	No Charge
70SSTAPER	70° Tapered Swing In/Swing Out Ftrsts	No Charge
ELRSS	Elevating Legrests	255.00

FOOTPLATE TYPE (PROX4S) (PAIR)

1651	Flip-Up Composite Footplates	No Charge
AT5543	Adjustable Angle Flip-Up Footplates (K0040)...	165.00 ²⁶
TUBULAR	Individual Tubular Footplates	59.00

FOOTPLATE TYPE (PROX4F70 & PROX4F80) (PAIR)

AT5545	One Pc. Adjustable Angle Flip Up Footplates.....	165.00 ⁸
--------	--	---------------------

FOOTRST/LEGREST ACCESSORIES

A1337	Calf Strap (each) (K0038)	24.00
IMPACTG	Impact Guards (pair)	51.00
ALPT	Longer Pivot and Slide Tube (pair).....	89.00 ²⁷
A9342	"H" Strap (each) (K0039)	51.00
HLOOP	Heel Loops w/o Ankle Straps (pair) (E0951)	51.00

Knee-To-Heel Range (PROX4S)		
Front Rigging Type	Min.	Max.
60°	14.00	18.75
70°	13.75	18.75
70° Tapered	14.50	19.50
Range adjusts in 1" increments		
Measurements taken from top of seat rail to back of footrests.		

Knee-To-Heel Range (PROX4F70 & PROX4F80)		
Front Rigging Type	Min.	Max.
Adjustable Angle Flip Up	14.00	18.00
1 Pc. Tubular Flip Up	14.00	18.00
Range adjusts in 1/4" increments		
Measurements taken from top of seat rail to back of footrests.		

FOOTNOTES

- ²⁴ Not available with URETH tires.
- ²⁵ Not available with D14, D15, or D16 seat depths.
- ²⁶ Doris/Plantar angle adjustable.
- ²⁷ Not available with AT5545 & AT 5540.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

Invacare® ProSPIN™ X4 Ultralightweight Active Performance Folder

Accessories and Options

SEATING & POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027**

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Basic Chair Includes

- Folding Curved Frame w/Swingaway Front
- Seat Width (W14" - W19")
- Seat Depth (D14" - D19")
- Swingaway Padded Arms
- Front Seat To Floor Height Selection (17"-21")
- Adjustable Height Back Selection w/o Push Handles
- Slip On Back Upholstery
- Powder Coated Frame Finish
- Caster Selection (SPLT4, AT902, AT901)
- Fork Type (BRG1 or BRG2)
- Fork Spacer Type (SPCR0, SPCR1, SPCR2)
- Wheel Options (20Sp, 22SP, 24SP, 20COM, 22COM, 24COM7)
- Urethane (treaded) Tire (URETH)
- Aluminum Anodized Handrims (ALUMPC)
- Axle Type (PAXLE or A1577)
- Wheel Lock Selections (Push or Pull to Lock)
- Camber Angle Selection
- Handrim Options (AHDRM or CHDRM)
- Swingaway Front Rigging
- Flip Up Composite Footplates

Available through Specials Department

All Specials are subject to final approval. Some special quotes are not compatible for all frame configurations. For your customization requests, please contact Customer Service for details and prices.

1. Oxygen Holder
2. 80° Swingaway Front Riggings

Specifications

• Frame Style	Folding w/Swingaway Front End
• Seat Widths	14" - 20" & 22"
• Seat Depths	14" - 22"
• Back Heights	10" up to 20"
•	(See Back Type for Range)
• Adjustable Angle Back	0°, 5°, 10°, 15°
• Camber Angles	0°, 3°, 6°, 9°
• Tubing O.D.	0° - 45°
• Weight of Chair	15 lbs./30 lbs**
• Weight Limit	250 lbs.
• Heavy Duty	350 lbs.
• Warranty	Lifetime on frames and crossbraces
•	

*16"x16", BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, no rear wheels, no footrests, no anti-tippers, no wheel locks, no arms

**16"x18", AT99, BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, SPIN24, PRIMO, ALUMPC, A1577, 70SSNTAPER, 1651, HLOOP, IMPACTG, 1671LG, WLPULL.



Invacare Corporation

USA

One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

KIK is a registered trademark of KIK Technology, Inc.

Spinergy is a registered trademark of Spinergy, Inc.

Frog Legs is registered trademark of Frog Legs, Inc.

Specifications and prices are subject to change without notification.

©2008 Invacare Corporation

Form No. 04-097