



Yes, you can.[®]

Price List and Order Form

Effective July 1, 2008

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

CHAIR TYPE ¹

- SOLARA2G Solara Tilt-in-Space wheelchair (E1161)..... 2,952.00

SEAT WIDTH RANGE

- SWR1216 12" – 16" Seat Width Range..... No Charge
- SWR1418 14" – 18" Seat Width Range..... No Charge
- SWR1620 16" – 20" Seat Width Range..... No Charge
- SWR1822 18" – 22" Seat Width Range..... No Charge

FRAME WIDTH

- W12 12" Wide..... No Charge
- W13 13" Wide..... No Charge
- W14 14" Wide..... No Charge
- W15 15" Wide..... No Charge
- W16 16" Wide..... No Charge
- W17 17" Wide..... No Charge
- W18 18" Wide..... No Charge
- W19 19" Wide..... No Charge
- W20 20" Wide (E2201)..... 300.00
- W21 21" Wide (E2201)..... 300.00
- W22 22" Wide (E2201)..... 300.00

SEAT FRAME LENGTH

- SDR1216 12" – 16" Seat Depth Range..... No Charge
- SDR1519 15" – 19" Seat Depth Range..... No Charge
- SDR1822 18" – 22" Seat Depth Range..... No Charge

SEAT DEPTH

- D12 12" Seat Depth..... No Charge
- D13 13" Seat Depth..... No Charge
- D14 14" Seat Depth..... No Charge
- D15 15" Seat Depth..... No Charge
- D16 16" Seat Depth..... No Charge
- D17 17" Seat Depth..... No Charge
- D18 18" Seat Depth..... No Charge
- D19 19" Seat Depth..... No Charge
- D20 20" Seat Depth (E2203)..... 300.00
- D21 21" Seat Depth (E2203)..... 300.00
- D22 22" Seat Depth (E2204)..... 300.00

CUSTOM PACKAGES ²

- HDPG Heavy Duty Package 375.00

FOOTNOTES

- ¹ A seating system must be purchased for use with this wheelchair. No seating system provided.
- ² Reinforced frame. Weight limit of 350 lbs.

SEAT STYLE ³

- FLSEAT Flush Mounted Aluminum Seat Pan (K0108)..... 95.00
- DPSEAT Drop Mounted Aluminum Seat Pan (K0108) 95.00 ⁴

TRANSPORT/VENT TRAY SELECTION

TRANSPORT TIE DOWN ^{5 6}

- TRRO Transport Ready Option (K0108)..... 250.00
- TRBKTS Transport Brackets..... 250.00

*TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*** Transport Tie Down option is not retrofittable to existing models and is not field serviceable.**

VENT TRAY OPTION ⁷

- AVENT Ventilator Tray (each) (E1029) 295.00 ⁸
- ATBOX Ventilator Battery Box (each) 97.00 ⁹

FOOTNOTES

- ³ FLSEAT or DPSEAT are solid seat pans and do not come with the chair; must be selected if needed.
- ⁴ DPSEAT, drop mounted seat pan accommodates 2" lower seat to floor heights. Use a cushion not less than 3" thick.
- ⁵ TRRO/TRBKTS not available with Heavy Duty Package or Reclining Back.
- ⁶ **Transport Tie-Down options not retrofittable to existing models, not field serviceable.**
- ⁷ Available with 18" and above Seat-to-Floor Heights only.
- ⁸ Removable from back canes. Only available on chair widths 16" or wider.
- ⁹ Non-Removable from back canes.

**Invacare® Solara™ 2G – Adult Tilt-in-Space Wheelchair
Accessories and Options**

ARM & BACK SELECTION

FRAME FINISH COLORS CONTINUED

ARM TYPE (PAIR) ¹⁰

- TAD T-Arm Adjustable Height, Desk Length (E0973) 225.00
- TAF T-Arm Adjustable Height, Full Length (E0973) 225.00
- FXDESK Conventional Fixed Ht, Desk Lgth Arms No Charge
- FXFULL Conventional Fixed Ht, Full Lgth Arms No Charge
- ADJDESK Conventional Adj Ht, Desk Lgth Arms (E0973) 195.00
- ADJFULL Conventional Adj Ht, Full Lgth Arms (E0973) 195.00
- FADJDESK Two-Point Flip-Back Adjustable Ht Desk (E0973) 245.00
- FADJFULL Two-Point Flip-Back Adjustable Ht Full (E0973) 245.00

T-ARM HEIGHT

- L8-12 Low Range 8" – 12" No Charge
- H10-16 Height Range 10" – 16" No Charge

Arm Heights

Adjustable-Height Dual-Point
8.5" – 12.5"
Conventional only.
Range Adjusts in 1" increments.
Measurements taken from top of seat rail to top of armrest.

Call Customer Service for specials for additional arm heights.

T-Arm
8" – 12" (0.5" increments)
10" – 16" (1" increments)
Measurements taken from top of seat rail to top of armrest.

BACK TYPE ^{11, 12}

- RECLINER Reclining Back (E1226) 810.00

BACK HEIGHT

- BH20 Back Height 20" No Charge
- BH24 Back Height 24" No Charge
- ABH1720 Adjustable Height Back Canes (17" – 20") No Charge ¹³

FOLD OPTION BACK

- FXBACK Non-Fold Down Back Canes **Standard**
- FLDBACK Fold-Down Back Canes No Charge

HEADREST & EXTENSION

- RECEXT Reclining Headrest Pillow and Extension (E0966) 75.00

FRAME FINISH COLORS

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal No Charge
- 121P Lolly Pop Blue No Charge
- 122P Cosmic Blue No Charge
- 125P Emerald Green No Charge
- 127P Grasshopper Green No Charge
- Touch Up Paint 20.00 ¹⁴

SEAT-TO-FLOOR HEIGHT

- STF145 Seat-To-Floor 14.5" No Charge
- STF15 Seat-To-Floor 15" No Charge
- STF155 Seat-To-Floor 15.5" No Charge
- STF16 Seat-To-Floor 16" No Charge
- STF165 Seat-To-Floor 16.5" No Charge
- STF17 Seat-To-Floor 17" No Charge
- STF175 Seat-To-Floor 17.5" No Charge
- STF18 Seat-To-Floor 18" No Charge
- STF185 Seat-To-Floor 18.5" No Charge
- STF19 Seat-To-Floor 19" No Charge

CASTER & REAR WHEEL SELECTIONS

CASTER OPTIONS (PAIR)

- 4 SPLT 4" Composite Urethane Light Up No Charge
- SSOFT 5" Semi-Pneumatic casters No Charge
- 5ALUM 5" Aluminum Urethane Casters 92.00
- AT902 6" Composite Urethane 3-Spoke Casters No Charge
- 1259 6" x 2" Semi-Pneumatic Casters 60.00
- 1228 8" x 1.75" Semi-Pneumatic Casters Lt Grey 60.00
- AT903 8" x 1.25" Composite Urethane Casters No Charge

Rear Wheel	Front And Rear Seat-to-Floor Heights ¹⁵				
	*Measured From Floor To Top Of Seat Rail				
	4" 4SPLT	5" SSOFT 5ALUM	6" AT902	6" 1259	8" 1228 AT903
12"	16.00"	16.00" 16.50"	16.00" 16.50" 17.00"	16.00" 16.50" 17.50" 18.00"	17.00" 17.50"
18"	14.50" 15.00" 15.50" 16.00"	14.50" 15.00" 15.50" 16.00"	15.50" 16.00"	16.00"	18.50" 19.00"
20"	15.00" 15.50" 16.00"	15.00" 15.50" 16.00" 16.50"	15.50" 16.00" 16.50" 17.00"	16.00" 16.50"	17.00"
22"	N/A	16.00" 16.50"	16.00" 16.50" 17.00"	16.00" 16.50" 17.50" 18.00"	17.00" 17.50"
24"	N/A	N/A	17.00"	17.50" 18.00"	17.00" 17.50" 18.50" 19.00"
All Dimensions Are Approximate					

FOOTNOTES

¹⁰ With **TRRO** the only arms available are the TAF or TAD.
¹¹ Comes with back upholstery and adjustable angle rigidizing handles bar. Must choose **RECEXT** when selecting **RECLINER**.
¹² Not available on 21", 22" seat width, **TRRO**, **TRBKTS** or a Heavy Duty Package. Standard with **AASH** adjustable Rigidizing Handle Bar.
¹³ Must choose the **AASH** when ordering the **ABH1720**. The **ABH1720** is not available with the 1483 Non-Adjustable Rigidizing Handle Bar.

FOOTNOTES

¹⁴ Touch Up Paint not available in all colors. .
¹⁵ Refer to Seat-to-Floor chart.

Invacare® Solara™ 2G – Adult Tilt-in-Space Wheelchair

Accessories and Options

WHEEL OPTIONS (PAIR)

- 12COM 12" Composite Wheels No Charge ¹⁶
- 18COM 18" Composite Wheels No Charge ¹⁶
- 20COM 20" Composite Wheels No Charge
- 22COM 22" Composite Wheels No Charge
- 24COM 24" Composite 7-Spoke Wheels..... No Charge

AXLE TYPE (PAIR)

- A1577 Quick Release Axles No Charge
- PAXLE Permanent Axles..... **Standard**

TIRE OPTIONS (PAIR)

- PNEU Pneumatic Tires 60.00
- PNEUFF Flat Free Inserts (E2213)..... 65.00
- URETH Urethane (Treaded) Tires No Charge

HANDRIM OPTIONS (PAIR)

- AHDRM Aluminum Handrims No Charge
- CHDRM Chrome Plated Handrims No Charge
- 1512 Plastic Coated Aluminum Handrims..... 72.00
- AT1225-8 Plastic Coated w/8 Vertical Projections..... 240.00 ¹⁷
- AT1224-8 Plastic Coated w/8 Oblique Projections..... 240.00 ¹⁷
- AT1225-12 Plastic Coated w/12 Vertical Projections..... 240.00 ¹⁷
- L8702 Omit Handrims Less 15.00

POSITIONING STRAPS

POSITIONING STRAPS (EACH)

- 1515 Push Button Buckle Seat Positioning Strap (E0978) 50.00
- 1321BK Hook & Loop seat Positioning Strap (E0978)..... 50.00
- 1311BK Airline Buckle Seat Positioning Strap (E0978)..... 50.00
- 9515 Push Button Chest Positioning Strap (E0960)..... 100.00
- 7321BK Hook & Loop Chest Positioning Strap (E0960)..... 100.00
- 7311BK Airline Buckle Chest Positioning Strap (E0960)..... 100.00

ACCESSORIES & WHEEL LOCKS

BACK ACCESSORIES

- ADJRGD Adjustable Width Back Rigidizing Bar 85.00
- 1483 Non-Adjustable Stroller Handles 85.00
- AASH Adj. Angle Stroller Handles (Outside Mount)..... 235.00 ¹⁸
- AASHIN Adj. Angle Stroller Handles (Inside Mount)..... 235.00 ¹⁸

MISCELLANEOUS ACCESSORIES

- TLL Tilt Lock Out Levers (pair) No Charge
- TLC Tilt Lock Out Stops (pair) 30.00 ¹⁹
- 1209 Tire Pump 35.00
- P750 Oxygen Holder(E2208) 128.00 ^{20,21,22}

WHEEL LOCKS

- PULLOCK Pull To Lock Wheel Locks (pair)..... No Charge
- PUSHLOCK Push To Lock Wheel Locks (pair) No Charge
- EXTLOCK Push To Lock Wheel Locks w/6" Bolt-On Ext(pr)(E0961) .. 65.00 ²³
- FOOTHUB Foot Operated Hublock 200.00 ²³
- ATHUB Cane Mount Hub Lock 50.00 ²³
- USEHUB Seat Rail Mount Hub Lock..... 50.00 ²³

FOOTNOTES

- ¹⁶ Not available with urethane tires.
- ¹⁷ Only available with 24COM7.
- ¹⁸ AASH release buttons are mounted on the outside and AASHIN release buttons are mounted on the inside.
- ¹⁹ Additional pair of tilt out stops, one pair standard with chair.
- ²⁰ Oxygen holder is designed to only carry one "D" size tank.
- ²¹ Not available with BH20, RECLINER, TRRO or TRBKTS.
- ²² Includes Adjustable Width Back Rigidizer Bar.
- ²³ Not available with vent tray.

FRONT RIGGING

FRONT RIGGING (PAIR)

- 70LIFTNT 70° Lift Off Non-Taper Footrests..... 95.00
- 70LIFTMFX 70° Lift Off MFX Footrests..... 95.00
- 90LIFTMFX 90° Lift Off MFX Footrests..... 95.00
- SWING1319 60° Swingaway Footrests No Charge
- SWING914 60° Swingaway Footrests MFX..... No Charge
- 70NTAPER 70° Swingaway Footrests No Charge
- 70MFX 70° Swingaway Footrests MFX..... No Charge
- 70STAPER 70° Tapered Footrests No Charge
- 90SWING 90° Swingaway Footrests (pair) No Charge ²⁴
- AT5044 Elevating Legrests (E0990) 255.00
- BILAT Bi-Lateral Contracture Footrests 120.00
- ATLA Center Pivot Articulating Legrests (K0053)..... 320.00 ²⁵
- 110CP Contracture Platform 280.00 ²⁵

Knee-To-Heel Range		
Front Rigging Type	Min.	Max
60°	13.50	17.50
60° MFX	9.00	12.50
70°	13.50	16.50
70° Tapered	12.75	16.25
70° MFX	5.75	11.75
90°	4.50	7.50
Elevating	14.00	20.00
Range adjust in 1" increments Measurement taken from top of seat rail to back of footrests		

FOOTPLATE TYPE (PAIR)

- 1651 Flip-Up Composite Footplates No Charge
- A1562 Aluminum Footplates 59.00
- A1564 Extra Large Aluminum Footplates (K0041) 59.00
- AT5543 Adjustable Angle Flip-Up Footplates (K0040)..... 165.00

ADJUSTABLE ANGLE FOOTPLATE RIGHT

(BI-LATERAL CONTRACTURE FOOTREST ONLY) (EACH)

- AAF325R Adj Angle Footplate Right 3.25"W (K0040) 83.00
- AAF425R Adj Angle Footplate Right 4.25"W (K0040) 83.00
- AAF525R Adj Angle Footplate Right 5.25"W (K0040) 83.00
- AAF625R Adj Angle Footplate Right 6.25"W (K0040) 83.00
- AAF725R Adj Angle Footplate Right 7.25"W (K0040) 83.00

ADJUSTABLE ANGLE FOOTPLATE LEFT (EACH)

- AAF325L Adj Angle Footplate Right 3.25"W (K0040) 83.00
- AAF425L Adj Angle Footplate Right 4.25"W (K0040) 83.00
- AAF525L Adj Angle Footplate Right 5.25"W (K0040) 83.00
- AAF625L Adj Angle Footplate Right 6.25"W (K0040) 83.00
- AAF725L Adj Angle Footplate Right 7.25"W (K0040) 83.00

FOOTREST/LEGREST ACCESSORIES (PAIR)

- 1337 Calf Strap (K0038) 24.00
- 9342 "H" Strap (K0039) 56.00
- IMPACTG Impact Guards (E0951) 51.00
- ALPT Longer Pivot And Slide Tubes (K0039)..... 84.00 ²⁶
- ALPT3 3" Extension For 90° Footrests..... 31.00

HEEL LOOP/TOE LOOP (PAIR)

- A1600BK Heel Loops With Ankle straps (E0951) 51.00
- HLOOP Heel Loops Without Ankle Straps (E0951) 51.00
- 1345 Toe Loops (E0952) 39.00

FOOTNOTES

- ²⁴ Not available with 8" casters.
- ²⁵ Not available on seat widths less than 16".
- ²⁶ Not available with 70STAPER.

Invacare® Solara™ 2G – Adult Tilt-in-Space Wheelchair Accessories and Options

SEATING

SEATING CUSHION SELECTION ²⁷

<input type="checkbox"/> EC	Absolute Cushion (E2601).....	104.00	²⁸
<input type="checkbox"/> ITSG	InTouch™ Stabilite™ Gentle Contour Cushion (E2607)	340.00	
<input type="checkbox"/> ITSM	InTouch Stabilite Max Contour Cushion (E2605)	340.00	
<input type="checkbox"/> ITFG	InTouch Flovair™ Gentle Contour Cushion (E2607)	360.00	
<input type="checkbox"/> ITFM	InTouch Flovair Max Contour Cushion (E2607).	360.00	

SEATING BACK SELECTION ²⁷

<input type="checkbox"/> PBR	Personal Back 10™ Reg. Height (E2615)	565.00
<input type="checkbox"/> PBT	Personal Back 10 Tall Height (E2615)	565.00
<input type="checkbox"/> PPR	Personal Back 10 Plus Reg. Height (E2615)	645.00
<input type="checkbox"/> PPT	Personal Back 10 Plus Tall Height (E2615)	645.00
<input type="checkbox"/> ITPR	InTouch Propel™ Back (E2611).....	438.00

**FOR ADDITIONAL SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

**INVACARE REHAB SEATING AND POSITIOING PRODUCTS
PRICE LIST FORM # 07-027.**

Basic Chair Includes:

- Solara Tilt-in-Space Frame
- Seat Width Range
- Seat Width (W15 – W19)
- Seat Depth (D15 – D19)
- Arm Type (FXDESK or FXFULL)
- Fixed Back Height Selection
- Fixed or Fold Down Back
- Frame Finish
- Caster Options (5SOFT, AT902 or AT903)
- Wheel Options
- Axle Type Selection
- URETH Urethane (treaded) Tires
- Handrim Options (AHDRM or CHDRM)
- Front Rigging Selection (70NTAPER, SWING1319, SWING914, 90SWING, 70STAPER, or 70MFX)
- 1651 Flip Up Composite Footplates
- TLC Tilt Lock Out Stops (1 pair)
- TLL Tilt Lock Out Levers (1 pair)
- Anti-tippers
- Push or Pull to Lock Wheel Locks

<ul style="list-style-type: none"> • Frame Style • Seat Widths • Seat Depths • Back Height • Adjustable Angle Back • Tilt Range • Camber Angles • Tubing O.D. • Weight Limit • Weight of Chair • Warranty 	<p>Tilt 14" – 22" 15" – 22" 20" and 24" (see Back Height Category) 0°, 10°, 20°, 30° -5° to 50° 0° 1" Reg Seat Rail 1" 250 lbs 28.5 lbs* 41.5 lbs** Lifetime on frames and crossbraces</p> <p><small>* Base chair weight without accessories ** 18" x 18", FXDESK, BH20, FLDBACK, AT902, 12COM, A1577, PNEU, 70NTAPER, 1651, TLC, Anti-tippers, Wheel Locks</small></p>
--	---

Available Through Our Specials Departments

Adjustable Clamp-on Sector Blocks (pair)	SM-031007103339
20" High Reclining Back	SM-050503023551
Telescoping I.V. Tube ²⁹	SM-990419014333
23" – 26" Wide	SM-990604021232
80° Footrest Hanger	SM-200411121756
1" to 4" Higher Arms ³⁰	SM-990609025414

FOOTNOTES

²⁷ For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.

²⁸ Available in even sizes only. See form 07-027 for additional information.

FOOTNOTES

²⁹ Not Available w/TRRO, TRBKTS.

³⁰ Available in 2-pt adjustable arms.



Invacare Corporation
USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications and prices are subject to change without notification.
©2008 Invacare Corporation