



Yes, you can.®

Price List and Order Form

Effective April 1, 2007 Revised July 1, 2008

For ease of ordering, contact customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

POWER CHAIR TYPE

- M61R M61 w/SureStep and Std Pwr Elev. Seat5,926.00
w/MK660ACC Controller & SPJ+ Joystick w/Actuator (K0830)
- M50SEAT Pronto M51 Seat Only.....740.00

SEAT SELECTION

- ASBA Seat w/ Adj. Seat Tilt & Back Angle..... No Charge
- M5SRB Van Seat Semi Recline with Solid Base.....70.00

ASBA SEAT

SEAT WIDTH

- W16 Seat Width 16" No Charge
- W17 Seat Width 17" No Charge
- W18 Seat Width 18" No Charge
- W19 Seat Width 19" No Charge
- W20 Seat Width 20" No Charge

SEAT DEPTH

- D16 Seat Depth 16".....**Standard**
- D17 Seat Depth 17"..... No Charge
- D18 Seat Depth 18"..... No Charge
- D19 Seat Depth 19"..... No Charge
- D20 Seat Depth 20".....529.00

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery..... No Charge

BACK SELECTIONS

SEAT BACK TYPE

- SECT Sectional Rear Uprights No Charge

SEAT BACK HEIGHT

- BH16 Back Height 16" No Charge
- BH17 Back Height 17" No Charge
- BH18 Back Height 18" No Charge
- BH19 Back Height 19" No Charge
- BH20 Back Height 20" No Charge

ARM SELECTIONS

ARM TYPE ¹

- 28F Desk Length Fixed Height Conv.**Standard**
- 29F Desk Length Adj Height Conv. (E0973)230.00
- 18F Full Length Fixed Height Conv. No Charge
- 19F Full Length Adj Height Conv. (E0973)230.00

FOOTNOTES

¹ Arms are flip-up and detachable.

VAN SEAT – WITH SOLID BASE

VAN UPHOLSTERY

- GV Grey Vinyl..... No Charge

VAN WIDTH AND DEPTH

- VW16 16" Wide Van Seat No Charge
- VD16 16" Depth No Charge
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW18 18" Wide Van Seat **Standard**
- VD16 16" Depth **Standard**
- VD17 17" Depth No Charge
- VD18 18" Depth No Charge
- VW20 20" Wide Van Seat329.00
- VD18 18" Depth No Charge
- VD19 19" Depth No Charge
- VD20 20" Depth329.00

VAN ARM TYPE

- VAH Full Length Adjustable Height (E0973)230.00
- VFH Full Length Fixed Height..... **Standard**
- VAHD Desk Length Adjustable Height (E0973).....230.00
- VFHD Desk Length Fixed Height No Charge

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap **Standard**
- 1311BK Airline Buckle Seat Positioning Strap50.00

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960)100.00
- 9515 Push Button Style Chest Pos. Strap (E0960)100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960)100.00

DRIVE WHEEL TYPE

DRIVE WHEEL TYPE

- DW10535 10½ "x 3½" Pneumatic Tires..... **Standard**
- DW10535FF 10½ "x 3½" Flat Free Tires No Charge

Invacare® Pronto™ M61™ with SureStep Accessories and Options

COLORS

SHROUD COLORS

- RD Candy Red Pearl..... No Charge
- BL Deep Blue Pearl..... No Charge

BATTERY OPTIONS

INSTALLATION OF BATTERIES

- UIINSTAL U1 Battery Installation (On Chair) (E2365)232.00

ACCESSORIES

MISCELLANEOUS ACCESSORIES

- SEATO2 Oxygen Holder (ASBA Seat Only) (E2208) 150.00 ²
- 1209 Tire Pump35.00
- VANO2 Oxygen Holder (Solid Base Van Seat) (E2208)... 150.00
- P725 Crutch and Cane Holder (E2207) 55.00 ²
- P770 Walker Holder99.00 ²
- 1831M4 Easy Remote Programmer E.R.P.316.00
- LESTR Lester 8 Amp Battery Charger (E2367).....320.00

FRONT RIGGING OPTIONS

FRONT RIGGING

FOOTBOARD

- 90FB Footboard**Standard**
- L90FB Omit Footboard Less 40.00

FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13½"– 19"165.00
- PH904A Swingaway ELR 16" – 19¾" (E0990).....165.00
- PHAL4A Smartleg® Art. Legrests 17¼" – 21½" (K0053)...225.00

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates..... 115.00
- 1450 Aluminum Footplates Deluxe No Charge
- 1650 Aluminum Footplates No Charge
- 1651 Flip-Up Composite Footplates No Charge
- AT5543 Adjustable Angle Flip Up Footplates 165.00

FOOTREST/LEGREST ACCESSORIES

- LPT2 Longer Pivot Slide Tube For 135083.00
- A1337 Calf Strap (K0038).....27.00
- IMPACTG Impact Guards51.00
- ALPT Longer Pivot And Slide Tube84.00
- ALPT4 Longer Pivot And Slide Tube 4"84.00

FOOTNOTES

² Only one option can be used at one time.

SEATING & POSITIONING

FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
See Price List Form 07-027

Specifications	ASBA	Van Seat w/Solid Base
DIMENSIONS		
Lengths (w./Footboard folded)	34"	34"
Lengths (w./Footboard down)	38" - 40"	38" - 40"
Widths (excluding Joystick)	24"	24"
Heights *	35" - 37"	35" – 37"
Widths	17" - 21" (seat pan to floor)	19" - 23"
Depths	16" – 20"	16", 18", 20"
Back Ht w/o Headrest	16" – 20"	19"
Seat to Floor Heights	17" – 19"	17" – 19"
Weight		
With Batteries	192 lbs.	198 lbs.
Without Batteries	145 lbs.	151 lbs.
Wheel and Caster Sizes		
Front Casters	6" x 2"	6" x 2"
Drive Wheels	10.5" x 3.5"	10.5" x 3.5"
Rear Casters	6" x 2"	6" x 2"
Maximum Speed	4.25 mph	4.25 mph
Incline Capability	9 deg	9 deg
Weight Capacity	300 lbs.	300 lbs.
Turning Radius		
Front (w/footboard down)	19.5"	19.5"
Rear	19.5"	19.5"
Ground Clearance	2.6"	2.6"
Warranties	See Owner's Manual	See Owner's Manual



Invacare Corporation

USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

www.invacare.com

All rights reserved. Trademarks are identified by the symbols™ and®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications and prices are subject to change without notification.
©2007 Invacare Corporation