

Invacare® Storm Series® Torque™ SE Power Wheelchair



Yes, you can.®

Price List and Order Form

Effective October 9, 2006 Revised November 8, 2007

For ease of ordering, contact Customer Service Toll Free at
800-333-6900 Fax: 800-678-4682
To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

CHAIR TYPE

- 3GTQSE Torque SE w/MK6i™ 5,288.00 ¹
- 3GTQSEBASE Torque SE Base Only w/MK6i 5,048.00 ¹
- 3GSEAT Series Seat Only 740.00 ²

*TRANSPORT TIE-DOWN

TRRO includes four factory installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol I Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lbs to 209 lbs.

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

TRRO FEATURES FOR 3G TORQUE SE

Features Required for TRRO	Optional Features for TRRO
<ul style="list-style-type: none"> • Torque SE Base & Seat • 300 lb user weight limit • Short or Long Base Frame • Adj Rehab Seat w/Adj Seat Tilt & Back Angle • Any Seat Width and Depth • Low Seat-to-Floor Height • Sect or Straight Rear Uprights • 18" - 24" ASBA Back Height • 90° or 95° Back Angle • Any Arm Type • 8" x 2.25" Semi Pneu Casters • Suspension Shocks • 14" x 3" Tires w/Foam Filled Insr • Rear Wheel Standard Location • Any Color Frame Finish • Center Pivot Style Bracket • 60° Swingaway, 70° Swingaway or 70° Tapered Footrests • Group 24 Battery Box 	<ul style="list-style-type: none"> • TrueTrack Upgrade • Fender Option • Flip-Up Composite Footplates, Aluminum Footplates-Deluxe, Adjustable Angle Flip Up Footplates • Calf Strap • Impact Guards • All PSR/PSF Joysticks • Left Hand Mounted Joystick • "T" Handle Flexible Joystick Ext • Straight Handle Flexible Joystick Ext • MK6i Remote Programmer • SPJ+ or MPJ™ + Joysticks • Reinforced Back Upholstery • Group 24 Batteries

Please note that all TRRO chairs are set at a 10° seat angle

FOOTNOTES

- ¹ Torque SE comes standard with 300 lb user weight capacity, 6.5 mph max speed. Performance 4-pole motor package. MK6i SPJ+ Joystick, MK690 Controller. Base only comes standard with MK690ACC Controller and SPJ+ w/PSS Joystick. Base only does not come standard with Seat Positioning Strap.
- ² When ordering a replacement seat and changing depth to 20" or deeper, order new arms.

TRANSPORT TIE-DOWN ³

- TRRO Transport Ready Option 300.00
- TRBKTS Wheelchair Transport Brackets 250.00

USER WEIGHT LIMIT

- U300 User Weight <300 lbs No Charge

FRAME TYPE

- SFS Short Base **Standard** ⁴
- LFS Long Base No Charge

ADJUSTABLE SEAT WIDTH RANGE ⁵

- 1216W 12" to 16" Seat Width Range No Charge
- 1620W 16" to 20" Seat Width Range No Charge
- 2022W 20" to 22" Seat Width Range No Charge

FRAME WIDTH ⁶

- W12 12" Seat Width 110.00 ⁵
- W13 13" Seat Width 110.00 ⁵
- W14 14" Seat Width 110.00 ⁵
- W15 15" Seat Width No Charge ⁵
- W16 16" Seat Width No Charge
- W17 17" Seat Width No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width No Charge
- W20 20" Seat Width 329.00
- W21 21" Seat Width 329.00
- W22 22" Seat Width 329.00

VAN SEAT WIDTH ⁴

- VW18 18" Wide No Charge
- VW20 20" Wide 329.00
- VW22 22" Wide 329.00

DRIVE TYPE OPTION

- TRUTK TrueTrack Upgrade 1,899.00 ⁷
- PERF Performance 4-Pole Motor No Charge

FOOTNOTES

- ³ Must select Group 24 battery boxes for TRRO or TRBKTS.
- ⁴ Not available on seat depths of 19" to 22". Manual Reclining Back is not available. Van Seat Depth set at the following: 18" wide set at 16" deep, 20" wide set at 18" deep, 22" wide set at 20" deep.
- ⁵ Jr. seat widths require a straight rear upright and have a user weight limit of 150 lbs. Manual reclining back not available. Additional charge may apply with widths/depths less than 16".
- ⁶ Contoura™ Seating available in even widths from 16" - 22".
- ⁷ TrueTrack driving technology, which enables driver to track straighter regardless of the terrain, includes Monroe® Spring Shock Suspension, <300 lb user weight limit and 7 mph max speed.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® Storm Series® Torque™ SE

Accessories and Options

SEAT SELECTION

- ADJASBA Adj Rehab Seat w/ Adj Seat Tilt & Back Angle.. **Standard**
- VAN Van Style Captain Seat..... 100.00
- ADJRECL Manual Recline for Adj ASBA (E1226)..... 744.00
- CTS Contoura Seating 395.00 ⁸
- OMCU Omit Contoura Seat Cushion..... Less 20.00

ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range..... No Charge ⁹
- 1619D 16" to 19" Seat Depth Range..... No Charge
- 1922D 19" to 22" Seat Depth Range..... No Charge

SEAT DEPTH SETTING ⁸

- D12 Seat Depth 12" (Short Base) 110.00 ¹⁰
- D13 Seat Depth 13" (Short Base) 110.00 ¹⁰
- D14 Seat Depth 14" (Short Base) 110.00 ¹⁰
- D15 Seat Depth 15" (Short Base) No Charge ¹⁰
- D16 Seat Depth 16" (Short Base) No Charge
- D17 Seat Depth 17" (Short Base) No Charge
- D18 Seat Depth 18" (Short Base) No Charge
- D19 Seat Depth 19" (Long Base)..... No Charge ¹¹
- D20 Seat Depth 20" (Long Base)..... 529.00 ¹¹
- D21 Seat Depth 21" (Long Base)..... 529.00 ¹¹
- D22 Seat Depth 22" (Long Base)..... 529.00 ¹¹

ADJASBA SEAT TO FLOOR HEIGHT ¹²

- LOW Seat To Floor – Low..... No Charge
- MED Seat To Floor – Med..... No Charge
- TALL Seat To Floor – Tall..... No Charge

VAN SEAT STYLE

- SEMIR Semi-Recline No Charge
- FULLR Full-Recline (E1226)..... 469.00

BACK SELECTION

ADJASBA BACK TYPE

- SECT Sectional Rear Uprights..... No Charge
- SRUP Straight Rear Uprights..... 59.00
- 9720P No Push Handles (Straight Canes)..... 82.00
- L8661 Omit Rear Uprights Less 20.00

RECLINER BACK HEIGHT ¹³

- BHR18 Back Height 18" No Charge
- BHR20 Back Height 20" No Charge
- BHR22 Back Height 22" No Charge
- BHR24 Back Height 24" No Charge
- BHR26 Back Height 26" No Charge

FOOTNOTES

⁸ Selected seat depth on Contoura assumes 1" of back cushion compression. For example: when choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please see page 5 for seat size offerings. Contoura seating available in depths of 16" to 22".

⁹ Must select 12" - 18" depth range (1218D) when ordering the 12" - 16" seat width range (1216W).

¹⁰ Manual reclining back not available.

¹¹ Must use long base when selecting 19" - 22" depth.

¹² Based on 5° tilt and 16" seat depth on ADJASBA add .25" for additional 2" of seat depth. See page 5 for nominal Seat-To-Floor Heights.

¹³ Used only with manual reclining seat selection.

BACK HEIGHT

ADJASBA

- BH12 Back Height 12" No Charge
- BH13 Back Height 13" No Charge
- BH14 Back Height 14" No Charge
- BH15 Back Height 15" No Charge
- BH16 Back Height 16" No Charge
- BH17 Back Height 17" No Charge
- BH18 Back Height 18" No Charge
- BH19 Back Height 19" No Charge
- BH20 Back Height 20" No Charge
- BH21 Back Height 21" No Charge
- BH22 Back Height 22" No Charge
- BH23 Back Height 23" No Charge
- BH24 Back Height 24" No Charge

CONTOURA

- TSBH25 Back Height 25" No Charge
- TSBH295 Back Height 29.5" No Charge

BACK ANGLE DEGREE

- 85BA 85° No Charge
- 90BA 90° No Charge
- 100BA 100° No Charge
- 105BA 105° No Charge

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
- RBP Reinforced Back Upholstery No Charge

ARM SELECTION

ARM TYPE ¹⁴

- 18F Full Length Fixed Hgt Conv Flip Up No Charge ¹⁵
- 28F Desk Length Fixed Hgt Conv Flip Up **Standard** ¹⁵
- 19F Full Length Adj Hgt Conv Flip Up (E973) 230.00
- 29F Desk Length Adj Hgt Conv Flip Up (E0973) 230.00

FOOTNOTES

¹⁴ All arms are flip-up and detachable.

¹⁵ Not available with seat widths 12" - 15" must use adjustable height.

Invacare® Storm Series® Torque™ SE Accessories and Options

MK6i ELECTRONICS

MK6i CONTROLLERS

- MK690 90 Amp Controller..... **Standard**
- MK6TT TrueTrack Controller..... No Charge

STANDARD DRIVER CONTROL

- SPJM6 SPJ+ Single Drive Joystick **Standard**

BASE ONLY SPJ+ OPTIONS ¹⁶

- SPJAP SPJ+ with Powered Seating Switch **Standard** ¹⁶
- SPJAPM Single Actuator Operation through SPJ+ (E2310) ..1,300.00 ¹⁷

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6 MPJ™+ Multiple Drive Joystick (E2377)..... 750.00
- PSFM6 Personal Joystick, Inductive on Front (E2399) . 1,130.00 ¹⁸
- PSRM6 Personal Joystick, Inductive on Rear (E2399) .. 1,130.00 ¹⁸
 - PBOD Push Button On/Off w/PB Drive Select No Charge ¹⁹
 - PBSS Push Button On/Off w/Speed Select No Charge ^{19,20}
 - OR On/Off Mounted Right..... No Charge ¹⁹
 - OL On/Off Mounted Left..... No Charge ¹⁹
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377)..... 750.00 ²¹

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

- 1500M6 Prop RIM Head Control w/Display (E2327) 2,875.00
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick w/Display (E2373)..... 1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325)..... 1,500.00 ²²
- 5018M6 Wafer Board & Digital Interface Box (E2322) . 1,550.00 ²²
- 5020M6 Mini Tash Jstck & Digital Interface Box (E2321) . 1,750.00 ²²

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326) .. 380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick – Less Std Mount..... 770.00
- 5018 Wafer Board 310.00 ²³
- 5020 Mini Tash Joystick 405.00 ²³
- PACM6 Proportional Attendant Control..... 850.00
- 1552M Digital Attendant Control..... 850.00 ²³

FOOTNOTES

- ¹⁶ Base Only Systems for Multiple Actuator Powered Seating Systems will require MPJ+, PSR+, PSF+, or MK6i Display & Alternative Driver Control. Includes Egg Switch.
- ¹⁷ Allows Actuator Operation through Driver Control.
- ¹⁸ TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer standard.
- ¹⁹ PSRM6 Only.
- ²⁰ Drive Select through mode switch, or mode switch and Joystick Commands.
- ²¹ Display Only, No Driver control. Allows choice from order form.
- ²² Must also use either a Multiple Drive Joystick or MK6i Display.
- ²³ Requires SNPM6.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract Jstck Mnt MPJ+ (E1028)..... 279.00 ²⁴
- QLA Quad Link Retract Jstck Mnt SPJ+ (E1028) 279.00 ²⁴
- QLM6 Quad Link Retract Jstck Mnt MPJ+ Van (E1028) . 230.00 ²⁵
- QL Quad Link Retract Jstck Mnt SPJ+ Van (E1028) 230.00 ²⁵
- ARM250 Stealth Height Adjustable Jstck Mount..... 210.00
- LEFT Left Handed Mounted Jstck..... No Charge
- RIMHW RIM Hardware for 1500 (E1028) 219.00
- PKG32669 Complete Bib Assembly (Tash/1558) 265.00
- GATMPJ6 Gatlin Midline for MPJ+ (E1028) 600.00
- GAT1812 Gatlin Midline for 1812 (E1028)..... 600.00
- RCM Rear Cane Mount for Compact Joystick 395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card 475.00
- 1813SD Pro Memory Card (USB Ready) 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1826 Chin Cup (E2324)..... 48.00
- A24VPS 24V MK6 Auxiliary Power Source 200.00
- AUX12M6 Auxiliary Module for 1 & 2..... 795.00
- AUX34M6 Auxiliary Module for 3 & 4..... 795.00
- MPY Monoport "Y" Cable for MPJ+ or Display 15.00 ²⁶
- 2PB Dual Push Button..... 147.00 ²⁷
- 2SW Dual Toggle 172.00 ²⁷

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black 65.00
- ASL304 Light Wobble Switch w/611 Mount..... 320.00
- ASL202J Fiber Optic Switch..... 610.00 ²⁸
- ASL208J Adjustable Proximity Switch..... 550.00 ²⁸
- ASL504JM6 Remote Emergency Stop Switch 535.00 ²⁸

FOOTNOTES

- ²⁴ QLAM6 and QLA include ARM250 standard.
- ²⁵ Only available with Van Seating.
- ²⁶ For use with expandable controller systems to provide additional switch sites.
- ²⁷ Option to replace Monoport "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- ²⁸ Must add MK6 Auxiliary Power Source.

Invacare® Storm Series® Torque™ SE

Accessories and Options

FRAME AND UPHOLSTERY COLOR

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue..... No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal..... No Charge
- 121P Lolly Pop Blue..... No Charge
- 122P Cosmic Blue No Charge
- 125P Emerald Green..... No Charge
- 127P Grasshopper Green No Charge

UPHOLSTERY TYPE

- U240 Black Nylon Upholstery..... **Standard**
- U67 Black Vinyl..... No Charge
- GC Grey Cloth..... No Charge ²⁹
- GV Grey Vinyl..... No Charge ²⁹
- TV Tan Vinyl..... No Charge ²⁹

CASTER AND WHEEL SELECTION

CASTER OPTIONS

- 1259 6" x 2" Semi Pneumatic Casters..... No Charge
- 1228-3 8" x 2.25" Semi Pneumatic Casters Lt Grey... **Standard**
- 1251-3 8" x 2" Pneumatic Casters No Charge
- 1220-3 9" x 2.75" Pneumatic Casters No Charge
- 1465-3 9" x 2.75" Foam Filled Casters No Charge

CASTER SUSPENSION ³⁰

- 9213-6 Shock Fork For 6" Casters(E1016) 119.00 ³¹
- 9213-8 Shock Fork For 8" Casters (E1016) 119.00 ³¹

WHEEL OPTION

- 1430-3 14" x 3" Pneumatic Tires..... No Charge
- 1431-3 14" x 3" Tires w/Foam Filled Inserts 250.00
- 1440 14" x 4" Pneumatic Tires..... 200.00 ^{32,33}
- 1441 14" x 4" Tires w/Foam Filled Inserts 320.00 ^{32,33}

FREE WHEEL HUB ³⁴

- FHUB Free Wheel Hub..... 125.00

FOOTNOTES

- ²⁹ Used only with Van Style Seat Selection.
- ³⁰ Not available with TRRO or TRBKTS.
- ³¹ Not available with short base. Not available with 9" casters.
- ³² 14" x 4" wheels increase overall width of the base by 2.5". Not available with TrueTrack (TRUTK).
- ³³ Does not include Storm Series designed mag wheel or fender.
- ³⁴ Free wheel hub increases overall width .75". Not needed or available with TRUTK TrueTrack package.

REAR WHEEL LOCATION ³⁵

- RWF2 Rear Wheel Forward 2" No Charge
- RWF1 Rear Wheel Forward 1" No Charge
- RWS Rear Wheel Standard **Standard**

FENDER OPTION ³⁶

- FDR Fender 320.00

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap..... 50.00 ³⁷
- 1515 Push Button Style Positioning Strap **Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Positioning Strap (E0960)..... 100.00
- 9515 Push Button Style Chest Positioning Strap (E0960).. 100.00
- 7321BK Hook & Loop Chest Positioning Strap (E0960)... 100.00

MISC ACCESSORIES

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks 45.00
- VANO2 O2 Holder – Van (E2208) 150.00
- SEATO2 Oxygen Holder (ASBA Seat Only) 150.00
- 9205-3 Suspension Shocks (E1016) 400.00 ³⁸
- 9205H-3 Suspension Shocks Heavy Duty 400.00 ³⁸

BATTERY SELECTION

BATTERY BOXES

- 22NF 22NF Battery Box..... No Charge
- GP24 Group 24 Battery Boxes No Charge

BATTERIES ³⁹

- 22NFBATTERY 22NF Battery (On Chair) (E2361)..... 300.00
- 24BATTERY 24Gel Battery (On Chair) (E2363) 410.00

HANGER BRACKET SELECTION

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge
- NHB No Hanger Bracket No Charge

FOOTNOTES

- ³⁵ This location changes the center-of-gravity and performance. Forward wheel position will make turning and climbing easier, but more prone to "wheely".
- ³⁶ Standard with TrueTrack (TRUTK) option.
- ³⁷ 1311BK not available in widths 12" - 15" wide.
- ³⁸ Not available with TRUTK package selected.
- ³⁹ Batteries will be installed by Invacare Corporation at no additional charge. Chair will be shipped complete w/batteries.

Invacare® Storm Series® Torque™ SE Accessories and Options

FRONT RIGGING OPTIONS

FRONT RIGGING

CENTER PIVOT

- 60HD 60° Swingaway Footrests No Charge
- 70NHD 70° Swingaway Footrests **Standard** ⁴⁰
- 70STAPER 70° Tapered Footrests..... No Charge
- AT5544 Elevating Legrests (E0990)..... 255.00
- 90SWING 90° Swingaway Footrests No Charge

PIN STYLE

- PAS4A Smartleg® Art. Legrests 17.25" – 21.5" (K0053). 320.00
- PW93 Swingaway Footrests 13.5" – 19" No Charge
- P904A Swingaway Legrests 16" – 19.75" (E0990)..... 255.00
- PAL4A Smartleg Art. Legrests 17.25" – 21.5" (K0053) .. 320.00 ⁴¹

POWER LEGS

- ELRPW Power Elev Legrest w/Articulation (E1010)..... 1,895.00 ^{42,43}

POWER LEGS CONTROL MODULE ⁴⁴

- 4WSB Four Way Switch Box No Charge ⁴⁵
- S4WSB Multiple Actuator Interface Box (E2311) 2,600.00 ⁴⁶
 - FWT 4 Way Toggle..... No Charge
 - QPB 4 Push Buttons No Charge

For all systems with two or more actuators you must choose 4WSB at No Charge or the S4WSB with an Up-Charge.

S4WSB Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

FWT or QPB is required when choosing the 4WSB or S4WSB.

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- 1450 Aluminum Footplates-Deluxe No Charge
- AT5043 Flip-Up Footplates No Charge
- S1263 One Piece Wooden Footboard..... No Charge
- AT5543 Adjustable Angle Flip Up Footplates..... 165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) 27.00
- 1600BK Heel Loops with Ankle Straps (pair) (E0951) 56.00
- IMPACTG Impact Guards..... 51.00
- A9342 "H" Heel Strap (K0039)..... 51.00
- LPT2 Longer Pivot Slide Tube For 1350 83.00
- ALPT Longer Pivot And Slide Tube 84.00
- ALPT3 3" Extension For 90° Footrest 31.00
- ALPT4 Longer Pivot And Slide Tube 4" 84.00

SEATING & POSITIONING
FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM 07-027

FOOTNOTES

- ⁴⁰ Standard with impact guards.
- ⁴¹ For seat widths greater than or equal to 16".
- ⁴² Independent powered legrests include elevating and articulating legrests with angle adjustable footplates. Extension range of 13" - 20".
- ⁴³ To operate power legs must choose: 4WSB - 4 Way Switch Box or S4WSB Multiple Actuator Interface Box. Must include with pick either FWT- 4-Way Toggle or QPB - 4 Push Buttons.
- ⁴⁴ These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories.
- ⁴⁵ Required with all multiple actuator systems. Provides a D-9 port for the 4-way toggle or 4 push buttons.
- ⁴⁶ Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control. Not compatible with SPJ+ joysticks.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® Storm Series® Torque™ SE Accessories and Options

Chair Features and Options
<ul style="list-style-type: none"> • Independent rear wheel suspension with new rubber element design. • 6.5 mph standard speed range. • Programmable performance to meet individuals needs. • 17.5" Seat-to-Floor Standard. • 19.75" and 21.5" Optional. • Seat location system allows for balanced weight distribution. • Continuous tilt seat 0° - 10°. • Dynamic anti-tippers. • Rigid Seat Pan. • Convenient rear battery access. • 22NF batteries (not supplied) in a single 22 NF box Single battery connector. • Arms are flip back/removable. • Adjustable depth rigging attachment 2" (1" increments). • Automatic dual mode battery charger (Lester Brand). • Compact size: 24.75" overall base width w/14" x 3" wheels. • Weight limit – 300 lbs.

Torque SE Specifications		
* Overall Width	24.75" 27.25"	w/14" x 3" wheel w/14" x 4" wheel
**Overall Length	29.5" 32.5"	Short Base Long Base
***Seat-to-Floor	17.5" Standard 19.75" Optional	Low Medium
Battery	22NF deep cycle; gel/sealed (Not supplied)	
Weight w/o batteries 149 lbs (40 – 45 lbs each battery)		
Shipping Weight 266 lb. (18" x 18" seat w/Group 24 batteries)		
*Measurement based on wheel base. Seat width is width picked plus 7.5" for the armrest and joystick.		
**Measured from headtube to anti-tipper.		
***Based on 5° tilt and 16" seat depth; add .25" for additional 2" of seat depth		
Note: All specifications and dimensions are approximate.		

Weight Limit 300 Pounds (Standard)
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

NOTE: All weight limits are inclusive of add-on equipment and seat systems.

Footrest Extension Range			
	AT5043	AT5543	1651
60HD	N/A	12.25" – 16.5"	13.5" – 17.5"
70MFX	6.75" – 12.5"	4.5" – 10"	N/A
70NHD	N/A	10.75" – 15.75"	11.75" – 17"
70STAPER	N/A	12.5" – 16.5"	13.25" – 17.25"
90SWING	4" – 9"	5" – 8"	N/A
Measured from Seat Pan to Back of Footplate			

MEMORY CARD - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

MEMORY CARD - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Manual Recline Specifications			
Seat Widths	14" - 22" (1" increments)		
Seat Depths	16" - 22" (1" increments) (Adjustable 16" - 19", Fixed 20" - 22")		
Back Heights	24" Standard, 18.5", 20", 22", 24", 26"		
Weight Limits	300 lbs. with Standard Cylinders		
Headrest	10" Height		
Recline Range	90° to 170° (Approximate)		
Overall Back Height	Low Seat	Med. Seat	High Seat
{w/o extension}			
Standard Base	41.5"	43.5"	45.25"
{Measured with 5° Tilt and 10" for Extension}			

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights	25" or 29.5"						



Invacare Corporation www.invacare.com

USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols™ and®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Therafin is a registered trademark of W/R Medical Electronics Co. Otto Bock is a trademark of Otto Bock Orthopaedische Industrien GMBH & Co. Freedom Designs is a registered trademark of Freedom Designs, Inc. Specifications and prices are subject to change without notification.
©2007 Invacare Corporation