

Invacare® PTO (Power Tilt Only) for Torque™ SE Series



Price List and Order Form Effective November 15, 2006 Revised April 9, 2007

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Yes, you can.®

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME

3G TORQUE SE

3GTQSE-PTO Torque™ SE for PTO w/MK6i™ (K0877).....5,048.00

POWERED SEATING

PTO-STM Power Tilt Only (E1002).....4,285.00

USER WEIGHT LIMITS

STDWT User Weight <250 lbs.....No Charge

SEAT SELECTION

ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle..Standard
 CTS Contoura™ Seating.....395.00
 OMCU Omit Contoura Seat Cushion.....Less 20.00

WIDTH AND DEPTH

ADJUSTABLE SEAT WIDTH RANGE

1216W 12" to 16" Seat Width Range.....No Charge
 1620W 16" to 20" Seat Width Range.....No Charge
 2022W 20" to 22" Seat Width Range.....No Charge

ADJASBA SEAT WIDTH SETTING

^{3,4} TRR12 Seat Width 12".....No Charge
^{3,4} TD13 Seat Width 13".....No Charge
^{3,4} TD14 Seat Width 14".....No Charge
^{3,4} TR15 Seat Width 15".....No Charge
 TR16 Seat Width 16".....No Charge
 TR17 Seat Width 17".....No Charge
 TR18 Seat Width 18".....No Charge
 TR19 Seat Width 19".....No Charge
 TR20 Seat Width 20".....No Charge
 TD21 Seat Width 21".....No Charge
 TR22 Seat Width 22".....No Charge

ADJUSTABLE SEAT DEPTH RANGE

1216D 12" to 16" Seat Depth Range.....No Charge
 1619D 16" to 19" Seat Depth Range.....No Charge
 1922D 19" to 22" Seat Depth Range.....No Charge

ADJASBA SEAT DEPTH SETTING

^{3,4} TAR12 Seat Depth 12" (Short Base).....No Charge
^{3,4} TAR13 Seat Depth 13" (Short Base).....No Charge
^{3,4} TAR14 Seat Depth 14" (Short Base).....No Charge
^{3,4} TAR15 Seat Depth 15" (Short Base).....No Charge
 TAR16 Seat Depth 16" (Short Base).....No Charge
 TAR17 Seat Depth 17" (Short Base).....No Charge
 TAR18 Seat Depth 18" (Short Base).....No Charge
 TAR19 Seat Depth 19" (Short Base).....No Charge
 TAR20 Seat Depth 20" (Long Base).....No Charge
 TAR21 Seat Depth 21" (Long Base).....No Charge
 TAR22 Seat Depth 22" (Long Base).....No Charge

ARMS

5 ARM TYPE

18F Full Length Fixed Hgt. Conv. Flip Up.....No Charge
 19F Full Length Adj. Hgt. Conv. Flip Up.....No Charge
 28F Desk Length Fixed Hgt. Conv. Flip Up.....No Charge
 29F Desk Length Adj. Hgt. Conv. Flip Up.....No Charge

BACK

TILT BACK OPTION (PTO)

CURV Curved Back.....No Charge
 POST Post Only.....No Charge

ADJASBA BACK HEIGHT

TSBH12 Back Height 12".....No Charge
 TSBH13 Back Height 13".....No Charge
 TSBH14 Back Height 14".....No Charge
 TSBH15 Back Height 15".....No Charge
 TSBH16 Back Height 16".....No Charge
 TSBH17 Back Height 17".....No Charge
 TSBH18 Back Height 18".....No Charge
 TSBH19 Back Height 19".....No Charge
 TSBH20 Back Height 20".....No Charge
⁸ TSBH21 Back Height 21".....No Charge
⁸ TSBH22 Back Height 22".....No Charge
⁸ TSBH23 Back Height 23".....No Charge
⁸ TSBH24 Back Height 24".....No Charge

CONTOURA BACK HEIGHT

TSBH25 Back Height 25".....No Charge
 TSBH295 Back Height 29.5".....No Charge

DEGREE OF TILT (PTO)

TR045 0° - 45° of Tilt.....No Charge
 TR550 5° - 50° of Tilt.....No Charge

FRAME TYPE

SFS Short Base.....No Charge
 LFS Long Base.....No Charge

FOOTNOTES

- Power Tilt Only (PTO) is a weight shifting tilt system up to 45° with a weight limit of up to 250 lbs.
- Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please refer to page 5 for seat size offerings.
- Not available on curved back w/headrest tilt back option, must use Post on Junior sizes.
- Additional costs may apply with ≤16" width/depth.
- All arms are flip-up and removable.
- Fixed height arms not available with 12"-15" seat widths or depths.
- Not available on 12" -15" Junior seat widths .
- Not available on post only.

NOTE: All specifications and dimensions are approximate.

Invacare® PTO (Power Tilt Only) for Torque™ SE Series

Accessories and Options

BACK (CONTINUED)

HEADREST TYPE (OPTIONAL)

- TAR127700 Curved Headrest (E0955).....210.00
- TAR127700HW Removable & Adjustable Hardware (E1028)....210.00
- TAR128000 Two-Step Headrest (E0955).....210.00
- TAR128000HW Removable & Adjustable Hardware (E1028)....210.00

TARSYS LATERALS

- TCSLM Clamp Swingaway Laterals Med Pair (E0956)....210.00
- TCSLMHW Swingaway Lateral Hardware (E1028).....210.00
- TCSLL Clamp Swingaway Laterals Large Pair (E0956)..210.00
- TCSLLHW Swingaway Lateral Hardware (E1028).....210.00

DRIVE TYPE OPTION

TRUETRACK UPGRADE

- TRUTK TrueTrack Upgrade.....1,899.00

WHEELS

CASTER OPTIONS

- 1259 6" X 2" Semi-Pneumatic CastersNo Charge
- 1228-3 8" X 2 1/4" Semi Pneu. Casters Lt Grey**Standard**
- 1251-3 8" X 2" Pneumatic CastersNo Charge
- 1220-3 9" X 2 3/4" Pneumatic CastersNo Charge
- 1465-3 9" X 2 3/4" Foam Filled CastersNo Charge

CASTER SUSPENSION

- 9213-8 Shock Fork For 8" Casters (E1016).....119.00
- 9213-6 Shock Fork For 6" Casters (E1016).....119.00

WHEEL OPTIONS

- 1430-3 14" X 3" Pneumatic Tires.....**Standard**
- 1431-3 14" X 3" Tire w/Foam Filled Inserts250.00
- 1440 14" X 4" Pneumatic Tires.....200.00
- 1441 14" X 4" Tires w/Foam Filled Inserts320.00

FREE WHEEL HUB

- FHUB Free Wheel Hub.....125.00

FENDER OPTION

- FDR Fender.....320.00

FRAME COLOR

FRAME FINISH CONTEMPORARY

- 24P Wet Black.....No Charge
- 60P Silver Metallic.....No Charge
- 71P Silver Vein.....No Charge
- 105P Midnight Blue.....No Charge
- 115P Black Prism.....No Charge

TOO HOT

- 30P Sunny Yellow.....No Charge
- 61P Electric Red.....No Charge
- 119P Bubblegum Pink.....No Charge
- 120P Grape MadnessNo Charge

WAY COOL

- 62P Electric Blue.....No Charge
- 104P Electric TealNo Charge
- 121P Lolly Pop Blue.....No Charge
- 122P Cosmic Blue.....No Charge
- 125P Emerald Green.....No Charge
- 127P Grasshopper Green.....No Charge

ELECTRONICS

CONTROLLERS

- MK690ACC 90 Amp Controllerw/Actuator Control.....**Standard**
- MK6TT True Track Controller.....No Charge

STANDARD DRIVER CONTROLS

- SPJAP SPJ™+ w/Power Seating Switch.....**Standard**

ACTUATOR CONTROLS THROUGH SPJ+

- SPJAPM SPJ+ w/Actuator Op thru Joystick (E2310).....1,300.00

EXPANDABLE CONTROLLER SYSTEMS

- MPJM6U MPJ™+ Multiple Drive Joystick (E2377).....750.00
- PSFM6 Personal Joystick, Inductive on Front (E2399).....1,130.00
- PSRM6 Personal Joystick, Inductive on Rear (E2399).....1,130.00
 - PBOD Push Button On/Off w/PB Drive Select.....No Charge
 - PBSS Push Button On/Off w/Speed Select.....No Charge
 - OR On/Off Mounted RightNo Charge
 - OL On/Off Mounted Left.....No Charge
- ASLRDYM6 MK6 Display Only
(No Driver Control) (E2377).....750.00

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

- 1500M6 Proportional RIM Head Cntrl w/ Display (E2327)2,875.00
 - SMHD Small Headrest.....No Charge
 - LGHD Large Headrest.....No Charge
- 1558MM6 Compact Joystick w/Display (E2373).....1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325).....1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322).....1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321).....1,750.00

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326).....380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount.....770.00
- 5018 Wafer Board.....310.00
- 5020 Mini Tash Joystick.....405.00
- PACM6 Proportional Attendant Control.....850.00
- 1552M Digital Attendant Control850.00

FOOTNOTES

1. Available back heights of 17" or greater.
2. TrueTrack driving technology, which enables driver to track straighter regardless of the terrain, includes Monroe® Spring Shock Suspension, <300 lb. User Weight Limit and 7 mph top speed.
3. Must pick long base.
4. 14" X 4" wheels increase overall width of the base by 2 1/2". No fender available. Does not include Storm Series mag wheel. Not available with TrueTrack Package.
5. Wheel will add 3/4" to overall width. Not available with TrueTrack Package. Does not include Storm Series wheel or fender.
6. Allows Actuator Operation through the driver Control and has a port for an attendant actuator switch..
7. Includes Monoport "Y" Cable for additional switch input port.
8. TOD-SS: Toggle On/Off Drive Select, Speed Potentiometer Standard.
9. Drive Select through mode switch, or mode switch and Joystick Commands.
10. PSRM6 Only.
11. Display Only, No Driver Control. Allows choice from order form.
12. Must also use either a Multiple Drive Joystick or MK6i Display.
13. Requires SNPM6 Sip-N-Puff Interface Box.

NOTE: All specifications and dimensions are approximate.

Invacare® PTO (Power Tilt Only) for Torque™ SE Series

Accessories and Options

ELECTRONICS (CONTINUE)

JOYSTICK MOUNTING HARDWARE

<input type="checkbox"/>	QLAM6	Quad Link Retract. Joystick Mnt MPJ+ (E1028).....	279.00
<input type="checkbox"/>	QLA	Quad Link Retract. Joystick Mnt SPJ+ (E1028).....	279.00
<input type="checkbox"/>	ARM250	Stealth Height Adjustable Joystick Mount (E1028) ..	210.00
<input type="checkbox"/>	LEFT	Left Handed Mounted Joystick.....	No Charge
<input type="checkbox"/>	RIMHW	RIM Hardware For 1500 (E1028)	219.00
<input type="checkbox"/>	PKG32669	Complete Bib Assembly (Tash/1558)	265.00
<input type="checkbox"/>	GATMPJ6	Gatlin Midline for MPJ (E1028)	600.00
<input type="checkbox"/>	GAT1812	Gatlin Midline for 1812 (E1028)	600.00
<input type="checkbox"/>	RCM	Rear Cane Mount for Compact Joystick.....	395.00

ELECTRONIC ACCESSORIES

<input type="checkbox"/>	I813M6	MK6i Programmer w/Pro Memory Card.....	475.00
<input type="checkbox"/>	I813SD	Pro Memory Card (USB Ready)	90.00
<input type="checkbox"/>	I560	"T" Handle Flexible Joystick Ext. (E2323).....	77.00
<input type="checkbox"/>	I561	Straight Handle Flexible Joystick Ext. (E2323).....	77.00
<input type="checkbox"/>	I826	Chin Cup (E2324)	48.00
<input type="checkbox"/>	A24VPS	24V MK6 Auxilliary Power Source	200.00
<input type="checkbox"/>	AUX12M6	Auxillary Module For 1 & 2	795.00
<input type="checkbox"/>	AUX34M6	Auxillary Module For 3 & 4	795.00

ASL OPTIONAL ACCESSORIES

<input type="checkbox"/>	EGSBLK	Egg Switch- Black	65.00
<input type="checkbox"/>	ASL304	Light Wobble Switch w/611 Mount.....	320.00
<input type="checkbox"/>	ASL202J	Fiber Optic Switch.....	610.00
<input type="checkbox"/>	ASL208J	Adjustable Proximity Switch	550.00
<input type="checkbox"/>	ASL504M6	24V Remote Emergency Stop Switch	535.00
<input type="checkbox"/>	ASL300A	Egg Switch w/611 Mount	300.00

POWER SEATING ACCESSORIES

<input type="checkbox"/>	SANODE	One Actuator Op. thru Driver Control (E2310) ..	1,300.00
--------------------------	--------	---	----------

POSITIONING STRAPS

SEAT POSITIONING STRAPS

<input type="checkbox"/>	I311BK	Airline Buckle Seat Positioning Strap.....	50.00
<input type="checkbox"/>	I515	Auto Style Buckle Seat Positioning Strap ...	Standard

CHEST POSITIONING STRAPS

<input type="checkbox"/>	I731BK	Airline Buckle Chest Positioning Strap (E0960).....	100.00
<input type="checkbox"/>	I9515	Auto Style Chest Positioning Strap (E0960)	100.00
<input type="checkbox"/>	I7321BK	Hook & Loop Chest Positioning Strap (E0960).....	100.00

ACCESSORIES

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	I303	Wheel Locks	45.00
<input type="checkbox"/>	I209	Tire Pump	35.00

BATTERIES

BATTERIES

<input type="checkbox"/>	22NFBATTERY	MK Battery Gel-2 Each (on chair) (E2361)	300.00
<input type="checkbox"/>	24BATTERY	24 Gel Battery (on chair) (E2363).....	410.00

FRONT RIGGING OPTIONS

FRONT RIGGINGS

PIN & SADDLE

<input type="checkbox"/>	PAS4A	Smartleg® Articulating Lgrsts 17 1/4"-21 1/2" (K0053) ...	320.00
<input type="checkbox"/>	PVW93	Swingaway Footrest. 13 1/2"-19"	No Charge
<input type="checkbox"/>	P904A	Swingaway Elevating Lgrsts 16"-19 3/4" (E0990) ...	303.00
<input type="checkbox"/>	PVW93ST	Footrests 14 3/4"-19 3/4"	165.00
<input type="checkbox"/>	PAL4A	Smartleg® Articulating Lgrsts 17 1/4"-21 1/2" (K0053) ...	320.00

CENTER PIVOT

<input type="checkbox"/>	70NHD	70° Swingaway Footrests	Standard
<input type="checkbox"/>	70STAPER	70° Tapered Footrests	No Charge
<input type="checkbox"/>	AT5544	Manual Elevating Legrests (E0990)	255.00

FOOTPLATE OPTIONS

<input type="checkbox"/>	NN163	Standard Footplate.....	No Charge
<input type="checkbox"/>	I651	Flip-Up Composite Footplates	No Charge
<input type="checkbox"/>	AT5543	Adjustable Angle Flip-Up Footplates	165.00

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/>	I820	Non-Removable Device	12.00
<input type="checkbox"/>	I337	Calf Strap (E0038)	27.00
<input type="checkbox"/>	I600BK	Heel Loops With Ankle Straps (E0951)	56.00
<input type="checkbox"/>	I345	Toe Loops (E0952)	39.00
<input type="checkbox"/>	I346	Heel Loops (E0951)	39.00
<input type="checkbox"/>	I9342	"H" Strap (K0039)	56.00
<input type="checkbox"/>	IMPACTG	Impact Guards.....	51.00
<input type="checkbox"/>	LPT2	Longer Pivot Slide Tube For I350	83.00
<input type="checkbox"/>	ALPT	Longer Pivot And Slide Tube.....	84.00

FOOTNOTES

1. QLAM6 and QLA include ARM250 standard.
2. Not available with SPJ+ Joystick. Must add MK6 Auxilliary Power Source.
3. Required for Tilt Operation through Driver Control - Multiple Drive Systems only. For Single Drive Systems choose SPJAPM. The Egg Switch with the 611 Mount does not come standard when ordering a SANODE.
4. Not available on 12" - 15" wide.
5. Not available with pneumatic tires.
6. Footrests are standard with heel loops.
7. Only available on 16" or wider.
8. 70STAPER not available with 10" casters and flip up composite footplates on 16", 17", and 20" seat depths, also not available with 9" casters and flip-up composite footplates on 19" and 20" widths, impact guards are standard.
9. Standard on 70NHD & AT5544.
10. Not available with 70STAPER.

NOTE: All specifications and dimensions are approximate.

SEATING & POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING
REHAB SEATING & POSITIONING PRODUCTS
SEE PRICE LIST FORM # 07-027.**

Invacare® PTO (Power Tilt Only) for Torque™ SE Series

Accessories and Options

Specifications for PTO System	
Seat Height (seat pan to floor) - Standard Frame (16" depth)	0° 19 1/2" 5° 20 1/2"
Range Tilt	45°
Overall length	
16" deep	41 1/2"
18" deep	42 1/2"
20" deep	43 1/2"
Switch Option Interface	Optional
Back Style	Back Pan w/ 1" vicous foam Back Post only
Back Height	12" - 24"
Fixed Adjustable Back	Standard
Armrest Style	Adjustable Height Desk or Full Length
<p>NOTE: Seat height measurement from front edge of seat pan to floor with 5° of seat angle. Seat height will vary slightly depending upon seat depth and Storm frame option.</p> <p>NOTE: All specifications are subject to change without notice.</p>	

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights	25" or 29.5"						

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Updates, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Back Height Specification Chart								
Back Heights	18"	19"	20"	21"	22"	23"	24"	25"
Tilt Only Systems								
- Back Cane Height	18"	19"	20"	20"	20"	20"	20"	N/A
- ABS Back Pan/Cushion Height	14 1/4"	15 1/4"	16 1/4"	17 1/4"	18 1/4"	19 1/4"	20 1/4"	N/A

Torque SE Specifications		
Drive Option	Rear Wheel	
Base Options	Short Standard Long Standard	
I Seat Option	Adjustable Back Angle	Standard
Seat Widths	14" - 22" 1" increments	Standard Base
Seat Depths	14" - 18" 1" increments	Standard Base
	19" - 22" 1" increments	Long Base
Wheels	14" x 3"	Standard Base
Casters	8" Semi-Pneumatic	Standard Base
	9" Pneumatic	Long Base
Overall Width	25"	Standard Base
	27 1/2"	14" x 4" wheels
Overall Length	41 1/2" (w/70° swingaway)	Short Base
	44 1/2" (w/70° swingaway)	Long Base
Battery	Group 24 Deep Cycle; (Not Supplied) or 22NF	
Turning Radius	29 1/2" to 33 1/2" RWD	Short Base
	32 1/2" to 35 1/2" RWD	Long Base
	Front 21" MWD	
Electronics	MK6i	
Weight Limits	250 lb PTO Junior Sizes (12"-15"w) 150 lb.	
Shipping Weight	230 lb. with batteries	
	190 lb. w/o batteries	
	* Based on 5° Tilt and 16" Seat Depth, Add 1/4" for additional 2" of Seat Depth	
<p>NOTE: All specifications and dimensions are approximate</p>		

Tech Tips
<p>Operating Powered Seating through the Driver's Control All chairs require adding or replacing one component to allow actuator operation through the driver control.</p> <p>SPJ+ with PSS: Replace joystick with SPJ+ with ACC (Actuator Control) -E2310</p> <p>Single Actuator Systems with MPJ+, PSR+, PSF+ MK6i Display: Add a SANODE (Single Actuator Node) to the system - E2310</p> <p>Multiple Actuator Systems: Replace the standard 4 Way Switch box with a Multiple Actuator Interface Box - E2311</p>

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Therafin is a registered trademark of W/R Medical Electronics Co.
Monroe is a registered trademark of Tenneco Automotive, Inc.
Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc.
Specifications and prices are subject to change without notification.
©2007 Invacare Corporation