

Invacare® TDX® SI Heavy Duty Power Wheelchair

Group 2 Heavy Duty; 450 lb. Weight Capacity



Price List and Order Form

Effective December 15, 2008 rev March 17, 2009

For ease of ordering, contact customer Service Toll Free at:
 800-333-6900 Fax: 800-678-4682
 To order Literature: 800-828-6282

****Recent updates are in RED**

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

STANDARD CONFIGURATIONS

For the best value, order the TDX-SI HD in Standard Configurations.
 Each Standard Configuration is pre-configured with the most popular options – you just choose the seat size and the frame color.

Group 2 HD - K0825 - Van Seat

TDXSIV-HD-S..... \$6,500

Standard Features:

- Semi-Recline Van Seat, Gray Vinyl
- 22NF Gel Batteries Installed on Chair (E2365)
- Full Length, Adjustable Height Arms (E0973)
- 14" Flat Free Tires
- Foot Platform
- Push-button Style Seat Positioning Strap
- 4-Pole, 650 Watt Motors
- MK690 90 Amp Controller
- 8 Amp Off-Board Charger

FRAME COLOR

- 61P Electric Red No Charge
 62P Electric Blue No Charge

VAN SEAT SIZE

- VW20 20" Wide Van Seat No Charge
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth No Charge
- VW22 22" Wide Van Seat No Charge
 VD18 18" Seat Depth No Charge
 VD19 19" Seat Depth No Charge
 VD20 20" Seat Depth No Charge

FOOTREST / FRONT RIGGINGS

- 90FB Foot Platform.....Standard
 70NHD Swingaway Footrests165.00
 AT5544 Manual Elevating Legrest (E0990).....255.00

Group 2 HD - K0824 - Van Seat w/ Solid Pan

TDXSI-HD-S..... \$6,500

Standard Features:

- Semi-Recline Van Seat w/ Solid Seat Pan, Gray Vinyl
- 22NF Gel Batteries Installed on Chair (E2365)
- Full Length, Adjustable Height Arms (E0973)
- 14" Flat Free Tires
- Foot Platform
- Push-button Style Seat Positioning Strap
- 4-Pole, 650 Watt Motors
- MK690 90 Amp Controller
- 8 Amp Off-Board Charger

FRAME COLOR

- 61P Electric Red.....No Charge
 62P Electric BlueNo Charge

VAN SEAT W/ SOLID SEAT PAN SIZE

- VW20 20" Wide Van Seat w/ Solid Seat Pan.....No Charge
 VD18 18" Seat DepthNo Charge
 VD19 19" Seat DepthNo Charge
 VD20 20" Seat DepthNo Charge
- VW22 22" Wide Van Seat w/ Solid Seat Pan.....No Charge
 VD18 18" Seat DepthNo Charge
 VD19 19" Seat DepthNo Charge
 VD20 20" Seat DepthNo Charge

FOOTREST / FRONT RIGGINGS

- 90FB Foot Platform.....Standard
 70NHD Swingaway Footrests165.00
 AT5544 Manual Elevating Legrest (E0990)255.00

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information go to www.invacare.com.

CUSTOM CONFIGURATIONS

For clients with more advanced seating and positioning needs, Custom Configurations allow you to customize the TDX-SI HD to meet the specific needs of each client.

CUSTOM CONFIG - FRAME SELECTION

POWER CHAIR TYPE

- TDXSIV-HD TDX-SI HD w/ Van Seat (K0825)5,970.00
- TDXSISEAT-HD TDX-SI HD Seat Only.....755.00

FRAME COLOR

- 61P Electric Red No Charge
- 62P Electric Blue No Charge

USER WEIGHT LIMIT

- U450 User weight < 450 lb.Standard

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets.....250.00

*** Transport Tie Down**

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

VAN SEAT SELECTIONS

SEAT SELECTION

- SEMIR Semi-Recline Van SeatStandard

UPHOLSTERY COLORS

- GV Gray Vinyl UpholsteryStandard
- GC Gray Cloth Upholstery No Charge
- TV Tan Vinyl Upholstery No Charge

VAN SEAT SIZE

- VW18 18" Wide Van Seat..... No Charge
 - VD16 16" Seat Depth No Charge
 - VD17 17" Seat Depth No Charge
 - VD18 18" Seat Depth No Charge
- VW20 20" Wide Van SeatStandard
 - VD18 18" Seat Depth No Charge
 - VD19 19" Seat Depth No Charge
 - VD20 20" Seat DepthStandard
- VW22 22" Wide Van Seat..... No Charge
 - VD18 18" Seat Depth No Charge
 - VD19 19" Seat Depth No Charge
 - VD20 20" Seat Depth No Charge
- VW24 24" Wide Van Seat.....469.00
 - VD20 20" Seat Depth No Charge
 - VD21 21" Seat Depth No Charge
 - VD22 22" Seat Depth No Charge

ARM TYPE (PAIR) ¹

- VAH Full Length, Adj. Height (E0973) 230.00
- VFH Full Length, Fixed Height Standard
- VAHD Desk Length, Adj. Height (E0973)..... 230.00
- VFHD Desk Length, Fixed Height.....No Charge

MIXED ARMS (EACH) ¹

- RVAHF Right Arm, Full Length, Adj. Height (E0973) 115.00
- RVFHF Right Arm, Full Length, Fixed HeightNo Charge
- RVAHD Right Arm, Desk Length, Adj. Height (E0973) ... 115.00
- RVFHD Right Arm, Desk Length, Fixed HeightNo Charge
- LVAHF Left Arm, Full Length, Adj. Height (E0973) 115.00
- LVFHF Left Arm, Full Length, Fixed HeightNo Charge
- LVAHD Left Arm, Desk Length, Adj. Height (E0973) 115.00
- LVFHD Left Arm, Desk Length, Fixed HeightNo Charge

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap..... Standard
- 1311BK Airline Buckle Seat Positioning Strap..... 50.00

MISCELLANEOUS

INSTALLATION OF BATTERIES

- 22NFBATTERY 22NF Battery (On-Chair)(1 pair)(E2361) . 300.00

BATTERY CHARGER OPTIONS

- 110CHARGER CTE 8-Amp Battery Charger Standard
- LESTR Lester 8-Amp Battery Charger 320.00
- L8900 Omit Battery ChargerLess 80.00

DRIVE WHEEL TYPE

- 1431-3 14" x 3" w/ Foam Filled Inserts Standard

ELECTRONICS

- MK690 MK6i 90-Amp Controller Standard
- SPJPLUS SPJ+ Single Drive Joystick Standard

JOYSTICK MOUNTING HARDWARE

- QL Quad-Link Swingaway Joystick Mount (E1028). 230.00
- LEFT Left handed Joystick MountNo Charge

ACCESSORIES

- VANO2 Oxygen Holder (E2208) 150.00
- HFO2-ML6 HomeFill O2 Holder ML6 Only 150.00
- 1813M6 MK6i Hand-Held Programmer..... 475.00
- 1303 Wheel Locks 45.00

FOOTNOTES

¹ Arms are flip up and adjustable width.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.**

FRONT RIGGINGS

FOOTBOARD

- 90FB Foot Platform.....Standard

CENTER PIVOT STYLE RIGGINGS

- 70NHD 70° Swingaway Footrest 13.5" – 17.5".....165.00 ²
- 70STAPER 70° Tapered Footrest 13" – 17"165.00 ³
- AT5544 Manual Elevating Legrest (E0990)255.00

PIN STYLE RIGGINGS

- PHW93 Swingaway Footrest 13.5" – 19"165.00 ²
- PHAL4A Smartleg Articulating ELR (K0053).....225.00

FOOTPLATE OPTIONS

- 1651 Composite Footplates..... No Charge
- 1350 Extra Large Aluminum Footplates115.00 ⁴
- AT5543 Adjustable Angle Flip-Up Footplates165.00 ⁵
- CTMTFB Footboard for Center Mount150.00

FOOTREST / LEGREST ACCESSORIES

- 1337 Calf Strap27.00
- 1600BK Heel Loops w/ Ankle Straps (Pair)56.00
- IMPACTG Impact Guards51.00
- LPT2 Longer Pivot Slide Tube for 1350.....83.00
- ALPT Longer Pivot and Slide Tube84.00
- ALPT4 Longer Pivot and Slide Tube 4".....84.00

CENTER MOUNT RIGGINGS

- CTMT-F Center Mount w/ Adjustable Knee Angle300.00 ⁶
- CTMT-FJ Junior Center Mount w/ Adj. Knee Angle300.00 ⁷
- SM Small Footplates (5.25" x 7.5")..... No Charge
- MED Medium Footplates (5.50" x 9.5"). No Charge
- LRG Large Footplates (5.90" x 11.5")..... No Charge
- CALF Calf Pads for Center Mount.200.00

5 STAR SERVICE PLAN


5 STAR SERVICE PLAN

- 5STAR1 Invacare 5 Star Service Plan - Guarantees predictable service costs for one year on Invacare Consumer Power and Scooters.150.00 ⁸

FOOTNOTES

- ² Standard with 1651 footplates.
- ³ Standard with impact guards (IMPACTG) and flip up footplates (AT5543).
- ⁴ Not available with AT5544.
- ⁵ One piece footboard for center mount riggings (CTMT-F).
- ⁶ Allows 60 to 90 degrees static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15".
- ⁷ Allows 60 to 90 degrees static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12".
- ⁸ Visit www.invacare.com/ipartnerH for more information.

Specification	Van Seat	Van Seat w/ Solid Seat Pan
Weight capacity	450 lb.	450 lb.
Maximum speed	5 MPH	5 MPH
Incline capability	6°	6°
Per-charge range	Up to 13 mi.	Up to 13 mi.
Turning radius	23"	23"
Ground clearance	3"	3"
Motor Package	4-pole, 650 watt	4-pole, 650 watt
Standard MK6i Electronics	MK690/SPJ+	MK690/SPJ+
Overall Dimensions		
Base length w/o front riggings	35"	35"
Base length w/ front riggings	38.5"	38.5"
Height	38"	38"
with headrest	45" - 47"	45" - 47"
Seat Dimensions		
Width	18" - 24"	18" - 22"
Depth	16" - 22"	16" - 20"
Back height	18"	18"
with headrest	25" - 27"	25" - 27"
Back angle	85° - 115°	85° - 115°
Seat-to-floor height		
Low	21"	21"
Weight		
With batteries	267 lb.	267 lb.
Without batteries	191 lb.	191 lb.
With transport brackets	Add 10 lb.	Add 10 lb.
Battery options	22NF	22NF
Wheel and caster size		
Drive wheels	14" x 3" foam filled	14" x 3" foam filled
Front/rear casters	6" x 2"	6" x 2"
Limited warranties		
Base frame	60 months	60 months
Seat frame	36 months	36 months
Electronics	12 months	12 months
Motor/gearbox	12 months	12 months


Yes, you can.®
 Invacare Corporation
 USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com
 All rights reserved. Trademarks are identified by the symbols TM and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
 Specifications and prices are subject to change without notification.
 ©2008 Invacare Corporation