





Yes, you can:

Invacare® Solace® Resolution Series Pressure-Reducing Mattresses




SRS1080
with raised
side rails (RSR)




SRS1080

This series is the “best” Invacare solution in foam pressure-reducing mattresses for residents at **high risk** of pressure ulcer development. Featuring two layers of the highest-density fortress cut foam, this mattress series offers pressure reduction and therapy in a non-powered mattress.


Invacare® Solace® Resolution Series Pressure-Reducing Mattresses




Body Interface Pressure Imaging Test



Test demonstrates contact measurement between key pressure points. Subject is 5'5" female, 110 lbs Simulation. Individual resident results will vary.



Fortress Cut Foam
provides more surface area, more relief



12° heel slope
protects vulnerable heels

Specifications

- Resident Risk Level : High
- Built-in Side Supports
- Cover : 70 Denier nylon top, 11oz. tri-laminate vinyl bottom, waterproof, vapor permeable
- Mattress Construction : Foam
- Limited Warranty : 5 years on mattress, 1 year on cover

SRS1000 Series

Width.....	36"
Available in 76," 80," 84" Length.....	*
Weight Capacity.....	350 lbs
Available with Raised Side Rails (RSR).....	In 80" Length only

Invacare Continuing Care, Inc. www.invacare-ccg.com info_ccg@invacare.com

USA
One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada
994 Hargrievie Road
London, Ontario
N6E 1P5 Canada
800-668-2337

For more information, contact your Invacare Continuing Care Sales Representative at 1.800.668.2337 or visit www.invacare-ccg.com.

© 2009 Invacare Continuing Care, Inc. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

Form 08.62.9/09



Yes, you can.