

Invacare® CareGuard™ Wheelchairs

It's a dilemma. You're being forced to find the lowest priced wheelchair for your customer. But at the same time, you want to make sure that what you're offering best fits their needs. And finding wheelchairs that meet both criteria hasn't been easy. That is, until now.

Introducing Invacare® CareGuard™ Wheelchairs. A complete line of high quality wheelchairs, priced to meet your cost competitive needs.



Invacare CareGuard wheelchairs range from standard, lightweight, high-strength lightweight and heavy duty. And they're all available through Total One Stop ShoppingSM. That means you get what you're looking for without dealing with a multitude of sales reps. You'll get financing services.* Plus, you'll get a commitment to quality and product innovation. What's more you'll get excellent customer service and strong supply lines, guaranteeing you prompt delivery. You see, with CareGuard Wheelchairs you'll get everything you expect from Invacare.

* Subject to approval by financial services.



**"Invacare says
'Yes, you can™', and
I couldn't have said
it better myself."**

Arnold Palmer
for Invacare



Yes, you can.™

Invacare® CareGuard™ K2000



Features

- Durable, low maintenance; chrome plated carbon steel frame
- Dual axle position for adult or hemi seat-to-floor heights
- Flame retardant, dual embossed black vinyl upholstery is mildew and bacteria resistant
- 24" solid rubber rear tires mounted on composite mag wheel
- Padded armrests
- 250-lb. weight limit
- Limited lifetime warranty on side frames and crossbraces

Model No. Description

Model No.	Description
K26628FR	16", Removable Desk Padded Arm, Swingaway Footrest, Black Vinyl
K26628LR	16", Removable Desk Padded Arm, Swingaway Legrest, Black Vinyl
K28608FR	18", Permanent Full Padded Arm, Swingaway Footrest, Black Vinyl
K28608LR	18", Permanent Full Padded Arm, Swingaway Legrest, Black Vinyl
K28628FR	18", Removable Desk Padded Arm, Swingaway Footrest, Black Vinyl
K28628LR	18", Removable Desk Padded Arm, Swingaway Legrest, Black Vinyl
K28618FR	18", Removable Full Padded Arm, Swingaway Footrest, Black Vinyl
K28618LR	18", Removable Full Padded Arm, Swingaway Legrest, Black Vinyl

Invacare® CareGuard™ K3000



Features

- Durable, low maintenance; chrome plated carbon steel frame
- Dual axle position for adult or hemi seat-to-floor heights
- Flame retardant, black nylon upholstery
- 24" solid rubber rear tires mounted on composite mag wheel
- Padded armrests
- 250-lb. weight limit
- Weighs less than 36 lb.
- Limited lifetime warranty on side frames and crossbraces

Model No. Description

Model No.	Description
K36628FR	16", Removable Desk Padded Arm, Swingaway Footrest, Black Nylon
K36628LR	16", Removable Desk Padded Arm, Swingaway Legrest, Black Nylon
K38628FR	18", Removable Desk Padded Arm, Swingaway Footrest, Black Nylon
K38628LR	18", Removable Desk Padded Arm, Swingaway Legrest, Black Nylon
K38618FR	18", Removable Full Padded Arm, Swingaway Footrest, Black Nylon
K38618LR	18", Removable Full Padded Arm, Swingaway Legrest, Black Nylon
K30628FR	20", Removable Desk Padded Arm, Swingaway Footrest, Black Nylon
K30628LR	20", Removable Desk Padded Arm, Swingaway Legrest, Black Nylon

Invacare® CareGuard™ K4000



Features

- Durable, low maintenance; chrome plated carbon steel frame
- Aluminum crossbraces are painted black
- Flip-back space saver arms
- Adjustable back height
- Dual axle position for adult or hemi seat-to-floor heights
- Flame retardant, black nylon upholstery
- 24" solid rubber rear tires mounted on composite mag wheel
- Padded armrests
- 250-lb. weight limit
- Weighs less than 34 lb.
- Limited lifetime warranty on side frames and crossbraces

Model No. Description

Model No.	Description
K46628FR	16", Removable Desk Padded Arm, Swingaway Footrest, Black Nylon
K46628LR	16", Removable Desk Padded Arm, Swingaway Legrest, Black Nylon
K48628FR	18", Removable Desk Padded Arm, Swingaway Footrest, Black Nylon
K48628LR	18", Removable Desk Padded Arm, Swingaway Legrest, Black Nylon
K48618FR	18", Removable Full Padded Arm, Swingaway Footrest, Black Nylon
K48618LR	18", Removable Full Padded Arm, Swingaway Legrest, Black Nylon

Invacare® CareGuard™ K7000



Features

- Dual axle position
- Heavy duty reinforced frame
- Flame retardant, dual embossed black vinyl upholstery is mildew and bacteria resistant
- Rear 5/8" axle
- 24" solid rubber rear tires mounted on composite mag wheel
- 8" x 2" semi pneumatic caster
- 400 lb. weight capacity
- Limited lifetime warranty on side frames and crossbraces

Model No. Description

Model No.	Description
K72228FR	22", Removable Desk Padded Arm, Swingaway Footrest, Black Vinyl
K72218FR	22" Removable Full Padded Arm, Swingaway Footrest, Black Vinyl
K72228LR	22", Removable Desk Padded Arm, Swingaway Legrest, Black Vinyl
K72218LR	22", Removable Full Padded Arm, Swingaway Legrest, Black Vinyl
K72428FR	24, Removable Desk Padded Arm, Swingaway Footrest, Black Vinyl
K72428LR	24, Removable Desk Padded Arm, Swingaway Legrest, Black Vinyl