

Invacare® Storm Series® 3G Ranger X™ Power Wheelchair

Price List and Order Form Effective October 9, 2006 Revised November 6, 2006



Yes, you can.

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone # _____ P.O.# _____ Date _____

Ship To Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTION

CHAIR TYPE

- 3GRX Ranger X Base & Seat w/MK6i™ Electronics7,569.00
- 3GRXBASE Ranger X Base Only w/MK6i Electronics6,829.00
- 3GSEAT Seat Only740.00

TRANSPORT TIE-DOWN

- TRRO Transport Ready Option.....300.00
- TRBKTS Wheelchair Transport Brackets.....250.00

* TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

TRRO FEATURES FOR 3G RANGER X

Features Required for TRRO	Optional Features for TRRO
<ul style="list-style-type: none"> • Storm Series Ranger X Base & Seat • 400 lb. user weight limit • Short or Long Base Frame • Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle • Any Seat Width • Any Seat Depth • Low Seat-to-Floor Height • Sectional or Straight Rear Uprights • 18" - 24" ASBA Back Height • 90° or 95° Back Angle • Any Arm Type • 8" x 2.25" Semi Pneumatic Casters • Rear Wheel Standard Location • 14" x 3" Tires with Foam Filled Inserts • Any Color Frame Finish • Black Nylon or Black Vinyl Upholstery • Center Pivot Style Bracket • 60° Swingaway, 70° Swingaway or 70° Tapered Footrests 	<ul style="list-style-type: none"> • TrueTrack Upgrade • Fender • Flip-Up Composite Footplates, Aluminum Footplates-Deluxe or Adjustable Angle Flip Up Footplates • Calf Strap • Impact Guards • MK6i MPJ™+, PSR+ / PSF Joysticks • Left Hand Mounted Joystick • "T" Handle Flexible Joystick Ext. • Straight Handle Flexible Joystick Ext. • MK6i Programmer • Reinforced Back Upholstery • Group 24 Batteries

Please note that all TRRO chairs are set at a 10° seat angle

SEAT SELECTION

3 USER WEIGHT LIMITS

- U400 User Weight <400 lb.....No Charge

FRAME TYPE

- SFS Short Base.....Standard
- LFS Long Base.....No Charge

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range.....No Charge
- 2022W 20" to 22" Seat Width Range.....No Charge

4 FRAME WIDTH

- W16 16" Wide.....No Charge
- W17 17" Wide.....No Charge
- W18 18" Wide.....No Charge
- W19 19" Wide.....No Charge
- W20 20" Wide.....329.00
- W21 21" Wide.....329.00
- W22 22" Wide.....329.00
- VW18 18" Wide.....No Charge
- VW20 20" Wide.....329.00
- VW22 22" Wide.....329.00

* TRUETRACK UPGRADE

- TRUTK TrueTrack Upgrade.....1,899.00

SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle...No Charge
- VAN Van Style Captain Seat100.00
- ADJRECL Manual Recline for Adjustable ASBA (E1226)744.00
- CTS Contoura Seating395.00
- OMCU Omit Seat Cushion.....Less 20.00

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range.....No Charge
- 1922D 19" to 22" Seat Depth RangeNo Charge

* TrueTrack upgrade currently not available with MK6i electronics. Please see special MK5 Electronic Options price sheet effective 10/1/06, for MK5 TrueTrack Controller.

FOOTNOTES

- RX comes standard with 400 lbs. user weight capacity, 5 mph maximum speed, Heavy Duty 4 Pole Motor package. SPJ+ Joystick, and MK690 controller. Base only does not come with a seat positioning strap. Base only standard w/MK690ACC controller.
- When ordering a replacement seat and changing depth to 20" or deeper, order new arms.
- U400 user weight uses 4-Pole Heavy Duty Motors with a maximum speed of 5 mph.
- Contoura™ Seating available in even widths 16"-22".
- Manual reclining back not available. Additional cost may apply with less than 16" width/depth.
- Selected seat depth on Contoura assumes 1" of back cushion compression. For example: When choosing an 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The Contoura back cushion is 2" thick. Please see page 5 for seat size offerings.

NOTE: All specifications and dimensions are approximate.

Invacare® 3G Ranger X™ Power Wheelchair

Accessories and Options

1 SEAT DEPTH

- D16 Seat Depth 16" (Short Base)No Charge
- D17 Seat Depth 17" (Short Base)No Charge
- D18 Seat Depth 18" (Short Base)No Charge
- D19 Seat Depth 19" (Long Base)No Charge
- D20 Seat Depth 20" (Long Base)529.00
- D21 Seat Depth 21" (Long Base)529.00
- D22 Seat Depth 22" (Long Base)529.00

2 SEAT-TO-FLOOR HEIGHT

- LOW Seat-To-Floor - LowNo Charge
- MED Seat-To-Floor - Med.....No Charge
- TALL Seat-To-Floor - Tall.....No Charge

3 VAN SEAT STYLE

- SEMIR Semi-Recline **Standard**
- FULLR Full-Recline (E1226)469.00

BACK SELECTION

ADJASBA BACK TYPE

- SECT Sectional Rear Uprights.....No Charge
- SRUP Straight Rear UprightsNo Charge
- 9720P No Push Handles (Straight Canes)82.00
- L8661 Omit Rear Uprights.....Less 20.00

4 BACK HEIGHT RECLINER

- BHR18 Back Height 18"No Charge
- BHR20 Back Height 20"No Charge
- BHR22 Back Height 22"No Charge
- BHR24 Back Height 24"No Charge
- BHR26 Back Height 26"No Charge

BACK HEIGHT

ASBA

- BH12 Back Height 12"No Charge
- BH13 Back Height 13"No Charge
- BH14 Back Height 14"No Charge
- BH15 Back Height 15"No Charge
- BH16 Back Height 16"No Charge
- BH17 Back Height 17"No Charge
- BH18 Back Height 18"No Charge
- BH19 Back Height 19"No Charge
- BH20 Back Height 20"No Charge
- BH21 Back Height 21"No Charge
- BH22 Back Height 22"No Charge
- BH23 Back Height 23"No Charge
- BH24 Back Height 24"No Charge

CONTOURA

- TSBH25 Back Height 25"No Charge
- TSBH295 Back Height 29.5"No Charge

BACK ANGLE DEGREE

- 85BA 95°No Charge
- 90BA 90°No Charge
- 95BA 95°No Charge
- 100BA 100°No Charge
- 105BA 105°No Charge

BACK ACCESSORIES

- L8601 Omit Back Upholstery.....Less 10.00
- RBP Reinforced Back Upholstery No Charge

ARM SELECTION

5 ARM TYPE

- 18F Full Length Fixed Hgt. Conv. Flip Up..... No Charge
- 28F Desk Length Fixed Hgt. Conv. Flip Up..... **Standard**
- 19F Full Length Adj. Hgt. Conv. Flip Up (E0973)230.00
- 29F Desk Length Adj. Hgt. Conv. Flip Up (E0973).....230.00

MK6i SELECTION

MK6i ELECTRONIC CONTROLLERS

- MK690 90 Amp Controller **Standard**
- * MK6TT True Track Controller No Charge

❖ STANDARD DRIVER CONTROL

- MPJM6 MJPJ+ Multiple Drive Joystick..... **Standard**

❖ ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- PSM6 Personal Joystick, Inductive on Front (E2399)750.00
- PSRM6 Personal Joystick, Inductive on Rear (E2399)750.00
- PBOB Push Button On/Off w/PB Drive Select...No Charge
- PBSS Push Button On/Off w/Speed SelectNo Charge
- OR On/Off Mounted RightNo Charge
- OL On/Off Mounted LeftNo Charge
- ASLRDYM6 MK6i Display Only (**No Driver Control**) No Charge

FOR YOUR ASL DRIVE SYSTEM NEEDS

PLEASE REFER TO THE FOLLOWING:

ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

- 1500M6 Prop RIM Head Control w/MK6i Display (E2327) 2,875.00
- SMHD Small HeadrestNo Charge
- LGHD Large HeadrestNo Charge
- 1558MM6 Compact Joystick w/Display (E2399) 1,130.00

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325) 1,500.00
- 5018M6 Wafer Board & Digital Interface Box (E2322)1,550.00
- 5020M6 Mini Tash Joystick & Digital Interface Box (E2321).....1,750.00

SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326) .380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount770.00
- 5018 Wafer Board310.00
- 5020 Mini Tash Joystick405.00
- PACM6 Proportional Attendant Control850.00
- 1552M Digital Attendant Control850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028)279.00
- QLM6 Quad Link Retract. Joystick Mnt MPJ+ Van (E1028) 230.00
- ARM250 Stealth Height Adjustable Joystick Mount (E1028) 210.00
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 (E1028)279.00
- PKG32669 Complete Bib Assembly (Tash/1558) 265.00
- GATMPJ6 Gatlin Midline for MPJ+ (E1028)600.00
- GAT1812 Gatlin Midline for 1812 (E1028)600.00
- RCM Rear Cane Mount for Compact Joystick395.00

FOOTNOTES

1. Contoura Seating available in depths 16"-22".
2. Based on 5° tilt and 16" seat depth on ASBA add 1/4" for additional 2" of seat depth.
3. Used only with Van style seat selection.
4. Used only with manual reclining seat selection.
5. All arms are flip up and detachable. Not available on 12-15" seat widths, must use adjustable height.
6. TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
7. Drive Select through mode switch, or mode switch and Joystick Command.
8. PSRM6 Only.
9. Display only, allows choice from order form.
10. Must also use either a Multiple Drive Joystick or MK6i Display.
11. Requires SNPM6 Sip-N-Puff Interface Box.
12. QLAM6 comes standard with ARM250.
13. Only available with Van Seating.

NOTE: All specifications and dimensions are approximate.

* TrueTrack upgrade currently not available with MK6i electronics. Please see special MK5 Electronic Options price sheet effective 10/1/06, for MK5 TrueTrack Controller.

❖ Includes Mono Port "Y" Cable to provide additional switch input port.

For the most current pricing information go to www.invacare.com.

Invacare® 3G Ranger X™ Power Wheelchair

Accessories and Options

ELECTRONIC ACCESSORIES

<input type="checkbox"/>	1813M6	MK6i Programmer w/Pro Memory Card.....	350.00
<input type="checkbox"/>	1813SD	Pro Memory Card (USB Ready).....	90.00
<input type="checkbox"/>	1560	"T" Handle Flexible Joystick Ext. (E2323)	77.00
<input type="checkbox"/>	1561	Straight Handle Flexible Joystick Ext. (E2323).....	77.00
<input type="checkbox"/>	MPJLHS	MPJ+ Added Left Hand Support Platform.....	49.00
<input type="checkbox"/>	MPJRHS	MPJ+ Added Right Hand Support Platform	49.00
<input type="checkbox"/>	1826	Chin Cup (E2324).....	48.00
<input type="checkbox"/>	A24VPS	24V MK6 Auxiliary Power Source.....	200.00
<input type="checkbox"/>	AUX12M6	Auxiliary Module For 1 & 2	795.00
<input type="checkbox"/>	AUX34M6	Auxiliary Module For 3 & 4	795.00

ASL OPTIONAL ACCESSORIES

<input type="checkbox"/>	EGSBLK	Egg Switch - Black	65.00
<input type="checkbox"/>	ASL304A	Light Wobble Switch w/611 Mount	320.00
<input type="checkbox"/>	ASL202J	Fiber Optic Switch.....	610.00
<input type="checkbox"/>	ASL208J	Adjustable Proximity Switch	550.00
<input type="checkbox"/>	ASL504M6	Remote Emergency Stop Switch.....	535.00

FRAME AND UPHOLSTERY COLOR

FRAME FINISH CONTEMPORARY

<input type="checkbox"/>	24P	Wet Black.....	No Charge
<input type="checkbox"/>	60P	Silver Metallic	No Charge
<input type="checkbox"/>	71P	Silver Vein.....	No Charge
<input type="checkbox"/>	105P	Midnight Blue.....	No Charge
<input type="checkbox"/>	115P	Black Prism.....	No Charge

TOO HOT

<input type="checkbox"/>	30P	Sunny Yellow.....	No Charge
<input type="checkbox"/>	61P	Electric Red.....	No Charge
<input type="checkbox"/>	116P	Sunkissed Orange.....	No Charge
<input type="checkbox"/>	117P	Tangerine Red.....	No Charge
<input type="checkbox"/>	118P	Fruit Fusion	No Charge
<input type="checkbox"/>	119P	Bubblegum Pink	No Charge
<input type="checkbox"/>	120P	Grape Madness	No Charge

WAY COOL

<input type="checkbox"/>	62P	Electric Blue.....	No Charge
<input type="checkbox"/>	104P	Electric Teal.....	No Charge
<input type="checkbox"/>	121P	Lolly Pop Blue	No Charge
<input type="checkbox"/>	122P	Cosmic Blue.....	No Charge
<input type="checkbox"/>	123P	Celestial Purple.....	No Charge
<input type="checkbox"/>	124P	Vivid Violet.....	No Charge
<input type="checkbox"/>	125P	Emerald Green	No Charge
<input type="checkbox"/>	126P	Shamrock Green.....	No Charge
<input type="checkbox"/>	127P	Grasshopper Green.....	No Charge

UPHOLSTERY COLORS

<input type="checkbox"/>	U240	Black Nylon Upholstery.....	Standard
<input type="checkbox"/>	U67	Black Vinyl.....	No Charge
<input type="checkbox"/>	GC	Grey Cloth.....	No Charge
<input type="checkbox"/>	GV	Grey Vinyl.....	No Charge
<input type="checkbox"/>	TV	Tan Vinyl.....	No Charge

CASTER AND WHEEL SELECTIONS

CASTER OPTIONS

<input type="checkbox"/>	1259	6" X 2" Semi-Pneumatic Casters.....	No Charge
<input type="checkbox"/>	1228-3	8" X 2 ¼" Semi Pneu. Casters Lt Grey.....	Standard
<input type="checkbox"/>	1251-3	8" X 2" Pneumatic Casters.....	No Charge
<input type="checkbox"/>	1220-3	9" X 2 ¾" Pneumatic Casters.....	No Charge
<input type="checkbox"/>	1465-3	9" X 2 ¾" Foam Filled Casters.....	No Charge

CASTER SUSPENSION

<input type="checkbox"/>	9213-6	Shock Fork For 6" Casters (E1016).....	119.00
<input type="checkbox"/>	9213-8	Shock Fork For 8" Casters (E1016).....	119.00

WHEEL OPTION

<input type="checkbox"/>	1430-3	14" X 3" Pneumatic Tires	Standard
<input type="checkbox"/>	1431-3	14" X 3" Tires w/Foam Filled Inserts.....	250.00
<input type="checkbox"/>	1440	14" X 4" Pneumatic Tires	200.00
<input type="checkbox"/>	1441	14" X 4" Tires w/Foam Filled Inserts	320.00

FREE WHEEL HUB

<input type="checkbox"/>	FHUB	Free Wheel Hub.....	25.00
--------------------------	------	---------------------	-------

REAR WHEEL LOCATION

<input type="checkbox"/>	RWF2	Rear Wheel Forward 2".....	No Charge
<input type="checkbox"/>	RWF1	Rear Wheel Forward 1".....	No Charge
<input type="checkbox"/>	RWS	Rear Wheel Standard	Standard

FENDER OPTION

<input type="checkbox"/>	L8603	Omit Fender	Less 50.00
--------------------------	-------	-------------------	------------

POSITIONING STRAPS

SEAT POSITIONING STRAPS

<input type="checkbox"/>	1311BK	Airline Buckle Seat Positioning Strap.....	50.00
<input type="checkbox"/>	1515	Auto Style Buckle Seat Positioning Strap	Standard

CHEST POSITIONING STRAPS

<input type="checkbox"/>	7311BK	Airline Buckle Chest Positioning Strap (E0960) ..	50.00
<input type="checkbox"/>	9515	Auto Style Chest Positioning Strap (E0960)	50.00
<input type="checkbox"/>	7321BK	Hook & Loop Chest Positioning Strap (E0960) ..	40.00

MISC ACCESSORIES

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	1303	Wheel Locks	45.00
<input type="checkbox"/>	VAN02	O ₂ Holder - Van (E2208)	150.00
<input type="checkbox"/>	SEATO2	Oxygen Holder (ASBA Seat Only) (E2208).....	150.00
<input type="checkbox"/>	9205-3	Suspension Shocks (E1016).....	400.00

BATTERY SELECTION

BATTERY BOXES

<input type="checkbox"/>	GP24	Group 24 Battery Boxes	No Charge
--------------------------	------	------------------------------	-----------

BATTERIES

<input type="checkbox"/>	24BATTERY	24 Gel Battery (on chair) (E2363).....	410.00
--------------------------	-----------	--	--------

FOOTNOTES

- 1 . Must add 24V Auxiliary Power Source.
- 2 . Not available with short base.
- 3 . Not available with 9" casters.
- 4 . 14" x 4" wheels increase overall width of the base by 2 1/2". No fender included. Not available with short base.
- 5 . Does not include StormSeries Mag Wheel. Free Wheel Hub increases overall width 3/4".
- 6 . Not available with pneumatic wheels.
- 7 . This location changes the center-of-gravity and performance. Forward wheel position will make turning and climbing easier, but more prone to "wheely".

NOTE: All specifications and dimensions are approximate.

Invacare® 3G Ranger X™ Power Wheelchair
Accessories and Options

HANGER BRACKET SELECTION

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket..... No Charge
- PSB Pin Style Bracket..... No Charge
- NHB No Hanger Bracket..... No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

CENTER PIVOT

- 60HD 60° Swingaway Footrests..... No Charge
- 70NHD 70° Swingaway Footrests..... **Standard**
- 70STAPER 70° Tapered Footrests..... No Charge
- AT5544 Elevating Legrests (E0990) 255.00

PIN STYLE

- PAS4A Smartleg® Art. Legrests 17 ¼" -21 ½" (K0053) 320.00
- PW93 Swingaway Footrest 13 ½" - 19"..... No Charge
- P904A Swingaway Elevating Legrests 16" - 19 ¾" (E0990) 255.00
- PAL4A Smartleg Art. Legrests 17 ¼" - 21 ½" (K0053) 320.00

POWER LEGS CONTROL MODULE

- 3 ELRPW Power ELR w/Articulation (pair) (E1010)..... 1,895.00
- 4 4WSB Four Way Switch Box..... No Charge
- S4WSB Multiple Actuator Interface Box (E2311)..... 2,600.00
 - FWT 4 Way Toggle..... No Charge
 - QPB 4 Push Buttons..... No Charge

FOOTPLATE OPTIONS

- 1651 Flip-Up Composite Footplates..... No Charge
- 1350 Extra Large Aluminum Footplates..... No Charge
- 1450 Aluminum Footplates-Deluxe..... No Charge
- S1263 One Piece Wooden Footboard..... No Charge
- AT5543 Adjustable Angle Flip Up Footplates 165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) 27.00
- IMPACTG Impact Guards 51.00
- A9342 "H" Heel Strap (K0039) 51.00
- LPT2 Longer Pivot Slide Tube For 1350 83.00
- ALPT Longer Pivot And Slide Tube..... 84.00
- ALPT3 3" Extension For 90° Footrest 31.00
- ALPT4 Longer Pivot And Slide Tube 4" 84.00

FOOTNOTES

1. For seat widths 16".
2. Independent powered legrests include elevating and articulating legrests with angle adjustable footplates. Extension range of 13"-20".
3. Required with all multiple actuator systems. Provides a D-9 port for the 4-way toggle or 4 push buttons.
4. Replaces the Four Way Switch Box and allows multiple actuator operation through the driver control.

NOTE: All specifications and dimensions are approximate.

SEATING & POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING**

INFINITY® SEAT AND BACK CUSHIONS
SEE PRICE LIST FORM # 01-351.

CONVENTIONAL PRODUCTS
SEE PRICE LIST FORM # 95-82.

**ABSOLUTE™, LAP TRAYS, AND
PERSONAL BACK PRODUCTS**
SEE PRICE LIST FORM # 01-207.

Invacare® 3G Ranger X™ Power Wheelchair

Accessories and Options

Chair Features and Options
<ul style="list-style-type: none"> ♦ Independent Monroe® shock suspension for clinical and functional comfort ♦ 5 MPH maximum speed - standard ♦ Basic switch option capability on integrated seat (ASBA) ♦ Programmable performance to meet individuals needs ♦ 17 1/2" seat-to-floor on standard, 19 3/4" & 21" optional ♦ Seat location system allows for balanced weigh distribution ♦ Continuous tilt seat 0° - 10° ♦ Dynamic anti tippers ♦ Rigid seat pan ♦ Convenient rear battery access ♦ Group 24 battery boxes ♦ Arms are flip back/removable ♦ Adjustable depth rigging attachment 2" {1" increments} ♦ Automatic dual mode battery charger (Lester brand) ♦ Compact size: 25" overall base width w/14" x 3" wheels ♦ Weight limit - Up to 400 lb. ♦ Optional long base

Ranger X Specifications		
* Overall Width	25" 27"	w/14" x 3" wheel w/14" x 4" wheel
** Overall Length	29.5" 32.5"	Short Long
***Seat-to-Floor	17 1/2" Standard 19 3/4" Optional 21" Optional	Low Medium Tall
Battery	Group 24 deep cycle; gel/sealed	
{Not supplied}		
Warranty	Lifetime limited on base frame 3 year limited seat frame 1 year limited on electronics 18 month limited warranty on gear motors	
Weight Limits	400 lb. inclusive of add-on seat systems and equipment	
Weight w/o batteries	153 lb.	
Shipping Weight	266 lb. {18" x 18" seat and Group 24 batteries}	
* Measurement based on wheel based on wheel base. Seat width is width picked o\plus 7.5" for armrest and joystick.		
**Measured from headtube to anti tipper.		
***Based on 5° tilt		
Note: All specifications and dimensions are approximate.		

Footrest Extension Range			
	AT5043	AT5543	1651
60HD	N/A	12 1/4" - 16 1/2"	13 1/2" - 17 1/2"
70MFX	6 3/4" - 12 1/2"	4 1/2" - 10"	N/A
70NHD	N/A	10 3/4" - 15 3/4"	11 3/4" - 17"
70STAPER	N/A	12 1/2" - 16 1/2"	13 1/4" - 17 1/4"
Measured from Seat Pan to Back of Footplate			

Weight Limit 400 Pounds
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.
NOTE: All weight limits are inclusive of add-on equipment and seat systems.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*	*	X	X	X	X	X
*Not available with Contoura Seat Cushion.							
Back Heights	25" or 29.5"						

Invacare Corporation www.invacare.com



USA
 One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

Invacare, Arrow, Infinity, Monroe, Smartleg, Storm Series, Yes, you can. and the Blue Medallion Design are registered trademarks and Absolute, Contoura, GB, MK6i, MPJ, SPJ, and Ranger X are trademarks of Invacare Corporation. Therafin is registered trademark of W/R Medical Eelectronics Co. Monroe is a registered trademark of Tenneco Automotive, Inc. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Specifications and prices are subject to change without notification.
 © 2006 Invacare Corporation