

Sleep Therapy



Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 03-031. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2008 Invacare Corporation, Elyria, Ohio. Invacare, the Medallion Design and Yes, you can. are registered trademarks of the Invacare Corporation. Polaris and Twilight are trademarks of Invacare Corporation. Information is subject to change without notification.

TABLE of CONTENTS

Description	Page Number
CPAP Polaris TR Accessories	15
CPAP Polaris TR System	13
ISP Polaris EX - ISP3000 & ISP4000	11
Nasal CPAP System - Accessories - Polaris/ Polaris LT	5
Nasal CPAP System - Polaris/ Polaris LT	3
Sleep Therapy Accessories	9
Sleep Therapy Masks	7

Nasal CPAP System Polaris/Polaris LT

Model ISP9800 - Polaris (Shown)

Model ISP9900 - Polaris LT

Model ISP9900AW

Mask and Headgear Models

ISP9800MKS

ISP9800MKL


ISP9900MKS

ISP9900MKL

ISP9800IQ

ISP9900IQ

ISP9900AWIQ




Nasal CPAP System Polaris/Polaris LT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	F	1082994-NLA	Cover with Lens, ISP9800 - Top	1	1
	A,G	1082429	Shroud Set, Polaris	1	1
1		1094568	Cover with Lens, ISP9900 - Top	1	1
2		1082995	Filter	3	1
	A,H	1082428	Filter, Fine	3	1
3	C	2001831	Cord, AC Power, 3 Prong	1	1
4	A,D	1104153	Cord, AC Power, 2 Prong	1	1
5		2002930	Bumper Feet, Rubber	4	1
6		1085420	Tube with Cuff, Smooth Bore, Clear, Packaged (6 FT. L)	1	1
	A,B	1100553	Tubing, Smooth Bore w/Cuff (12" L)	1	1
7	E	1119677	Case, Carrying	1	1
	A	1075698-NLA	Insert, Carton	1	1
	A	1075699	Carton	1	2
	A	1075696	Owner's Manual	1	1
	A	1082997	Instructions, System Set-Up	1	1
	A	1082999	Fuses, T 500MA	2	1
	A	ISP9800IQ	ISP9800 with IQ50106H Mask Assembly	1	1
	A	ISP9900IQ	ISP9900 with IQ50106H Mask Assembly	1	1
NOTE: A - Not Shown. B - To be used with ISP9700-CPAP Poseidon Humidifier. C - Before 11/3/00. D - After 11/2/00. E - Use with 9900AW. F - Individual item no longer available, refer to Part #1082429 Shroud Set with Lens for replacement. G - Includes Top, Middle, Bottom Shrouds and Lens. H - Filter is white in color and located on inside of unit.					

Nasal CPAP System - Accessories

Polaris/ Polaris LT




Nasal CPAP System - Accessories

Polaris/ Polaris LT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	ISP9910	Kit, Manometer	1	1
2		1092792	Manometer	1	1
3		1123914	Package, Pressure Line Adapters	10	1
4		1140387	Tube, Silicon (10-1/2")	1	1
5	B	16501-NLA	Mirage Mask Assembly, Standard Size	1	1
5	B	16502-NLA	Mirage Mask Assembly, Large Size	1	1
6		16504-NLA	Tubing, Swivel	10	1
7		16526-NLA	Forehead Assembly	1	1
8		16535	Arm, Adjustable	5	1
9		16514	Pad, Forehead	10	1
10		16512	Plug, Port	10	1
11		16505	Mask, Frame	1	1
12		16506	Cushion, Nasal - Standard	1	1
12		16507	Cushion, Nasal - Large	1	1
13		16510-NLA	Vent, Air, Bumped	10	1
14		IQ50160LS-NLA	Mask with Tubing and Instructions	1	1
	C	IQ50160S	IQ Mask with Tubing	1	1
	C	IQ55215	Headgear, Sleepnet One-Size Strap	1	1
	C	IQ55121	Headgear, Sleep Net - Medium Cap	1	1
	C	IQ55122	Headgear, Sleep Net - Large Cap	1	1
15		IQ55352	Headgear, Sleep Net - Large Strap	1	1
	C	1082428	Filter, Fine	3	1
	C	IQ55024	Tubing, Sleep Net Elbow with Swivel Adaptor	1	1
	C	1085421	Tubing, Smooth Bore, Clear (6' L)	1	1
	C,D	1100553	Tubing, Smooth Bore w/Cuff (12" L)	1	1
NOTE: A - Includes Items 2-4. B - Mirage Mask Assembly Can Be Used With Other Manufactures of CPAP and Bilevel Units. C - Not Shown D - To be used with ISP9700-CPAP Poseidon Humidifier.					

Sleep Therapy Masks


Twilight Full Face Mask
Shown for
Reference Only



Mask Shown for
Reference Only




Twilight Mask Shown
for Reference Only



8



9







10



Sleep Therapy Masks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1133279	Hose, Swivel, Twilight Series	1	1
2		1133284	Plug, Port - Twilight Series	1	1
3		16902	Pads, Forehead, Pink	1	1
4		16917	Forehead Arm Set (Includes: T-Bar, Forehead Pads, and Forehead Arms)	1	1
5		16906	Axial, Swivel	1	1
6		16901	Plug, Port	10	1
7		1133280	Ring, Retainer - Twilight II	1	1
8		1122734	Kit, Chin Support - Twilight	1	1
			Strap, Chin Assembly - Twilight	1	1
			Instruction Sheet	1	1
9		1122733	Kit, Headgear - Twilight Series	1	1
			Mask, Headgear Assembly	1	1
			Clip, Headgear Mask	1	2
			Instruction Sheet	1	1
10	A	16015	Support, Chin	1	1
	B	IQ55046	Headgear, Sleepnet	1	1
NOTE: A - Used With: ISP9800, 21001, 21031 and 21042. B - Not Shown.					

Sleep Therapy Accessories




Sleep Therapy Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		14909	Air Filter Assembly (Includes: Filter Holders, Primary, and Secondary Filters)	1	1
2	A	14911	Connector, Air Tubing	10	1
3	A	16015	Support, Chin	1	1
4		14928	Bag, Travel, CPAP	1	1
4		15907	Bag, Travel, VPAP	1	1
5		1085420	Tube with Cuff, Smooth Bore, Clear, Packaged (6 FT. L)	1	1
	B	21801	Operating Manual - Sullivan V, HM	1	1
	B	21802	Operating Manual - Sullivan Comfort	1	1
	B	21805	Clinical Manual - Sullivan V, HM	1	1
	B	21807	Clinical Manual - Sullivan Comfort	1	1
	B	1123914	Package, Pressure Line Adapters	10	1
	B	IRC840	Adapter, AC, 115V	1	1
NOTE: A - Used With: ISP9800, 21001, 21031 and 21042. B - Not Shown.					

ISP Polaris EX

ISP3000 & ISP4000




ISP Polaris EX

ISP3000 & ISP4000

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,G	ISP4000S	Humidifier, Heated Polaris EX, Silver	1	1
1	C	1132517	Cap, Reservoir	1	1
2	C	1132522	Reservoir, Water	1	1
3	C,H	1136073	Elbow and Reservoir Assembly	1	1
4	C	1117424	Elbow	1	1
5	C,J	1117423-NLA	Plate, Bottom, Assembly	1	1
5	C,K	1132523	Bottom Plate Assembly	1	1
6	C	1132521	O-Ring	1	1
7	B,I	1136074	Humidifier Reservoir Assembly	1	1
8	C	1120165	Coupling, Hose	1	1
9	C	11008X001	Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	1
10	B	1124211	Replacement Filters, ISP3000	3	1
11	B	1124213	Door, Access with Screw	1	1
12	B	1023281	Screw, Sheet Metal (#8)	1	1
13	B	1119677	Case, Carrying	1	1
	A,D	1085420	Tube with Cuff, Smooth Bore, Clear, Packaged (6 FT. L)	1	1
	A,C	1119768	Insert, Bottom Humidifier	1	2
	A,C	1119769	Carton, Shipper	1	1
	A,B	1119765	Insert, Bottom	1	1
	A,B	1119766	Insert, Top	1	1
	A,B	1119767	Carton, Shipper	1	1
	A	1110559	Owner's Manual, ISP3000	1	1
	A	1118358	Owner's Manual, ISP4000	1	1
14	A,B,L	1141773	Kit, Filter Ultra Fine	1	1
			Filter, Ultra Fine CPAP	3	1
			Instructions, Ultra Fine Filter	1	1
NOTE: A - Not Shown. B - For use on ISP3000. C - For use on ISP4000. D - Includes Part #1075695 Smooth Bore Tube with Cuff enclosed in Polybag. G - Includes Silver Humidifier Base, Silver Humidifier Shroud, Hot Plate Assembly, Hose Coupling Kit, Shipper Carton and Owners Manual. H - Includes Water Reservoir, Cap and Elbow. I - Includes Reservoir, Cap, Bottom Plate, O-Ring and Elbow. J - Plate used with old style Water Reservoir, Before 8/17/05. K - Used with new Water Reservoir, after 8/16/05. L - Can be used in addition to CPAP Filter, for additional filtration.					

CPAP Polaris TR System

For Technical and Service Issues with the Polaris TR please contact the Manufacturer at 866-406-0595.
This product is distributed by Invacare Corporation.



CPAP Polaris TR System

For Technical and Service Issues with the Polaris TR please contact the Manufacturer at 866-406-0595.
This product is distributed by Invacare Corporation.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		ISP3510	Polaris TR CPAP Assembly	1	1
2	C	1145255	Lid, Humidifier - Polaris CPAP TR	1	1
3	A	ISP4510	Heated Humidifier Assembly, Polaris TR	1	1
4		1145254	Cell, Battery, Polaris CPAP TR	1	1
5	D	1145260	Charger, Battery - Polaris CPAP TR	1	1
	B	1145452	Patient Manual - Polaris CPAP TR	1	1
	B	1145453	Clinician Manual - Polaris CPAP TR	1	1
	B	1145454	Interface Manual, Polaris CPAP TR	1	1
	B	1145475	Compliance Quick Guide - Polaris CPAP TR	1	1
NOTE: A - Includes Humidifier Shell Assembly and Humidifier Lid. B - Not Shown. C - Includes the Humidifier Lid Only. D - Includes Items 4 and 5.					

CPAP Polaris TR Accessories

For Technical and Service Issues with the Polaris TR please contact the Manufacturer at 866-406-0595.
This product is distributed by Invacare Corporation.

1 	2 	3 
4 	5 	6 
7 	8 	9 
10  Item 10 contains the Shell of each item only. These are non-functional units for display only.		

CPAP Polaris TR Accessories

For Technical and Service Issues with the Polaris TR please contact the Manufactuer at 866-406-0595.
This product is distributed by Invacare Corporation.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1145472	Demo Card	1	1
2		1145458	Docking Station - Polaris CPAP TR	1	1
3		1145460	Reader, Card	1	1
4		1145459	Data Card - Polaris CPAP TR	1	1
5	A	1145257-NLA	Kit, Filter	1	1
6		1145259	Bag, Travel - Polaris CPAP TR	1	1
7	B	ISP2510M	Seal, Nasal - Twilight NP (Medium)	1	1
7	B	ISP2510L-NLA	Seal, Nasal - Twilight NP (Large)	1	1
7	B	ISP2510XL	Seal, Nasal - Twilight NP (X-Large)	1	1
8		1145473	Pouch, Carry - Twilight NP	1	1
9		1145261	Headgear Assembly - Twilight NP	1	1
10	E,H	1145455	Shell, CPAP Humidifier Battery - Polaris TR	1	1
	C	1145256	Power Supply - Polaris CPAP TR	1	1
	C,D	1145461	Software Package - Polaris CPAP TR	1	1
	C	1145474	Manual, Compliance System - Polaris CPAP TR	1	1
	C,F	ISPPCS-B	Compliance System B - PolarisTR	1	1
	C,G	ISPPCS	Kit, Compliance - Polaris TR	1	1
NOTE: A - Includes (2) Replacement Filters. B - Note which size you need before placing order. C - Not Shown. D - Includes Installation Disk. E - Includes the shell only of the Polaris CPAP TR, Humidifier Shell and Battery Shell. These are non-functional units for display use only. F - Includes Items 3, Part # 1145461 Software Package, and Part #1145474 Compliance Manual. G - Includes ISPPCS-B, Item 2, 4 and Manual. H - These items are no longer available seperately, you must purchase these together.					

Invacare Corporation

www.invacare.com



Yes, you can. TM

One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Florida Technical Support
407-321-5630

Customer Service
800-333-6900

© 2008 Invacare Corporation, Elyria, Ohio.
All rights reserved. Trademarks are
identified by the symbols TM and ®. All
trademarks are owned by or licensed to
Invacare Corporation unless otherwise
noted. Information is subject to change
without notification. (Revised 10/17/08)

Form No. 03-031