



Frames




**1. MVP
Style**

Rear




**2. Jr.
Xtra**

Front



Swingaway

**3. MVP, Xtra,
Jr.**



Fixed

**4. Style (After 3/92)
MVP**

Frames

#1 Rear Frames - MVP/Style

When replacing rear frames for an MVP or Style Chair models built before 7/5/00, 2 frames must be ordered along with new rear canes. Armrest hardware may also be needed if the chair has 2 point arms.

Description	Part Number
Rear Frame (Curved) Before 7/5/00 Right See Note(s) A	1024766-NLA
Rear Frame (Curved) Before 7/5/00 Left See Note(s) A	1024767-NLA
Rear Frame (Curved) After 7/4/00 Right/Left	8881094943
NOTE: A - When replacing rear frames for an MVP or Style Chair models built before 7/5/00, 2 frames must be ordered along with new rear canes. Armrest hardware may also be needed if the chair has 2 point arms. Specify frame finish when ordering.	

#2 Rear Frames - Xtra/Jr.

Description	Part Number
Rear Frame (Straight)	8881018543
Specify frame finish when ordering.	

#3 Swingaway Front Frames - Xtra/Jr./MVP

Due to aesthetics, swingaway front frames for an MVP, Xtra, or Jr. chair built Before 7/5/00 must be replaced in pairs.

Right/Left	12" Deep	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep
Right	1018171-NLA See Note(s) A	8881018533-NL A See Note(s) A	8881018535-NL A	8881018537	8881018539-NL A	8881018541
Left	8881018172-NL A See Note(s) A	8881018534-NL A See Note(s) A	8881018536-NL A	8881018538	8881018540-NL A	8881018542
NOTE: A - Jr. chairs use Swingaway Front Frames for both Swingaway and Fixed Front chairs. See 90 Degree Fixed Supports category. Specify frame finish when ordering.						

#4 Fixed Front Frames - MVP/Style

Front Frame Type	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep
60 Degree	8881068512-NLA	8881068513-NLA	8881068514-NLA	8881068515-NLA	8881068516-NLA
80 Degree	8881038779-NLA	8881038780-NLA	8881038781-NLA	8881038782-NLA	8881038783-NLA
Specify frame finish when ordering.					