


Back Upholstery



Back Upholstery

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B,C	8881073550-NL A	NLA - Upholstery, Back with Foam Assembly	1	1
1	A,B,D	8881075128-NL A	NLA - Upholstery, Back with Foam Assembly	1	1
1	A,B,E	8881076323-NL A	NLA - Backrest Assembly	1	1
2	A,C	8881073654-NL A	NLA - Upholstery with Foam	1	1
2	A,D	8881075127-NL A	NLA - Upholstery with Foam	1	1
2	A,E	8881076322-NL A	NLA - Upholstery with Foam	1	1
3		1073655-NLA	NLA - Liner, Mesh	1	1
4		1105393-NLA	NLA - Kit, Hardware (Includes: 4 Screws and 4 Washers)	4	1
NOTE: A - Specify upholstery color, see chart on page titled "Packaging & Literature" B - Includes items 2-4 C - CP200 D - CP300 ,CP400, CP500 E - CP300, CP400, CP500, CAL133 Uph.					