


Seat and Ottoman Upholstery



Seat and Ottoman Upholstery

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B,C	8881073551-NL A	NLA - Seat/ Ottoman Upholstery Assembly	1	1
1	A,B,D	8881075137-NL A	NLA - Seat/ Ottoman Upholstery Assembly	1	1
1	A,B,E	8881076332-NL A	NLA - Seat/ Ottoman Upholstery Assembly	1	1
2	C	1074000-NLA	NLA - Cushion Assembly with Bag	1	1
2	E	1073657-NLA	NLA - Cushion with Bag	1	1
3		1073661-NLA	NLA - Liner, Ottoman Mesh	1	1
4	A,C	8881073660-NL A	NLA - Upholstery, Seat/Ottoman w/ Foam	1	1
4	A,D	8881075136-NL A	NLA - Upholstery, Seat/Ottoman w/ Foam	1	1
4	A,E	8881076331-NL A	NLA - Upholstery, Seat/Ottoman w/ Foam	1	1
NOTE: A - Specify upholstery color, see color chart on page titled "Packaging & Literature" B - Includes items 2-4 C - CP200 D - CP300, CP400, CP500 E - CP300, CP400, CP500, CAL133 Uph.					