



Clamp-On Style Wheel Locks

Tracer SX Recliner



Clamp-On Style Wheel Locks Tracer SX Recliner

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1080332-NLA	NLA - Wheel Lock Assembly, Black - Right	1	1
1	A	1080333-NLA	NLA - Wheel Lock Assembly, Black - Left	1	1
2		60320X000	Locknut (1/4-28)	1	3
3		60275X000	Locknut (1/4-20)	1	1
4		1010398	Washer, Nylon (1/4 x 1/2 x 1/16")	1	7
5		1080330	Bracket, Wheel Lock, Black - Right	1	1
5	C	1080331-NLA	NLA - Bracket, Wheel Lock, Black - Left	1	1
6		40157X016-NLA	NLA - Washer (1/4 x 1/2 x 1/16")	1	4
7		1015443-NLA	NLA - Shoe	1	1
8		1016592-NLA	NLA - Link, Chrome	1	1
9		1016588-NLA	NLA - Handle, Push, Adult, Chrome - Right	1	1
9		1016589-NLA	NLA - Handle, Push, Adult, Chrome - Left	1	1
10		47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
11		60274X012-NLA	NLA - Screw, Hex Head (1/4-28 x 3/4")	1	3
12		1007017-NLA	NLA - Screw, Hex Head (1/4-20 x 1-1/4")	1	1
13	B	1030987	Hardware, Wheel Lock Attaching	1	1
14			Bolt, Lug (5/16-24 x 1-3/8")	1	2
15			Locknut (5/16-24)	1	2

NOTE: A - Includes items 2-12.
B - Includes items 14 and 15.
C - Chrome is No Longer Available. Order item #1 1080333 in Black. Recommended to order in quantity of 2 to match.