


Seat - Pail - Splash Shield



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1017773-NLA	NLA - Seat, Padded Commode (16")	1	1
1		1123718	Padded Seat Assembly, Commode (16")	1	1
1	B	1028042-NLA	NLA - Seat, Padded Commode (18")	1	1
2		1025378	Hinge, Clip Seat	1	8
3		7930344-NLA	NLA - Package, Flat Head C-Sunk Screw (#10-32 x 1")	25	1
4		1029572-NLA	NLA - Pad, Rubber, Black	1	8
5	C	6319	Shield, Splash	12	1
6	C	6317	Package, Pail w/ Lid w/o Handle	12	1

NOTE: A - Model 6490.
 B - Model 6690.
 C - Accessory-12 Pack.