


Amputee Attachment - 9357L

DLX, SX, Tracer Titan



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1057223-NLA	NLA - Amputee Attachment (Pair)	1	1
2			Screw, Hex Head (1/4-20 x 1-1/2")	1	2
3			NLA - Screw, Hex Head (7/16-20 x 2-1/2")	1	2
4			Locknut (1/4-20)	1	2
5			Locknut (7/16-20)	1	2
6			NLA - Bracket, Amputee, Chrome - Right Adult/Left Hemi	1	1
6			NLA - Bracket, Amputee, Chrome - Left Adult/Right Hemi	1	1
			NLA - Instruction Sheet	1	1
7	B	1023618	Insert, Guide Tube, Black (1" Long)	1	1

NOTE: A - Includes items 2-6 and an instruction sheet.
 B - Not included in the Amputee Attachment hardware package.