


Footboard Assembly

After 11/05/01



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1099714-NLA	NLA - Footboard Assembly Includes: 1 ea. Footboard with Rubber Pad 1 ea. Adjustment Tube, Footboard 2 ea. Washer (8.2 x 15.8 x 3) 1 ea. Safety Pin 1 ea. Cable Tie 1 ea. Nut (M8 x 1.25 x 13) 1 ea. Screw, Hex Head Cap (M8 x 1.25 x 25) 2 ea. Plug Buttons, Ball 1 ea. Screw, Hex Head Cap (M8 x 1.25 x 60)	1	1
2		1099715-NLA	NLA - Pin, Safety	4	1